# Optima Turn

## User Manual

Be well, Live well  
**Wellell**

# CONTENTS

|                                  |     |
|----------------------------------|-----|
| <b>ENGLISH</b>                   |     |
| USER MANUAL .....                | 1   |
| <b>ESPAÑOL</b>                   |     |
| MANUAL DE INSTRUCCIONES .....    | 22  |
| <b>PORTUGUÊS</b>                 |     |
| INSTRUÇÕES DE UTILIZAÇÃO .....   | 51  |
| <b>FRANÇAIS</b>                  |     |
| INSTRUCTIONS D'UTILISATION ..... | 78  |
| <b>DEUTSCH</b>                   |     |
| GEBRAUCHSANLEITUNG.....          | 104 |
| <b>NEDERLANDS</b>                |     |
| MANUALE DI ISTRUZIONI.....       | 130 |
| <b>ITALIANO</b>                  |     |
| GEbruikersHANDLEIDING .....      | 155 |

MODEL NO.: 9P-085500

PLEASE READ THE MANUAL BEFORE USE

# IMPORTANT SAFEGUARDS

## READ ALL INSTRUCTIONS BEFORE USING

### **DANGER- To reduce the risk of electrocution:**

- (1) Always unplug this product immediately after using.
- (2) Do not use while bathing.
- (3) Do not place or store this product where it can fall or be pulled into a tub or sink.
- (4) Do not place in or drop into water or other liquid.
- (5) Do not reach for a product that has fallen into water. Unplug immediately.

### **WARNING- To reduce the risk of burns, electrocution, fire, or injury to persons:**

- (1) Evaluate patients for entrapment risk according to protocol and monitor patients appropriately.
- (2) The intended operator who can touch/operate the device should be someone other than the patient. Never touch the patient and the device at the same time during use.
- (3) The product may be used with caution for patients with spinal injury, but suggested to consult with physician before use. However, it should not be used for patients with unstable spinal fractures.
- (4) Close supervision is necessary when this product is used on or near children. Electrical burns or choking accident may result from a child swallowing a small part detached from the device.
- (5) Use this product only for its intended use as described in this manual. Do not use other mattress not recommended by the manufacturer.
- (6) Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the product to your supplier or Wellell for examination and repair.
- (7) Keep the cord away from heated surfaces.
- (8) Never block any air openings of this product or place it on soft surfaces, such as a bed or couch, where openings may be blocked. Keep the air opening free of lint, hair, and other similar particles.
- (9) Never drop or insert any object into any opening or hose.
- (10) Do not modify this equipment without authorization of the manufacturer.
- (11) Mattress covers have passed skin sensitization and skin irritation test. However, If you suspect that you may have had or are having an allergic reaction, please consult a physician immediately.
- (12) Do not leave long lengths of tubing around the top of your bed. It could lead to strangulation.
- (13) Third conductor in the POWER SUPPLY CORD is only a functional earth.
- (14) The Operator cannot operate the instrument and touch the patient at the same time.
- (15) If the RTC battery needs to be replaced, please contact local representative or place of purchase.

-  **CAUTION:** If there is a possibility of electro-magnetic interference with mobile phones, please increase the distance (3.3m) between devices or turn off the mobile phone.

#### NOTE, CAUTION AND WARNING STATEMENTS:

-  **NOTE:** Indicate some tips.
-  **CAUTION:** Indicate correct operating or maintenance procedures in order to prevent damage to or destruction of the equipment or other property.
-  **WARNING:** Calls attention to a potential danger that requires correct procedures or practices in order to prevent personal injury.

#### SYMBOLS

|                                                                                                 |                                                                                                                                                                                                                                                                                                                                                                                                                   |                                                                                     |                                                              |
|-------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------|
|                | Authorized representative in the European community.                                                                                                                                                                                                                                                                                                                                                              |    | Temperature Limitation                                       |
|                | Manufacturer                                                                                                                                                                                                                                                                                                                                                                                                      |    | Dry clean, Any Solvent Except Trichloroethylene              |
|                | Complies with standards protecting against electric shock for type BF equipment.                                                                                                                                                                                                                                                                                                                                  |    | Do Not Iron                                                  |
|                | Consult operating instructions for use                                                                                                                                                                                                                                                                                                                                                                            |    | Tumble Dry, Normal, Low Heat                                 |
|               | Attention, should read the instructions.                                                                                                                                                                                                                                                                                                                                                                          |   | Do Not Tumble Dry                                            |
|              | Attention, should read the instructions.                                                                                                                                                                                                                                                                                                                                                                          |  | Do Not Bleach                                                |
|              | Functional earth. (for UL only)                                                                                                                                                                                                                                                                                                                                                                                   |  | Do Not Dry Clean                                             |
|              | Class II                                                                                                                                                                                                                                                                                                                                                                                                          |  | Machine wash, regular / normal, 95 degrees C (203 degrees F) |
| <b>IP21</b>  | Protected against solid foreign objects of 12.5 mm and greater; Protection against vertically falling water drops                                                                                                                                                                                                                                                                                                 |                                                                                     |                                                              |
|              | Attention – Observe proper Disposal of Electrical & Electronic Equipment (WEEE): This product should be handed over to an appropriate collection point for the recycling of electrical and electronic equipment. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product. |                                                                                     |                                                              |

# 1. Introduction

This manual should be used for initial set up of the system and for reference purposes.

## 1.1. General Information

The system is a high quality mattress system suitable for treatment and prevention of pressure injury.

The system has been tested and successfully approved to the following standards:

IEC/EN 60601-1  
IEC/EN 60601-1-2  
IEC/EN 61000-3-2 Class A  
IEC/EN 61000-3-3  
CISPR 11 Group 1, Class B

## EMC Warning Statement

This equipment has been tested and found to comply with the limits for medical devices to the EN 60601-1-2. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device.
- Increase the separation between the equipment.
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) are connected.
- Consult the manufacturer or field service technician for help.

## 1.2. Intended Use

This product is intended to help and reduce the incidence of pressure injury while optimizing patient comfort. It also provides following purposes:

- for long term home care of patients suffering from pressure injury.
- for pain management as prescribed by a physician.

The product can only be operated by personnel who are qualified to perform general nursing procedures and have received adequate training in knowledge of prevention and treatment of pressure injury.


**NOTE: Equipment is not suitable for use in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide.**

## 2. Product Description

Unpack the box to check if there are any damage which may have occurred during shipment. Please contact the local representative or place of purchase.

### 2.1. Pump and Mattress System


- (1) CPR
- (2) Pump unit
- (3) **Mattress system**
- (4) Quick connector
- (5) Turning-cell ventilation connectors
- (6) Bolster control knob

### 2.2. Pump Unit


#### Front

- (1) Quick connector port
- (2) Front panel


#### Rear

- (1) Power switch
- (2) Air filter
- (3) Fuses
- (4) Mounting brackets

## 2.3. Front Panel

### (1) Power Button

The **Power** button is to turn on and turn off pump.

### (2) Skin Check Button

- Allows caregivers to check a patient quickly and thoroughly for pressure injuries. Press it once to start the routine: mattress tilts left to a selected angle → stays for 1 minute → returns to supine position → tilts right to a selected angle → stays for 1 minute → returns to supine position, routine complete.
- If a patient is already in a tilted position when Skin Check starts, the system will start the routine from that tilted position.
- Press the button again to stop the routine. Once stopped, the routine is canceled, and the patient will return to supine position.
- LED indications
  - I. Flashing green: turning cells inflating/deflating.
  - II. Solid green: mattress holding in tilted position.
  - III. Off: routine complete or manually stopped.

### (3) Panel Lock Button

User interface will be locked automatically 3 minutes after last operation or by pressing the Panel Lock button. To unlock the interface, press and hold the Panel Lock button for 1 second.

**NOTE: Always unlock user interface before turning off the pump.**

#### (4) Max Firm Button

The surface will go to max firm automatically when the power is turned on.

**NOTE:** The LED light of Max Firm mode will be on constantly until the mattress is reached to the maximum pressure setting.

**NOTE:** This button is to make sure the mattress can reach maximal operating pressure. This mode will last for 20 minutes and then automatically goes into ALT mode with previous settings. This function can be used during patient transfers or nursing procedures for better support. To disable this function, press the Max Firm button again.

#### (5) Mode Button

##### A. Alternating (ALT)

The air cells inside the mattress deflate and inflate periodically to avoid prolonged pressure in certain areas of patient.

##### B. Continuous Low Pressure (CLP, Static)

Its LED indicator is steady-on in constantly non-alternating mode. The system will redistribute body mass over a greater surface area at a constant low pressure.

##### C. Seat Inflation

Press the **Mode** button to enter the desired mode in ALT or CLP mode for additional supports to the patient during upright position without bottoming out.

#### (6) Turning Mode Button

There are three turning modes:  **LEFT TURN**,  **RIGHT TURN** and  **LEFT-SUPINE-RIGHT TURN**. The turning side is determined from the point of view of a caregiver standing at the pump facing the patient. The top air cells are non-alternating during the turning process and will change to alternating mode when it reaches the position of a desired angle. The same also applies when it comes back to supine position.

Different dwell times can be used in combination with these turning modes. Find more details in **(9) Time Button**.

#### (7) Turning Angle Button (Up to 30–40 degree turning under Full-turn Mode, subject 75 kg with 45.7 cm shoulder width)

-  Light-turn,  Half-turn,  Full-turn

#### (8) Pressure Adjust Buttons

- Press the buttons to get the correct setting of the pressure level based on patient's weight.

 **NOTE: When the power is on (i.e., the surface is inflated), it will automatically go to "Max Firm" mode to hasten inflation. Once the system is ready, the system will enter ALT mode.**

- The mattress pressure can be manually adjusted for patients' comfort using  for a firmer setting and  for a softer setting.

 **NOTE: The caregivers should check if the pressure is suitable for the patient by sliding one hand beneath the air cells at the level of the patient's buttocks. Always leave at least 1 inch space between patient and the static cell to prevent bottoming out.**

- Weight ranges:

|    |                                                                                     |              |              |
|----|-------------------------------------------------------------------------------------|--------------|--------------|
| 10 |    | 161 – 180 kg | 355 – 397 lb |
| 9  |    | 141 – 160 kg | 311 – 353 lb |
| 8  |    | 121 – 140 kg | 267 – 309 lb |
| 7  |    | 106 – 120 kg | 234 – 264 lb |
| 6  |    | 91 – 105 kg  | 201 – 231 lb |
| 5  |    | 81 – 90 kg   | 179 – 198 lb |
| 4  |    | 71 – 80 kg   | 157 – 180 lb |
| 3  |    | 61 – 70 kg   | 134 – 154 lb |
| 2  |    | 51 – 60 kg   | 112 – 132 lb |
| 1  |  | 40 – 50 kg   | 88 – 110 lb  |

**(9) Time Button** 

- While in **Alternating Mode**, press the button to select among 10, 15, 20 and 25 minutes **cycle time**.
- While in **Turning Mode**, press the button to select among 0, 1, 2, 3 and 4 hours (00, 1h, 2h, 3h, 4h) **dwell time**. A dwell time is the approximate time of a patient being tilted in a position or in supine position before being tilted again.

Example 1: When selecting 1h dwell time in Left Turn mode, the cycle will be: mattress tilts left to a selected angle → holds position then returns to supine position: 1 hour → stays in supine position then tilts left to a selected angle: 1 hour.

Example 2: When selecting 2h dwell time in Left-Supine-Right Turn mode, the cycle will be: mattress tilts left to a selected angle → holds position then returns to supine position: 2 hours → stays in supine position then tilts right to a selected

angle: 2 hours → holds position then returns to supine position: 2 hours → stays in supine position then tilts left to a selected angle: 2 hours.

**NOTE: Cycle time and dwell time duration can be selected based upon patient comfort and desired outcome.**

**NOTE: Whenever a patient is in a tilted or supine position in a cycle, ALT mode (20 mins) is turned on.**

- While in **Turning Mode**, press the button to select “00” will turn on single-sided manual turning mode, which means the caregivers can manually control Left/Right turn mode and return patients back to supine position.

**NOTE: For safety reasons, the system will automatically return the patient to supine position after being tilted for 4 hours.**

- To cancel an ongoing cycle or stop the manual turning mode, press either the Alternating Mode button or the Max Firm button. Doing this will return the patient back to supine position.

#### (10) Alarm Mute Button

When there is an ongoing alarmed situation, the number display will show corresponding error code with buzzer sounded every few seconds. The buzzer can be temporarily suspended by pressing the **Alarm Mute Button**. Alarm buzzer will sound repeatedly every 3 minutes until the alarmed issue is solved. For low-pressure alarm while turning, the repeated alarm is every 8 minutes.

#### A. PFA (Power Failure Alarm)

The LED light and buzzer will activate when the power failure happens. To disable both buzzer and LED, press the mute button.

**NOTE: The pump can seal the air inside the mattress, which temporarily keeps the patients from bottoming out during power failure. Patients should be moved if the power failure issue cannot be solved soon.**

#### B. Pressure Anomaly Indicator

When the Pressure Anomaly LED lights up, the pressure inside the air mattress is either below or above the normal range. Please refer to **10. Troubleshooting**.

**NOTE: The pump can seal the air inside the mattress, which temporarily keeps the patients from bottoming out during low pressure condition. Patient should be moved if the power failure issue cannot be solved soon.**


#### C. Service Alarm

The Service Alarm indicator lights up when an issue occurs related to the system. Please contact your local representative or place of purchase for technical support.

## 3. Installation

### 3.1. Pump & Mattress Installation

- ⚠ **NOTE: Unpack the box. Please contact your support personnel or package vendor if there are any damage occurred during shipment.**
- ⚠ **NOTE: USB function is provided for use by distributors or authorized technicians.**

- (1) Place the mattress on top of bed frame with the footprint logo pointing towards the foot end. Secure the mattress firmly by fixing the straps to bed frame to keep the mattress from shifting. The securing straps are located at the base of the mattress.

- ⚠ **NOTE: Make sure that ALL the Heel Relief Connectors and emergency CPR knob are closed and locked in position.**

- ⚠ **NOTE: Make sure the bed side rails are raised.**

- (2) Hang the pump onto bed rail (foot-end) and adjust hangers to best upright position of the pump or place the pump on flat surface.
- (3) Connect air hose connectors of air mattress to pump unit. When a "click" sound is felt or heard, the connection is completed and secured.

- ⚠ **NOTE: Check and ensure the air hoses are not kinked or tucked under mattress.**

- (4) Plug the power cord into electrical outlet.

- ⚠ **NOTE: Make sure the pump unit is suitable for the local power voltage.**

- ⚠ **NOTE: The plug is also served to disconnect the device.**

- (5) Turn the power switch ON.

**CAUTION: The pump can only be applied to the mattress recommended by the manufacturer. Do not use it for any other purpose. (Applied part: air mattress) Do not position the equipment so that it is difficult to operate the disconnecting device.**

- (6) The cable management tape is provided along the side of the mattress to manage power cable. Please make sure the wire has no signs of damage. After installation, extra length of the power cord, if any, should be neatly managed to prevent stumble accident and should be clear of moving bed mechanisms or other possible entrapment area.

**NOTE: The EQUIPMENT should be firmly placed at position where caregivers can access easily.**

## 4. Operation

**NOTE: Always read the operating instruction before use.**

### 4.1. General Operation

- (1) Switch the power on, the mattress starts the initial inflation. The Max Firm light flashes until the mattress is fully inflated, which takes approximately 30 min.
- (2) When the inflation is completed, the system will enter ALT mode with prior settings.
- (3) Use the Pressure Adjust buttons to get the correct setting of therapeutic pressure level based on patient weight.

**NOTE: Pressure levels with the weight guide on the panel is only an indication and clinical judgement should always be taken into consideration.**

- (4) For lateral turning, choose the Left-Supine-Right Turn button to get standard turning which helps the patient to turn to both sides and back to supine position, or either one direction in Left turn or Right turn mode.


**NOTE: Do the hand check to check out if the pressure reduction is adequate and the patient has not bottomed out.**

**NOTE: For patient comfort, always use the head pillow when patient is under turning mode.**


## 4.2. Emergency CPR Operation

⚠ **NOTE:** When there is a need to perform CPR on the patient, quickly turn the CPR knob to release air from mattress. CPR knob is located at the right-head side of the mattress. The quick connector on the pump unit can be disconnected for even faster deflation.

## 4.3. Emergency Patient Reposition

The turning-cell ventilation connectors are used to manually return the patient to a level supine position in the event of a power failure. They are hidden inside the air hose cover. Pull down the cover to find the two connectors labeled L and R: "L" for the left side and "R" for the right side. Simply pull apart the raised side connector to deflate and then re-connect to stop.

## 4.4. Heel Relief Function

The last five cells from the foot-end have heel-relief connectors. Carers can easily disconnect the connector, and heels can be suspended in air to eliminate pressure exerted on the heels.

For permanent off-loading/pressure relief:


- Select **no more than two** cells directly under the area that desired to relieve.

⚠ **NOTE:** These one to two cells can be left permanently deflated.

⚠ **NOTE:** Deflating more than two adjacent cells may affect the support of the patient during the normal alternating cycle.

## 4.5. Transportation


**NOTE:** When there is a need of patient transportation, pump with the power cord can be disconnected from the mattress during transport. For optimum comfort, change the system to CLP Mode and wait for 5 minutes before disconnecting is recommended. Cap on the quick connector as shown on the left to keep the air inside the mattress. The mattress remains inflated for at least 8 hours depending on the setting.

## 4.6. Patient Transferring

Located at the side near the foot end, the **Bolster Control Knob** is used to adjust the height of the side bolster when the patient needs to be moved.

- To deflate the side bolster: turn the knob counterclockwise to OPEN the ventilation.
- To inflate the side bolster: turn the knob clockwise to CLOSE the ventilation.

**NOTE:** Switching to Max Firm mode before transferring patients is recommended.

**NOTE:** Make sure the Bolster Control Knob is properly closed after each patient transfer.

## 5. Cleaning

It is important to follow the cleaning procedures to avoid cross contamination. Be sure to clean the surface in a dry and dust free environment. Wipe down the pump unit with a damp cloth pre-soaked with a mild detergent. Avoid contact with dust and proximity to dusty areas. Make sure that any cleaning agents you use will not harm or corrode the plastic casing on the pump unit. If your doctor or medical facilities have other special cleaning instruction, please follow the professional instruction.

**CAUTION: Do not immerse or soak pump unit.**

**WARNING: Do not remove the casing of the pump to avoid electrical shock. All disassembly or repair should be done by authorized technicians.**

**CAUTION: The pump does not need oil lubrication; please do not disassemble the system.**

|                                                               |                                                                                                                                                                                                                                                                                                                                                                                                                            |
|---------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Cover Material: Stretch<br/>Silver<sup>+</sup> Stretch</b> |      |
|---------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

Wipe down the mattress with 70 % alcohol solution or chlorine solution concentration must be less than or equal to 1,000 ppm. After disinfecting the mattress, use water and a soft cloth to wipe off all traces of disinfectant and make sure the mattress is dry before storage. Also, the mattress top cover can be completely removed for laundry with water temperature indicated on the laundry label; however, it is recommended that user should check with local policy or regulation to determine the time/temperature ratio required achieving thermal disinfection. The cover may also be cleaned using sodium hypochlorite diluted in water. After cleaning, please avoid dust and proximity to dusty areas and all parts should be air dried thoroughly before use.

**CAUTION: Do not use phenolic based products for cleaning.**

**CAUTION: After cleaning, dry the mattress without direct exposure of sunlight.**

**CAUTION: Do not use ozone or UV disinfection as these methods can damage the surface coating.**

## 6. Storage

- (1) To store the mattress, lay the mattress out flat and upsides down.
- (2) Roll from the head end towards the foot end with CPR valve open.
- (3) Foot-end strap can then be stretched around the rolled mattress to prevent unrolling.

**NOTE: Do not kink, crease or stack the mattresses and do not store the system in direct sunlight, high temperature or moisture area.**


## 7. Maintenance

### 7.1. General

- (1) Check main power cord and plug if there is abrasions or excessive wear.
- (2) Check mattress cover for signs of wear or damage. Ensure mattress cover and tubes are stubbed together correctly.
- (3) Check airflow from the opening of the tubes in the quick connector. The airflow should alternately blow every half of the cycle time between two tubes at the lower right corner and one of the other two tubes at the lower left corner will blow to give desired turning angle when it is under turning mode.
- (4) Check the air hoses if there is kink or break. For replacement, please contact local representative or place of purchase.

### 7.2. Fuse Replacement

- (1) Disconnect the plug from mains power when a blown fuse is suspected.
- (2) Remove the cover of the fuse holder by means of a flat screwdriver.
- (3) Insert a new fuse of the correct rating in and replace the cover of the fuse holder back. The fuse should be rated as T1AL/250V type and VDE approved.

### 7.3. Air Filter Replacement

- (1) Open the air filter cap located at the back of pump.
- (2) Replace with a clean filter and put the cover lid back. The filter is reusable and can be washed gently with a mild detergent and water. Dry the filter before use.
- (3) Suggested to check and replace air filter regularly every six months.

## 8. Servicing

The products are intended to offer safe and reliable operation when use or installed according to the instructions provided by Wellell. Wellell recommends that the system be inspected and serviced by authorized technicians if there are any signs of wear or concerns with device function and indication on products. Otherwise, service and inspection of the devices generally should not be required.

## 9. Limited Warranty

Wellell Inc. ("Wellell") provides a limited manufacturer warranty on new, original Wellell products and any replacement part fitted by Wellell in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from the date of purchase. This warranty does not cover damage caused by accident, misuse, abuse, alteration and other defects not related to material or workmanship. To exercise your rights under this warranty, please contact your local authorized Wellell dealer.

## 10. Troubleshooting

### **Q1 Power is not ON**

- Check if the plug is connected to the mains power.
- Check for a blown fuse.

### **Q2 Pressure Anomaly Alarm is on**

- Check if the Quick Connector is tightly secured.
- Check if all tubing connections along mattress are secured.
- Check if the air hoses are kinked or broken.
- Check if the CPR knob is tightly closed.
- Check if the bolster control knob is tightly closed.
- Check if the turning-cell ventilation connectors are properly connected and have no leakage.

### **Q3 Power Failure Alarm is on**

- Check the switch on/off position.
- Check if there is a power outage.
- Check if the power cord is connected properly.

### **Q4 Patient is bottoming out**

- Pressure setting might be inadequate for the patient, increase comfort range 1 level higher and wait for a few minutes for best comfort.

### **Q5 Mattress form is loose**

- Check if all the snap buttons or straps of mattress are all securely fastened.
- Check if the mattress is fixed to the bed frame by straps.

### **Q6 No air released from some air outlets of the air tube connector**

- This is in alternating mode. The air is released through outlets by turns at each cycle.

### **Q7 Possible system failure**

- Stop using the system and move patient onto another support surface.
- Call customer service for further assistance.

If the problems you have are not solved by the above information, please contact local representative or place of purchase.

## 11. Error Code Table

| <b>Error Code No.</b> | <b>Error Type</b> | <b>Description</b>                  | <b>Note</b>                                                                             |
|-----------------------|-------------------|-------------------------------------|-----------------------------------------------------------------------------------------|
| Er00                  | Low Pressure      | Initial inflation overtime          | Check if all tubing connections along mattress are secured.                             |
| Er01–03               | Low Pressure      | Alternating cells                   | Check if the air hoses are kinked or broken.                                            |
| Er04 & 05             | Low Pressure      | Turning cells                       | Check if the turning-cell ventilation connectors are properly connected and no leakage. |
| Er06                  | Low Pressure      | Head, bolster cells                 | Check if the bolster control knob is tightly closed.                                    |
| Er07                  | Service           | Compressor malfunction              | Service required                                                                        |
| Er08–12               | Service           | System malfunction                  | Service required                                                                        |
| Er13                  | Service           | Vent valve malfunction              | Service required                                                                        |
| Er16                  | NA                | Pressure control system malfunction | Service required                                                                        |
| Er17                  | NA                | System access control malfunction   | Service required                                                                        |
| Er18                  | NA                | Mattress disconnected               | Check if the Quick Connector is tightly secured.                                        |
| Er19                  | NA                | Wrong mattress detected             | Make sure the correct mattress is connected.                                            |

## 12. Technical Specifications

| Item                                                 |                      | Specification                                                                                                                                                              |
|------------------------------------------------------|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Pump</b>                                          |                      |                                                                                                                                                                            |
| Power Supply (Note: See rating label on the product) |                      | AC 230V 50 Hz, 0.25A (for 230V system)                                                                                                                                     |
| Fuse Rating                                          |                      | T1AL, 250V                                                                                                                                                                 |
| Cycle Time (Alternating)                             |                      | 10, 15, 20, 25 min.                                                                                                                                                        |
| Dwell Time (Turning)                                 |                      | 1, 2, 3, 4 hours                                                                                                                                                           |
| Dimensions (L x W x H)                               |                      | 34.1 x 16.5 x 26.0 cm                                                                                                                                                      |
| Weight                                               |                      | 4.5 kg / 9.9 lb                                                                                                                                                            |
| Environment                                          | Atmospheric Pressure | 700 hPa to 1060 hPa                                                                                                                                                        |
|                                                      | Temperature          | Operation: 10°C to 40°C (50°F to 104°F)<br>Storage: -15°C to 50°C (5°F to 122°F)<br>Shipping: -15°C to 70°C (5°F to 158°F)                                                 |
|                                                      | Humidity             | Operation: 10% to 90% non-condensing<br>Storage: 10% to 90% non-condensing<br>Shipping: 10% to 90% non-condensing                                                          |
| Classification                                       |                      | Class II , Type BF, IP21<br><b>Applied Part: Air Mattress</b><br>⚠ <b>Not suitable for use in the presence of a flammable anesthetic mixture (No AP or APG protection)</b> |
| Service Life (Design Life)                           |                      | System: 7 years                                                                                                                                                            |
| <b>Mattress</b>                                      |                      |                                                                                                                                                                            |
| Model                                                |                      | 5" cells & 10" side bolster                                                                                                                                                |
| Dimensions (L x W x H)                               |                      | 200 x 85/90 x 13 cm / 78.7 x 33.5/35.4 x 5 in.                                                                                                                             |
| Weight                                               |                      | 12.5 kg / 27.5 lb                                                                                                                                                          |
| Max. Support Weight                                  |                      | 180 kg / 397 lb                                                                                                                                                            |

### NOTE:

- Consult the distributor or EU representative for further technical documents.
- The specification is also suitable for other areas operating with same power supply.
- The manufacturer reserves the right to modify the specifications without notice.

# Appendix A: EMC Information

## Guidance and Manufacturer’s Declaration- Electromagnetic Emissions:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

| Emissions Test                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            | Compliance | Electromagnetic Environment-Guidance                                                                                                                                           |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| RF emissions<br>CISPR 11                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | Group1     | The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment |
| RF emissions<br>CISPR 11                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | Class B    | The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network            |
| Harmonic emissions<br>IEC61000-3-2                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | Class A    |                                                                                                                                                                                |
| Voltage fluctuations / Flicker emissions<br>IEC61000-3-3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | Complies   |                                                                                                                                                                                |
| <div><div></div><div><b>Warning:</b></div></div> <div><div>1.</div><div>The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used.</div></div> <div><div>2.</div><div>Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.</div></div> <div><div>3.</div><div>Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Pump, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.</div></div> |            |                                                                                                                                                                                |

## Guidance and Manufacturer’s Declaration- Electromagnetic Immunity:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

| Basic EMC standard      | Immunity Test Levels                         |                             | Compliance Levels | Electromagnetic Environment-Guidance |
|-------------------------|----------------------------------------------|-----------------------------|-------------------|--------------------------------------|
|                         | Professional healthcare facility environment | HOME HEALTHCARE ENVIRONMENT |                   |                                      |
| Electrostatic Discharge | ±8 kV contact<br>± 2, 4, 8, 15 kV air        | T                           | ±8 kV contact     | Floors should be wood, concrete or   |

|                                                                                                       |                                                                                                                                                                                               |                                                                                                                      |                                                                                       |                                                                                                                                                                                                                                                                               |
|-------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| (ESD)<br>IEC61000-4-2                                                                                 |                                                                                                                                                                                               |                                                                                                                      | ± 2, 4, 8, 15 kV air                                                                  | ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.                                                                                                                                                                   |
| Electrical fast transient/ burst<br>IEC61000-4-4                                                      | ±2 kV, 100 kHz for power supply line<br>±1 kV, 100 kHz for input/output line                                                                                                                  |                                                                                                                      | ±2 kV, 100 kHz for power supply line                                                  | Mains power quality should be that of atypical commercial or hospital environment                                                                                                                                                                                             |
| Surge<br>IEC61000-4-5                                                                                 | ± 0.5, 1 kV (line to line)<br>± 0.5, 1, 2 kV (line to earth)                                                                                                                                  |                                                                                                                      | ± 0.5, 1 kV (line to line)                                                            | Mains power quality should be that of atypical commercial or hospital environment.                                                                                                                                                                                            |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br>IEC61000-4-11 | Voltage Dips:<br>i) 100% reduction for 0.5 period,<br>ii) 100% reduction for 1 period,<br>iii) 30% reduction for 25/30 period,<br>Voltage Interruptions:<br>100% reduction for 250/300 period |                                                                                                                      | 230 Vac                                                                               | Mains power quality should be that of a typical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60Hz) magnetic field<br>IEC61000-4-8                                              | 30 A/m                                                                                                                                                                                        | 30 A/m                                                                                                               | 30 A/m                                                                                | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.                                                                                                                                     |
| Conducted RF<br>IEC 61000-4-6                                                                         | 3 Vrms<br>0,15 MHz – 80 MHz<br>6 Vrms in ISM bands between 0,15 MHz and 80 MHz<br>80 % AM at 1 kHz                                                                                            | 3 Vrms<br>0,15 MHz – 80 MHz<br>6 Vrms in ISM and amateur radio bands between 0,15 MHz and 80 MHz<br>80 % AM at 1 kHz | 3 Vrms<br>0,15 MHz – 80 MHz<br>6 Vrms in ISM and amateur radio bands between 0,15 MHz | Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there commended separation distance calculated from the equation applicable                                                                       |

|                                                                                                                                                                                                                                                                                                                                |                                                                                                                                          |                                                                                                                                           |                                                                                                                                           |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                                                                                                                                                                                                                                                                                                |                                                                                                                                          | kHz                                                                                                                                       | and 80 MHz<br>80 % AM<br>at 1 kHz                                                                                                         | to the frequency of the transmitter.<br><br>Recommended separation distance<br>$d = \sqrt{P}$ 150kHz to 80MHz<br>$d = 0.6\sqrt{P}$ 80MHz to 800MHz<br>$d = 1.2\sqrt{P}$ 800 MHz to 2.7GHz<br><br>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). <sup>b</sup><br><br>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency ranged.<br><br>Interference may occur in the vicinity of equipment marked with the following symbol:<br> |
| Radiated RF EM Fields<br>IEC61000-4-3                                                                                                                                                                                                                                                                                          | 3 V/m<br>80 MHz - 2.7 GHz<br>80 % AM at 1 kHz<br><br>380 - 5800 MHz, 9-28 V/m, Pulse modulation or FM $\pm 5$ kHz deviation, (1kHz sine) | 10 V/m<br>80 MHz - 2.7 GHz<br>80 % AM at 1 kHz<br><br>380 - 5800 MHz, 9-28 V/m, Pulse modulation or FM $\pm 5$ kHz deviation, (1kHz sine) | 10 V/m<br>80 MHz - 2.7 GHz<br>80 % AM at 1 kHz<br><br>380 - 5800 MHz, 9-28 V/m, Pulse modulation or FM $\pm 5$ kHz deviation, (1kHz sine) |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| NOTE 1: $U_T$ is the a.c. mains voltage prior to the application of the test level<br>NOTE 2: At 80 MHz and 800 MHz, the higher frequency range applies.<br>NOTE 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people |                                                                                                                                          |                                                                                                                                           |                                                                                                                                           |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |

- a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

**Recommended separation distances between portable and mobile RF communications equipment and this device:**

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment

| Rated maximum<br>output power<br>of transmitter<br>W | Separation distance according to frequency of transmitter m |                                        |                                         |
|------------------------------------------------------|-------------------------------------------------------------|----------------------------------------|-----------------------------------------|
|                                                      | 150 kHz to 80 MHz<br>$d = \sqrt{P}$                         | 80 MHz to 800 MHz<br>$d = 0.6\sqrt{P}$ | 800 MHz to 2.7 GHz<br>$d = 1.2\sqrt{P}$ |
| 0.01                                                 | 0.1                                                         | 0.06                                   | 0.12                                    |
| 0.1                                                  | 0.31                                                        | 0.19                                   | 0.38                                    |
| 1                                                    | 1                                                           | 0.6                                    | 1.2                                     |
| 10                                                   | 3.1                                                         | 1.9                                    | 3.8                                     |
| 100                                                  | 10                                                          | 6                                      | 12                                      |

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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All rights reserved.