EN	ANAESTHESIA BREATHING SYSTEMS
BG	АНЕСТЕТИЧНИ ДИХАТЕЛНИ СИСТЕМИ
CS	ANESTEZIOLOGICKÉ DÝCHACÍ SYSTÉMY
DA	ANÆSTESI-RESPIRATIONSSYSTEMER
DE	FLEXICARE ANÄSTHESIE-ATEMSYSTEME
EL	ΑΝΑΙΣΘΗΤΙΚΑ ΑΝΑΠΝΕΥΣΤΙΚΑ ΣΥΣΤΗΜΑΤΑ
ES	SISTEMAS RESPIRATORIOS PARA USO ANESTÉSICO
ET	ANESTEESIAHINGAMISSÜSTEEMID
FI	ANESTESIAHENGITYSJÄRJESTELMÄT
FR	SYSTÈMES RESPIRATOIRES D'ANESTHÉSIE
HR	SUSTAVI ZA DISANJE TIJEKOM ANESTEZIJE
HU	ANESZTÉZIÁS LÉGZŐRENDSZEREK
ID	SISTEM PERNAPASAN ANESTESI
IT	CIRCUITI RESPIRATORI PER ANESTESIA
JP	麻酔用ブリージングシステム
KO	마취 호흡 시스템
LT	ANESTEZIJOS KVĖPAVIMO SISTEMOS
LV	ANESTĒZIJAS ELPINĀŠANAS SISTĒMAS
NL	BEADEMINGSSYSTEMEN VOOR ANESTHESIE
NO	ANESTETISKE PUSTESYSTEMER
PL	UKŁADY ODDECHOWE DO ZNIECZULANIA
BT-BR	SISTEMAS DE RESPIRAÇÃO EM ANESTESIA
PT	SISTEMAS DE RESPIRAÇÃO SOB ANESTESIA
RO	SISTEME RESPIRATORII ANESTEZICE
RU	АНЕСТЕЗИОЛОГИЧЕСКИЕ ДЫХАТЕЛЬНЫЕ КОНТУРЫ
SR	ANESTEZIJSKI DISAJNI SISTEMI
SK	ANESTÉZIOLOGICKÉ DÝCHACIE SYSTÉMY
SV	ANDNINGSSYSTEM FÖR ANESTESI
TR	ANESTEZİ NEFES ALMA SİSTEMLERİ
UA	СИСТЕМИ ІНГАЛЯЦІЙНОЇ АНЕСТЕЗІЇ
VI	CÁC HỆ THỐNG KHÍ THỞ GÂY MÊ
ZH	麻醉呼吸回路

SYMBOLS GLOSSARY

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These instructions are supplied with and apply to all Flexicare® Adult, Paediatric and Infant/Neonatal Anaesthesia Breathing Systems, including Argentic™ & Nitrous Oxide/Oxygen Systems. Please see individual package labels for full product details.

Intended Use

Flexicare® Anaesthesia Breathing Systems are used for Adult, Paediatric and Infant/Neonatal patients undergoing general anaesthesia, and are intended for both ventilated and spontaneously breathing patients within hospital operating theatres, anaesthetic rooms and recovery (PACU), during emergency and standard anaesthesia practice. The user should be familiar and trained in the selection of the appropriate Anaesthesia Breathing System, the use and the care of the product, and should be aware of any clinical contraindications.

Precautions

Do not use if packaging is open, or packaging or contents are damaged, or after the expiry date.

- Before use, visually inspect all connectors to confirm they are securely in place. Immediately before each procedure, carry out a functional test to ensure there is no occlusion or leakage in the system.

Connecting the Breathing System

These instructions should be followed before each use in association with hospital policy and local guidelines:

1. Do not open packaging until immediately prior to use.
2. Before each use, visually inspect all components, ensure all connections are secure and carry out functional tests on the assembled breathing system to ensure there are no leakages or occlusions. Co-axial Bain and Ventiflow™ systems are supplied with a red lumen patency tester to isolate the inner tube for leak testing.
3. Remove the red patient end safety cap before connection to the patient airway device, and replace whenever the Breathing System is not in use to prevent contamination, damage or entry of foreign objects.
4. Use a push and twist action to connect components, and a twist and pull action to disconnect.

Warnings/Cautions

- Adult Anaesthesia Breathing Systems are single use. However, unless otherwise dictated by local hospital policy and local guidelines, Adult Anaesthesia Breathing Systems can be used on multiple patients provided a new Flexicare® or other appropriately validated* bacterial/viral breathing filter is used at the patient end of the system for each individual patient. In this case, Flexicare® Anaesthesia Breathing Systems are validated for a period of use up to 7 days or 10 operating lists, whichever is shorter.
- Paediatric and Infant/Neonatal Anaesthesia Breathing Systems are single use. However, hospital policy and local guidelines may permit repeat use under certain conditions.
- The removal of any component, other than a patient end filter, from any assembled breathing system renders it single use, and further use of that system is not permitted.
- The use of a suitable Anaesthetic Gas Scavenging System is recommended to limit user exposure to expired and vented gases.
- Do not use a Heat Moisture Exchanger (HME), Heat Moisture Exchanger with Filter (HMEF) or Pleated Paper Filter between a patient and a humidification device.
- Contamination: Regardless of the designated use, all the components of the Breathing System should be discarded immediately after use where there is actual or suspected infection that might be transmissible via the Breathing System or where visible contamination is present (e.g. blood or other fluids).
- This product may not be decontaminated in any way. This includes rinsing, washing or decontamination using gas, heat, steam, or boiling water.
- Do not reuse this device. Reuse may result in hazards such as, but not limited to; cross contamination, damage to device and/or suboptimal performance.
- After use, dispose of in accordance with local policies and / or as advised by your local healthcare provider.

Any serious incident that occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State of the user and/or patient.

Protective Filters

- Anaesthesia machines may be further protected from contamination by a validated* bacterial/viral filter at the machine end of the system, fitted to either the inspiratory, expiratory or both machine ports.
 - Appropriately validated* bacterial/viral filters & HMEFs may also be used at the patient end of the system.
 - All filters should be used in accordance with the manufacturer's instructions.
- *A validated bacterial/viral filter meets BS EN ISO 23328 Breathing system filters for anaesthetic and respiratory use. Parts 1 & 2.

Adjustable Pressure Limiting (APL) Valves

- Check the APL Valve for correct function prior to each use and adjust as clinically required during the procedure.
- Adult (blue dial): Maximum opening pressure 60cm H₂O (5.9kPa).
- Paediatric (yellow dial): Maximum opening pressure 35cm H₂O (3.4kPa)
- Do not use Flexicare® APL Valves within MRI scanners or similar equipment as the valve contains metallic components. It may be possible to use a Breathing System with an APL valve located a suitable remote distance from the scanning field, eg next to the anaesthetic machine.

Record Keeping & Expiry Date

- Each Anaesthesia Breathing System is supplied with an adhesive label for record keeping and traceability. Complete the date, time, change date, and lot number at the time of opening and attach the label to the Breathing System at the anaesthetic machine end.
- Do not use after the end of the expiry date printed on the packaging or after the change date noted on the adhesive label.

Argentic™ Anti-Microbial Anaesthesia Breathing Systems

- The tubing of Flexicare® Argentic™ Anaesthesia Breathing Systems are impregnated with silver ions to reduce microbial contamination of the tube outer surfaces and reduce the risk of cross-contamination and Hospital Acquired Infections (HAIs).

Storage

Store in a cool, dry place out of direct sunlight.

Technical Data

Please contact Flexicare® for any additional technical data.

Visit www.flexicare.com/downloads for symbols guide and information on the risks associated with the re-use of single use devices.

All information in this leaflet is correct at the time of printing.

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**Flexicare Medical Limited**
**Flexicare Deutschland GmbH**