Hygrovent HF55


UK
CA
0120

CE
1639


STERILE EO


EN - BACTERIAL VIRAL HME FILTER


Hygrovent


Hygrovent S


Hygrovent S
Angled

HF55

Hygrovent Child
Straight

Hygrovent Child
Angled

210mm

EN - BACTERIAL VIRAL HME FILTER

INTENDED USE

This device is intended to provide a medium through which gases may pass to reduce the transmission of bacteria and viruses to and from an adult or pediatric patient in anesthesia and critical care and/or provide humidification to those gases, in a pre-hospital, hospital and homecare environments. For use with anaesthesia, artificial respiration and other types of assisted ventilation. The device is intended on being used on a single patient up to 24hrs and is single use. It is to be used by trained users only.

DESCRIPTION

Electrostatic filters use a bipolar, permanently charged filter media. Electrostatic filters are recommended for applications where high efficiency and low pressure drop are required. Hygrovent filters incorporate an HME (Heat Moisture Exchange) medium which uses the patients passive humidification to humidify the patients inspired gases as dry gases can cause adverse effects to the intubated patient.

INDICATION

No specific indication.

CONTRAINDICATION

- Do not use HMEF with patients producing fulminating frothy secretions within their airways, or patients with conditions such as haemoptysis or pulmonary edema.
- Do not use in conjunction with conventional humidifiers and active humidification.

CLINICAL BENEFIT

HME Filter is able to protect the patient against the risk of infection by having electrostatically charged media which traps bacteria and viruses. They are also able to protect the patient from damaging the endothelial tissues of the airway due to inadequate humidity, by conserving a portion of the patient's exhaled heat and moisture, and condition inspired gas by warming and humidifying it.

RESIDUAL RISKS

The device should be discarded after a maximum usage of 24 hours. Reuse may result in hazards such as, but not limited to; cross contamination, damage to device and/or suboptimal performance.

PRECAUTIONS

- Do not use if packaging is opened or damaged.
- Before use, visually inspect all connections to confirm they are secure. Immediately before each procedure, carry out a functional test of the breathing system to ensure there is no occlusion or leakage.
- Do not use after the expiry date on the packaging.
- Do not mount the filter with the use of lubricants and other fluids.

WARNINGS AND CAUTIONS

- The HME filter should be replaced immediately if contaminated by secretions or haemoptysis, this can result in significantly increased resistance.
- Do not use between the patient and a humidification device or nebulizing device.
- Consider total system dead space before use.
- Replace the filter immediately if it becomes soiled, ventilation is compromised, or significant pressure change is observed.
- Do not attempt to 'dry out' a breathing system or a breathing system HME filter.
- Do not use with inflammable anesthetic gases.
- Monitor the pressure-drop across the filter if it is being used in combination with medicine nebulizers.
- Leakage around the endotracheal tube can cut down the effective humidification.
- Do not decontaminate / reprocess this device in any way. This includes rinsing, washing or decontamination using gas, heat, steam, or boiling water.
- Do not reuse this device. Reuse may result in hazards such as, but not limited to; cross contamination, damage to device and/or suboptimal performance.
- Do not add moisture to the HME prior to use.
- Do not use filter longer than 24hrs.
- Monitoring port should only be used for respiratory monitoring. Inappropriate connection of monitoring port to other devices may result in potential hazard.
- Ensure accurate HMEF connection on both patient and machine end to reduce humidity build-up at the device.
- Expert clinical judgment must be used in assessing the patient's humidification requirements.

- Maintenance of normal breathing pressure is important. This will prevent any condensation from being forced through the filter.

INSTRUCTIONS FOR USE

1. This device is single use. If used in combination with a reusable Breathing system, a new filter must be used for each patient for a maximum duration of 24hrs. A new Breathing System Filter must be used following this period of use. Follow manufacturer's instructions and local guidelines / hospital policy on the reuse of breathing systems.
2. Securely attach all connections using a 'push and twist' action to ensure a secure leak free connection.
3. Leak test the breathing system with all accessories attached prior to connecting to patient.
4. The HME filter is to be located at the patient / breathing system interface (patient Y piece) and it must be replaced after each patient.
5. If the luer-lock connection is not being used, check that the cap has been tightened properly.

Any serious incident that occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State of the user and/or patient.

DISPOSAL CONSIDERATIONS

After use, dispose of in accordance with local policies and / or as advised by your local healthcare provider.

STORAGE

Store in a cool, dry place out of direct sunlight.

EN - BACTERIAL VIRAL HME FILTER

PERFORMANCE CHARACTERISTICS

Aspects	HF55 HMEF	Hygrovent HMEF	Hygrovent S HMEF	Hygrovent S Angled HMEF	Hygrovent Child Angled HMEF	Hygrovent Child Straight HMEF
Bacterial Filtration Efficiency (%)	99.999	99.9999	99.999		99.99	
Viral Filtration Efficiency (%)	99.99	99.999	99.99		99.9	
Recommended Range of Tidal Volume (ml)	270-1000	475-1000	270-1000		70-500	
Pressure Drop at Gas Flow Rate (cmH ₂ O)	Inspiratory (Before Usage): 1.35@30L/min 3.13@60L/min 5.68@90L/min	Inspiratory (Before Usage): 0.72@30L/min 1.72@60L/min 3.09@90L/min	Inspiratory (Before Usage): 1.24@30L/min 2.85@60L/min 4.90@90L/min	Inspiratory (Before Usage): 1.17@30L/min 2.63@60L/min 4.37@90L/min	Inspiratory (Before Usage): 2.52@15L/min 5.58@30L/min 9.32@45L/min	Inspiratory (Before Usage): 2.46@15L/min 5.51@30L/min 9.15@45L/min
	Expiratory (Before Usage): 1.29@30L/min 2.87@60L/min 4.92@90L/min	Expiratory (Before Usage): 0.67@30L/min 1.44@60L/min 2.60@90L/min	Expiratory (Before Usage): 1.14@30L/min 2.52@60L/min 4.18@90L/min	Expiratory (Before Usage): 1.17@30L/min 2.57@60L/min 4.33@90L/min	Expiratory (Before Usage): 2.42@15L/min 5.18@30L/min 8.59@45L/min	Expiratory (Before Usage): 2.44@15L/min 5.35@30L/min 8.91@45L/min
	Inspiratory (After 24 hrs of Usage): 1.43@30L/min 3.34@60L/min 5.65@90L/min	Inspiratory (After 24 hrs of Usage): 0.66@30L/min 1.97@60L/min 1.31@90L/min	Inspiratory (After 24 hrs of Usage): 1.44@30L/min 2.74@60L/min 4.36@90L/min	Inspiratory (After 24 hrs of Usage): 1.23@30L/min 2.71@60L/min 4.50@90L/min	Inspiratory (After 24 hrs of Usage): 2.46@15L/min 5.36@30L/min 8.99@45L/min	Inspiratory (After 24 hrs of Usage): 2.41@15L/min 5.54@30L/min 9.17@45L/min
	Expiratory (After 24 hrs of Usage): 1.41@30L/min 2.95@60L/min 4.84@90L/min	Expiratory (After 24 hrs of Usage): 0.72@30L/min 1.53@60L/min 2.66@90L/min	Expiratory (After 24 hrs of Usage): 1.53@30L/min 2.45@60L/min 4.83@90L/min	Expiratory (After 24 hrs of Usage): 1.21@30L/min 2.69@60L/min 4.26@90L/min	Expiratory (After 24 hrs of Usage): 2.28@15L/min 5.22@30L/min 8.62@45L/min	Expiratory (After 24 hrs of Usage): 2.35@15L/min 5.30@30L/min 8.71@45L/min
Internal Volume (ml)	54.32	95.8	53.5	63.7	18.52	13.86
Particle Filtration Efficiency (%)	Before Usage: 99.350 After 24 hrs of Usage: 80.080	Before Usage: 99.800 After 24 hrs of Usage: 99.246	Before Usage: 99.450 After 24 hrs of Usage: 99.549		Before Usage: 80.940 After 24 hrs of Usage: 70.580	
Gas Leakage (ml/min)	Before Usage: 1.0@70cmH ₂ O After 24 Hours of Usage: 0.0@70cmH ₂ O	Before Usage: 0.5@70cmH ₂ O After 24 Hours of Usage: 0.0@70cmH ₂ O		Before Usage: 1.0@70cmH ₂ O After 24 Hours of Usage: 0.0@70cmH ₂ O	Before Usage: 0.0@70cmH ₂ O After 24 Hours of Usage: 0.0@70cmH ₂ O	Before Usage: 0.5@70cmH ₂ O After 24 Hours of Usage: 0.0@70cmH ₂ O
Moisture Output (mg/L)	Minimum Tidal Volume: 32.20 Operating Range (At 500ml): 28.95 Maximum Tidal Volume: 30.40	Minimum Tidal Volume: 33.17 Operating Range (At 500ml): 31.96 Maximum Tidal Volume: 32.70	Minimum Tidal Volume: 32.49 Operating Range (At 500ml): 32.15 Maximum Tidal Volume: 32.52		Minimum Tidal Volume: 32.90 Operating Range (At 250ml): 29.13 Maximum Tidal Volume: 29.81	
Moisture Loss (mg/L)	Minimum Tidal Volume: 6.32 Operating Range (At 500ml): 10.13 Maximum Tidal Volume: 8.42	Minimum Tidal Volume: 5.19 Operating Range (At 500ml): 6.60 Maximum Tidal Volume: 5.21	Minimum Tidal Volume: 5.97 Operating Range (At 500ml): 6.38 Maximum Tidal Volume: 6.39		Minimum Tidal Volume: 5.50 Operating Range (At 250ml): 10.87 Maximum Tidal Volume: 5.12	
Compliance (ml/cmH ₂ O)	0.057	0.085	0.044	0.043	0.014	0.015

Note: Results stated above are based on average value

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Flexicare Medical Switzerland GmbH
c/o BK-Services AG Suurstoffi 22, 6343
Rotkreuz Switzerland

ThermoShield

UK
CA
0120

CE
1639


STERILE EO


EN - BACTERIAL VIRAL HME FILTER


ThermoShield


ThermoShield
Angled

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ThermoShield
(Tethered Cap Port)

ThermoShield Angled
(Tethered Cap Port)


ThermoShield Mini
(Tethered Cap Port)

ThermoShield
(LuerSafe Port)

ThermoShield Mini
(LuerSafe Port)

ThermoShield Angled
(LuerSafe Port)

210mm

EN - BACTERIAL VIRAL HME FILTER

INTENDED USE

This device is intended to provide a medium through which gases may pass to reduce the transmission of bacteria and viruses to and from an adult or paediatric patient in anaesthesia and critical care and/or provide humidification to those gases, in a pre-hospital, hospital and homecare environments. For use with anaesthesia, artificial respiration, and other types of assisted ventilation. The device is intended on being used on a single patient up to 24hrs and is single use. It is to be used by trained users only.

DESCRIPTION

Electrostatic filters use a bipolar, permanently charged filter media. Electrostatic filters are recommended for applications where high efficiency and low-pressure drop are required. ThermoShield filters incorporate an HME (Heat Moisture Exchange) medium which uses the patient's passive humidification to humidify the patients inspired gases as dry gases can cause adverse effects to the intubated patient.

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WARNINGS AND CAUTIONS

- The HME filter should be replaced immediately if contaminated by secretions or haemoptysis, this can result in significantly increased resistance.
- Replace the HME filter immediately if it becomes soiled, ventilation is compromised, or significant pressure change is observed.
- Do not use between the patient and a humidification device or

nebulizing device.

- Consider total system dead space before use.
- Do not attempt to 'dry out' a breathing system or a breathing system HME filter.
- Do not use with inflammable anaesthetic gases.
- Do not reuse this device. Reuse may result in hazards such as, but not limited to; cross contamination, damage to device and/or suboptimal performance.
- Do not use filter longer than 24hrs.
- Do not add moisture to the HME prior to use.
- Do not decontaminate / reprocess this device in any way. This includes rinsing, washing or decontamination using gas, heat, steam, or boiling water.
- Monitor the pressure-drop across the HME filter if it is being used in combination with medicine nebulizers.
- Leakage around the endotracheal tube can cut down the effective humidification.
- Monitoring port should only be used for respiratory monitoring. Inappropriate connection of the monitoring port to other devices may result in potential hazard.
- Ensure accurate HMEF connection on both patient and machine end to reduce humidity build-up.
- Expert clinical judgment must be used in assessing the patient's humidification requirements.
- Maintenance of normal breathing pressure is important. This will prevent any condensation from being forced through the filter.

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EN - BACTERIAL VIRAL HME FILTER

PERFORMANCE CHARACTERISTICS

Aspects	ThermoShield Bacterial / Viral HME Filter	ThermoShield Angled Bacterial / Viral HME Filter	ThermoShield Mini Bacterial / Viral HME Filter	ThermoShield Neonatal Bacterial/Viral HME Filter
Bacterial Filtration Efficiency (%)	99.99		99.999	99.999
Viral Filtration Efficiency (%)	99.9			99.99
Recommended Range of Tidal Volume (ml)	300-1000		135-500	55 - 500
Pressure Drop at Gas Flow Rate (cmH ₂ O)	Inspiratory (Before Usage): 0.98@30L/min 2.24@60L/min 3.96@90L/min Expiratory (Before Usage): 0.92@30L/min 2.11@60L/min 3.67@90L/min Inspiratory (After 24 hrs of Usage): 0.96@30L/min 2.21@60L/min 3.94@90L/min Expiratory (After 24 hrs of Usage): 0.95@30L/min 2.11@60L/min 3.67@90L/min	Inspiratory (Before Usage): 0.99@30L/min 2.31@60L/min 4.05@90L/min Expiratory (Before Usage): 0.92@30L/min 1.79@60L/min 2.90@90L/min Inspiratory (After 24 hrs of Usage): 1.05@30L/min 2.40@60L/min 4.22@90L/min Expiratory (After 24 hrs of Usage): 1.09@30L/min 2.05@60L/min 3.68@90L/min	Inspiratory (Before Usage): 0.89@15L/min 1.66@30L/min 2.76@45L/min Expiratory (Before Usage): 0.86@15L/min 1.51@30L/min 2.34@45L/min Inspiratory (After 24 hrs of Usage): 0.98@15L/min 1.61@30L/min 3.22@45L/min Expiratory (After 24 hrs of Usage): 1.05@15L/min 2.05@30L/min 2.68@45L/min	Inspiratory (Before Usage): 1.88@15L/min 4.97@30L/min 9.36@45L/min Expiratory (Before Usage): 1.95@15L/min 5.05@30L/min 9.42@45L/min Inspiratory (After 24 hrs of Usage): 1.77@15L/min 4.87@30L/min 8.78@45L/min Expiratory (After 24 hrs of Usage): 1.83@15L/min 4.97@30L/min 8.87@45L/min
Internal Volume (ml)	57.2	65.31	27.12	14.59
Particle Filtration Efficiency (%)	Before Usage: 98.542 After 24 hrs of Usage: 98.722	Before Usage: 98.542 After 24 hrs of Usage: 98.722	Before Usage: 96.372 After 24 hrs of Usage: 90.194	Before Usage: 78.333 After 24 hrs of Usage: 72.300
Gas Leakage (ml/min)	Before Usage: 0.5@70cmH ₂ O After 24 Hours of Usage: 0.0@70cmH ₂ O	Before Usage: 0.5@70cmH ₂ O After 24 Hours of Usage: 0.0@70cmH ₂ O	Before Usage: 0.5@70cmH ₂ O After 24 Hours of Usage: 0.0@70cmH ₂ O	Before Usage: 1.4@70cmH ₂ O After 24 Hours of Usage: 1.5@70cmH ₂ O
Moisture Output (mg/L)	Minimum Tidal Volume: 31.78 Operating Range (At 500ml): 30.9 Maximum Tidal Volume: 30.79	Minimum Tidal Volume: 31.78 Operating Range (At 500ml): 30.9 Maximum Tidal Volume: 30.79	Minimum Tidal Volume: 33.60 Operating Range (At 250ml): 31.9 Maximum Tidal Volume: 32.19	Minimum Tidal Volume: 35.39 Operating Range (At 250ml): 29.43 Maximum Tidal Volume: 23.47
Moisture Loss (mg/L)	Minimum Tidal Volume: 6.81 Operating Range (At 500ml): 7.84 Maximum Tidal Volume: 7.96	Minimum Tidal Volume: 6.81 Operating Range (At 500ml): 7.84 Maximum Tidal Volume: 7.96	Minimum Tidal Volume: 4.69 Operating Range (At 250ml): 6.67 Maximum Tidal Volume: 6.34	Minimum Tidal Volume: 2.58 Operating Range (At 250ml): 9.56 Maximum Tidal Volume: 16.53
Compliance (ml/cmH ₂ O)	0.056	0.071	0.029	0.037

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www.flexicare.com/symbols

Flexicare Medical Limited
Cynon Valley Business Park,
Mountain Ash, CF45 4ER, UK
+44 (0)1443 474647
enquiries@flexicare.com
www.flexicare.com

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Flexicare Deutschland GmbH