Voluson™ Performance Series Quick Card

Extraordinary vision

- | | | |
|---|--|--|
| 1 Power on/off | 9 Archive/exam review key | 19 Menu exit key |
| 2 USB connections | 10 Report key | 20 Arrow key |
| 3 F1/help key | 11 PW mode and gain control/
X and Y axis [†] | 21 Annotation/ABC key |
| 4 F2/utilities key | 12 M mode and gain control/Z axis [†] | 22 Clear key |
| 5 Rotary and paddle switch controls | 13 CF mode and gain control/
parallel shift [†] | 23 Measure key |
| 6 Acoustic power output/
Doppler volume | 14 2D mode and gain control | 24 3D/4D [*] key |
| 7 New patient/end exam | 15 HD/pan zoom | 25 Print/save keys |
| 8 Probe/program key | 16 Depth Control | 26 Trackball and trackball keys |
| | 17 Focal zone position control | 27 Freeze/run |
| | 18 Navigation rotary control | 28 TGC controls |

[†] 3D Control

^{*} Software option, may not be available in all countries

Entering patient information

1. Press the New Patient key (7).¹
2. Enter patient information. If a Patient ID is not entered, the system can assign an ID number. Patient ID is required to store images. If LMP, GA or EDD is not given, no growth information will be displayed (for growth trending on subsequent exams, the same ID number must be used).
3. Select Start from the top trackball key (26) or from the menu.
4. To end the exam, select New Patient key (7) then select End Exam from display menu.

Selecting probe and presets

1. Press the Probe/Program key (8). All connected probes are shown on the screen.
2. Press the corresponding keys to choose the probe and selections.

Measurements

1. Press the Measurement key (23).
2. Both generic measurements and calculations for the active application (e.g. OB) will appear on the screen.
3. To perform generic measurements, move and place calipers with trackball and right or left trackball key (26). Adjust endpoints with top trackball key as needed.
4. To perform calculations, select the desired calculation, then move and place calipers with trackball and right or left trackball keys.
5. To access the patient report, press Report (10).

3D/4D mode

1. Press the 3D/4D key (24).
2. From the display menu choose 3D or 4D scan mode and appropriate setting.
3. Position transducer in the center of the structure to be imaged.
4. Adjust ROI (region of interest) box position and size with trackball and top trackball key (26) to include area of interest. If performing a surface rendering, the green line should be just above the area of interest. A fluid interface is helpful when performing surface rendering.
5. Start acquisition using right trackball (26) or Freeze key (27). The acquisition is automatic; do not move the transducer during the acquisition.
6. Save volume with appropriate Print key (25).

Color mode

1. To activate press the CF knob (13).
2. To adjust gain, turn CF knob.
3. Adjust Color ROI position and size with trackball and top trackball key (26).
4. To switch from CF to PD or HD-Flow,[™] select the appropriate choice on the display menu.
5. PRF (pulse repetition frequency) and WMF (wall motion filter) can be adjusted with appropriate toggle keys (5).

PW mode

1. To activate press PW knob (11).
2. Adjust PW gate position and size with trackball and top trackball key (26).
3. Press left trackball key (26) to start spectral Doppler.
4. To adjust spectral gain, turn the PW knob.

Note:

To deactivate all modes, press the 2D mode (14).

Image export

1. Insert media.
2. Press Archive (9) on console.
3. Select desired patient and exam using right or left trackball key (Set).
4. Select Export from menu.
5. Choose format type:²
 - For patient use choose 'PC/MAC' format
 - For later use on Voluson or 4DView choose '.4dv uncompressed'
 - Many other selections available
6. Name file if other than the default provided is desired.
7. Choose anonymization if desired.
8. Save.

For more information please consult eIFU (Electronic Instructions For Use) via the F1 key (3).

1. Key numbers can be found on front side.
2. Format types: PC/MAC format files (.jpg, .mp4), Voluson format compressed (.4dv), Voluson format uncompressed (.4dv), JPEG files (.jpg), Bitmap files (.bmp), TIFF files (.tiff), AVI files (.avi), MPEG files (.mpg), Volume files/Raw Data (.vol, .raw), DICOM[®] files (.dcm), and DICOM files with DICOMDIR.

© 2020 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram, Voluson and HD-Flow are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. DICOM is a registered trademark of the National Electrical Manufacturers Association. GE Medical Systems, Inc., doing business as GE Healthcare.

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
The screenshot shows the GE Healthcare Voluson P6/P8 software interface. At the top, there are fields for Patient ID (E60003-14-07-22-1), Last Name, First Name, and Middle Name. Below these are fields for DOB (mm/dd/yyyy), Age, and Sex (female). A search button is also present. The main area contains a grid of buttons for different exam types: OB, GYN, ABD, CARD, TR, VAS, DM F, PED, CEPH, and MSK. The OB button is highlighted. Below the grid, there are fields for LMP, DOC, EDD, GA, Gravida, Para, AB, Ectopic, Fetus #, Perf. Physician, Ref. Physician, Sonographer, Exam Type, Exam Comment, Access #, Custom, Scan Assistant (None), and Indication. At the bottom, there is a table with columns: Exam Date, Exam Time, Exam Type, M, Img, Application, and Exam Size.

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