Quick Card

Voluson™ E10

- | | | |
|----------------------------------|--|--|
| 1 Power on/off | 9 Multi-function control knobs | 18 Menu exit key |
| 2 F1/help key | 10 Acoustic power output | 19 Arrow key |
| 3 Utilities key | 11 PW mode and gain control/X axis* | 20 Annotation/ABC key |
| 4 Report key | 12 M mode and gain control/Y axis* | 21 Clear key |
| 5 New patient | 13 PD mode and gain control/Z axis* | 22 Measure/Calc key |
| 6 Probe/program key | 14 CF mode and gain control/Parallel Shift* | 23 3D & 4D keys |
| 7 Archive/exam review key | 15 2D mode and gain control | 24 Auto Optimize |
| 8 End exam | 16 Push/Turn – HD/Pan Zoom
Up/Down – Depth Control | 25 Print, Store and Save keys |
| | 17 Turn – Sector Width
Up/Down – Focal Position
Push – Focal Number | 26 Freeze/run |
| | | 27 Trackball and trackball keys |

Voluson E10

Entering patient data

1. Press the "New patient" key (5)
2. Enter patient information. Patient ID is required to store images. Growth chart information will only be displayed if LMP, EDD or GA is input
3. Select "Start" from the top trackball key (27), touch panel or menu
4. To End exam, select End Exam key (8)

Probe/program selection

1. Press the "Probe/program" key (6)
All connected probes are shown on the screen
2. Press the corresponding keys to choose the probe and program selections

Calculation measurements

1. Press the "Measure/Calc" key (22). The calculation package for the active application (e.g. OB) will appear on the touch panel¹
2. Select the desired calculation
3. The calipers can be set by "Right or Left trackball" key (27)
4. For patient report, press "Report" key (4)

Generic measurements

1. Press the Measure/Calc key twice for generic measurements (22)¹
2. Select the desired measurement from the touch panel (e.g. generic dist). Entering and setting of calipers is done with the "Right or Left trackball" key (27)
Endpoints can be adjusted with "Top trackball" key; calipers are set with "Right or Left trackball" key

3D/4D mode

1. Press the "3D or 4D" key (23)
2. If choosing 4D, select the desired 4D mode on the top of the touch panel.
The right trackball key (27) or "Freeze" key (26) will begin the acquisition

CF mode (Color Doppler mode)

1. Press the "CF mode" key (14)
2. The CFM box appears in the active 2D image
3. The 'Top trackball' key (27) switches between CF box position and size
The 'CF mode' key (14) is also the gain control.
Activate Power Doppler or HD-flow in the same way with the PD key (13).
To switch between PD and HD-flow use touchpanel control.

PW mode (Pulsed Wave Doppler mode)

1. Press the "PW mode" key (11)
2. Adjust the PW cursor position and gate size with the trackball on the 2D image. The top trackball key (27) switches between gate position and gate width
3. Press the left trackball key (27) to start the spectral Doppler
4. The "PW mode" key (11) is also the PW spectrum gain control

Note:

To deactivate all modes, press the "2D mode" key (15).

Important note:

An Electronic User Manual (EUM) is integrated on the Voluson* Expert Series. Press the F1 (2) key on the keyboard to activate the EUM.

¹Measure/Calc Button can be configured in the Utility menu to either go into Generic or Calculations first.

© 2014 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.

GE, the GE monogram and Voluson are trademarks of General Electric Company or one of its subsidiaries.

GE Healthcare, a division of General Electric Company.

GE Healthcare
9900 Innovation Drive
Wauwatosa, WI 53226
U.S.A.
www.voluson.gehealthcare.com

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
The screenshot shows a patient information form with the following fields:

- Patient ID: MED10907-10-12-21-1
- Last Name: [Empty]
- First Name: [Empty]
- Middle Name: [Empty]
- 2nd Patient ID: [Empty]
- DOB (mm/dd/yyyy): [Empty]
- Age: [Empty]
- Sex: female
- Search button
- Tabs: ABD, OB, GYN, CARD, URO, VAS, NEURO, SMP, PED, MSK
- LMP: [Empty]
- Gravida: [Empty]
- Perf. Physician: [Empty]
- DOC: [Empty]
- Para: [Empty]
- Ref. Physician: [Empty]
- EDO: [Empty]
- AB: [Empty]
- Sonographer: [Empty]
- GA: [Empty]
- Ectopic: [Empty]
- Exam Type: [Empty]
- Fetus #: 1
- Exam Comment: [Empty]
- Access #: [Empty]
- Custom: [Empty]
- Scan Assistant: None
- Indication: [Empty]
- Table with columns: Exam Date, Exam Time, Exam Type, M, Img, Application, Exam Size

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The screenshot displays the IDS MED software interface. On the left is a sidebar with navigation options: Patient Archive, Archive, Image History, Exam Review, Current Patient, Worklist, Details, Case Entries, Worksheet, Home Page, and Exit. The main area contains a form for patient information, including Patient ID (MED10907-10-12-21-1), Last Name, First Name, Middle Name, DOB, Age, Sex (Female), and a Search button. Below this is a section for exam settings with tabs for ABD, OB, GYN, CARD, URO, VAS, NEURO, SM P, PED, and PDR. The OB tab is selected, showing fields for LMP, POC, EDD, GA, Provider, Para, AB, Ectopic, Fetus # (1), Ref. Physician, Ref. Physician, Sonographer, Exam Type, Exam Comment, Custom, and Indications. At the bottom, there is a table with columns: Exam Date, Exam Time, Exam Type, H, Img, Application, and Exam Size. Below the table is a 'Start Exam' button and a 'Set' button. On the right, there is a 'Run' section showing 'Run: 155 sec' and a progress bar.

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