

HOGY Surgical Gowns

Reuse Prohibited

Product name: Surgical Gowns

Size(s): M, L, LL and XL

Fabric: "Surrem" progress usual SMS is constituted 5 layers (SSMMS)



SG-SMS K Series

Intended use:

This gown is intended for use by medical personnel to protect both the wearer and the patient from the transfer of microorganisms, fluids and blood. The gown is appropriate for use in major procedures.

Performance and Standards:

The product complies with EN 13795-1:2019 Standard Performance

Contraindications

Allergic reaction to the materials used in gowns.

How to wear:

- 1) Before using the product, please check to ensure there is no damage to the carton or packaging. If the product or the packaging is damaged or wet, please dispose appropriately.
- 2) Carefully read the labelling on the front of the packaging paying careful attention to the

harmonised symbols, expiry date and any warnings. If the expiry date has passed, please dispose of the gown appropriately.

- 3) Wash hands thoroughly.
- 4) Open the packaging and remove the gown.
- 5) Don the gown by placing both arms through the sleeves with the open part of the gown facing backwards.
- 6) Secure the hook and loop fastening at the neck.
- 7) Tie the waist belt at the back using a secure bow or knot and ensure that the gown is overlapping so that your body is not exposed at the rear.
- 8) Don gloves (not supplied) and ensure that they are secured above the seam between the gown cuff and sleeve.
- 9) After use, carefully dispose of the gown by having an assistant undo the neck closure and the waist belt. Slide the gown forwards away from your body directly into the appropriate clinical waste bin ensuring that at no point your skin or clothing comes into contact with the outer layer of the gown.
- 10) Dispose of the gown as clinical waste in accordance with your organisation's guidelines.
- 11) Wash hands thoroughly after use

Warnings & Contraindications:

- 1) Do not use this product if you are allergic to polypropylene, polyester or cotton.
 - 2) Do not use this product after the expiry date has passed.
 - 3) Do not use this product if the packaging is opened, wet or visibly damaged.
 - 4) Do not expose product to flames, heat or high intensity light source.
- ISOL/8 Disposable SMS
Medical Gown Instructions for Use Version 003

Shelf Life

Three years from the date of manufacture.

Handling and Storage Conditions

Store in a dry, cool, ventilated environment away from flammable and corrosive materials. Protect from direct sunlight and moisture.

Transport conditions

During transport, avoid damage to the outer packaging, direct sunlight and rain, snow and other forms of precipitation.

Disposal Instructions:

Please dispose of gowns as clinically contaminated waste as per your organization's guidelines for waste disposal.

Sterilization:

Product is sterilized by EO.

<Contact information (etc.)>

[Manufacturer/distributor and manufacturer name (etc.)]

<Manufacturer/distributor>

PT. Hogy Medical
Address: Infinia Park No. 45 Blok A
57, Jalan Dr. Saharjo, Manggarai,
Tebet, Kota Jakarta Selatan,
Daerah Khusus Ibukota Jakarta
Phone: (021) 83797001

Department<Manufacturer>

HOGY MEDICAL CO.,LTD.