FasTurn Syringes / Seringues FasTurn / Jeringas FasTurn

Instructions for Use / Mode d'emploi / Instrucciones de Uso

Injector Injecteur Inyector	Catalog Number Référence catalogue Número de Catálogo	Product Name Nom du Produit Nombre del Producto
MEDRAD® Mark V Plus	150-FT	150 mL FasTurn Sterile Disposable Syringe / Seringue jetable stérile FasTurn de 150 mL / Jeringa desechable estéril FasTurn de 150 mL
	150-FT-Q	Sterile Disposable Syringe 150 mL & QFT / Seringue jetable stérile de 150 mL et canule de remplissage rapide / Jeringa desechable estéril de 150 mL y tubo de llenado rápido
MEDRAD® Mark V MEDRAD® Mark V ProVis	150-FT	150 mL FasTurn Sterile Disposable Syringe / Seringue jetable stérile FasTurn de 150 mL / Jeringa desechable estéril FasTurn de 150 mL
	150-FT-Q	Sterile Disposable Syringe 150 mL & QFT / Seringue jetable stérile de 150 mL et canule de remplissage rapide / Jeringa desechable estéril de 150 mL y tubo de llenado rápido
	200-FT	200 mL FasTurn Sterile Disposable Syringe / Seringue jetable stérile FasTurn de 200 mL / Jeringa desechable estéril FasTurn de 200 mL
	200-FT-Q	200 mL FasTurn Sterile Disposable Syringe & QFT / Seringue jetable stérile FasTurn de 200 mL et canule de remplissage rapide / Jeringa desechable estéril FasTurn de 200 mL y tubo de llenado rápido
MEDRAD® CT 100/200/300 Series Injectors Injecteurs séries 100/ 200/300 Serie de inyectoros 100/ 200/300	200-FT	200 mL FasTurn Sterile Disposable Syringe / Seringue jetable stérile FasTurn de 200 mL / Jeringa desechable estéril FasTurn de 200 mL
	200-FT-Q	200 mL FasTurn Sterile Disposable Syringe & QFT / Seringue jetable stérile FasTurn de 200 mL et canule de remplissage rapide / Jeringa desechable estéril FasTurn de 200 mL y tubo de llenado rápido

Bayer Medical Care Inc.

1 Bayer Drive
Indianola, PA 15051-0780
USA
Phone: +1 (412) 767-2400
+1 (800) 633-7231
Fax: +1 (412) 767-4120

Bayer Medical Care B.V.

Horsterweg 24
6199 AC Maastricht Airport
The Netherlands
Phone: +31 (0) 43-3585601
Fax: +31 (0) 43-3656598

Importer BMC B.V.:

Bayer Medical Care B.V.
Avenue Céramique 27
6221 KV Maastricht,
The Netherlands

© 2011-2014, 2017, 2019 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer.

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Bayer, the Bayer Cross, MEDRAD, Mark V ProVis, MEDRAD Mark V ProVis, FluiDots, MEDRAD FluiDots, and Quik-Fit Syringe are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries. Other trademarks and company names mentioned herein are properties of their respective owners and are used herein solely for informational purposes. No relationship or endorsement should be inferred or implied.

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A glossary of the symbols used on the listed sterile disposable products can be found on page 2 of these Instructions for Use.

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Introduction: Read the information contained in this section. Understanding the information will assist you in operating the device in a safe manner.

Important Safety Notice: This device is intended to be used by individuals with adequate training and experience in diagnostic imaging studies.

Indications for use: The contents of this package are intended to be used in the delivery of contrast media. These devices are indicated for single-use only with injectors from Bayer.

Contraindications: These devices are not intended for multiple-use, drug infusion, chemotherapy, or any other use for which the device is not indicated.

Restricted Sale: Federal (USA) law restricts these devices for sale by or on the order of a physician.

Report any serious incident that has occurred in relation to this device to Bayer (radiology.bayer.com/contact) and to your local European competent authority (or, where applicable, to the appropriate regulatory authority of the country in which the incident has occurred).

Warnings

For devices labeled for single use, please note: This product is intended for single use only. Do not resterilize, reprocess or reuse. The disposable devices have been designed and validated for single use only. Re-use of the single use disposable devices pose risks of device failure and risks to the patient. Potential device failure includes significant component deterioration with extended use, component malfunction, and system failure. Potential risks to the patient include injury due to device malfunction or infection as the device has not been validated to be cleaned or re-sterilized.

Air embolization can cause death or serious injury to the patient. Do not connect a patient to the injector until all trapped air has been cleared from the syringe and fluid path. Carefully read the instructions for loading and the use of MEDRAD® FluiDots indicators to reduce the chance of air embolism.

Biological contamination can result from reusing disposable items or failure to follow aseptic technique. Properly discard disposable items after use.

Syringe sterility will be compromised, and patient infection may result, if the plunger is removed from the syringe. Do not remove the plunger to fill the syringe.

Bacterial contamination can occur if syringes are used to store contrast media. Use loaded syringes immediately. Do not store loaded syringes for later use. Discard unused loaded syringes.

Do not use if sterile package is opened or damaged. Patient or operator injury may result if package is opened or damaged, or if damaged components are used. Visually inspect contents and package before each use.

Patient or Operator injury may result from contrast media leaks or tubing ruptures. Ensure that the fluid path is open; do not exceed pressures identified on the front of the package. Use of greater pressures or occlusions in the fluid path may result in leaks or ruptures to the connector tube.

Cautions

Component damage or leaks may occur if not installed properly. Ensure all connections are secure; do not overtighten. This will help minimize leaks, disconnection, and component damage.

Refer to the Injector Operation Manual for further instructions.

MEDRAD Mark IV, MEDRAD Mark V, MEDRAD Mark V Plus, MEDRAD Mark V Provis, or MEDRAD CT-100/200/300 Series Injectors:

1. Ensure that the piston on the head is fully retracted. To retract the piston(s), press the Reverse Load Switch on the injector head.
2. Set the Mechanical (safety) Stop to zero in

the full forward position (applies to MEDRAD Mark IV and MEDRAD CT-100 Series Injectors only). Position the empty Pressure Jacket in the top of the turret. If there is a syringe in the top pressure jacket, remove and discard.

3. Remove the syringe from its package and slide it in the top Pressure Jacket (the flat side of the flange must face upward). This positions the nut with its cutout facing up for easier hub attachment.

4. Fully rotate the turret to place the syringe plunger in front of the piston.

Note: The turret must be fully rotated to the closed position to ensure piston/syringe plunger alignment to prevent contrast leakage.

The following are general guidelines only. For additional instructions, refer to the Operation Manual for your Liebel Flarsheim injector.

Liebel Flarsheim 6000 and 6000 CT Injectors:

1. Ensure that the piston on the head is fully retracted.
2. Unlock pressure jacket and rotate downward.
3. Remove syringe from its package and insert syringe into pressure jacket, making sure to position alignment key with pressure jacket slot.
4. Rotate pressure jacket plate upward and latch by turning knob fully clockwise.
5. Advance piston to engage syringe plunger.

Liebel Flarsheim Angiomat Injector:

1. Press the LOAD button on the console. Press the LOAD button on the head to fully retract the plunger.
2. Set the mechanical stop at zero.
3. Open the turret door so the pressure jacket is ready to accept a syringe. If there is a syringe in the pressure jacket, remove and discard.
4. Remove syringe from its package and insert syringe into pressure jacket.
5. Close and secure the turret.
6. Press the UNLOAD button to advance the piston forward.

Syringe Filling:

Note: One person should fill and arm the injector. If you must change the injector operator, be certain that the second operator verifies that the syringe is properly filled and that any air is eliminated.

1. Tilt the injector head up.
 2. Press and hold the forward motion control until the piston is fully advanced.
 3. Remove the protective covers from the syringe luer tip and attach a sterile filling device onto the tip of the syringe (Quick Fill Tube or other filling device).
 4. Transfer the contrast media into the syringe using a Quick Fill Tube. Raise the contrast bottle until the Quick Fill Tube is fully inserted. The end of the Quick Fill Tube should be near the bottom of the contrast bottle.
- Note:** Use a Quick Fill Tube or equivalent device to reduce the volume and size of air bubbles drawn into the syringe during filling. It is more difficult to remove the air bubbles if you use smaller diameter tubes, such as a catheter-over-needle, or a tube longer than 10 in. (25cm.)
5. Aspirate contrast media into the syringe.
 6. Press and hold the REVERSE piston motion key until the syringe is filled with the volume of contrast media desired.
 7. Use the FluiDots indicators to help you detect the presence of contrast media in the syringe. Verify that the FluiDots indicators are round in the filled portion of the syringe. The rounded shape of the FluiDots indicators varies according to the type of contrast media, but an oblong shape indicates the presence of air. Rounded FluiDots indicators do not indicate the total absence of air bubbles in the syringe tip

5. Aspirate contrast media into the syringe.
6. Press and hold the REVERSE piston motion key until the syringe is filled with the volume of contrast media desired.
7. Use the FluiDots indicators to help you detect the presence of contrast media in the syringe. Verify that the FluiDots indicators are round in the filled portion of the syringe. The rounded shape of the FluiDots indicators varies according to the type of contrast media, but an oblong shape indicates the presence of air. Rounded FluiDots indicators do not indicate the total absence of air bubbles in the syringe tip

Expelling Air and Attaching the Connector Tube:

After you fill the syringe, you must expel all air from the syringe and then attach the connector tubing. To do so:

1. Remove the filling device and expel air from the syringe by manually advancing

the piston so that the contrast media is in the tip of the syringe.

Gently tap the base of the pressure jacket to dislodge the air bubbles. Expel remaining air.

2. Attach the connector tube to the syringe:
 - a. **FasTurn** - Place the connector tube hub in the syringe nut's cutout. Turn the syringe nut clockwise approximately 2 1/4 turns until secure while holding the hub stationary.
 - b. **Push-N-Turn** - Place the connector tube hub in the syringe nut's cutout. Push the syringe nut towards the injector head until it snaps; turn the nut clockwise approximately 1 turn until secure.

NOTE: Position the hubs so that when the syringe nut is secure, the wings are across the cutout, not with one wing in the cutout.

3. Verify that the tubing is not kinked or obstructed.
4. Use the manual knob to advance the syringe plunger to provide a slow flow of contrast media through the connector tubing.

An absence of flow is an obvious indication of air in the fluid path

5. Prime the connector tubing and make certain that the air is expelled.
6. Verify that all of the air has been expelled from the syringe and fluid path.
7. Attach the connector tube to the vascular entry device.
8. Tilt the injector head down.

Removing a Syringe:

To remove the syringe when the procedure has been completed, fully retract the piston and rotate the turret counter-clockwise. Remove the syringe from the pressure jacket. Discard all used components.

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