 bringing technology to life		WORK INSTRUCTION FM-SOP-20002-04 Rev. 100		Doc. Number AS10031622-001		Revision 100	
WI Type: ARTWORK SPECIFICATION				Effective Date: 25-SEP-2024		Page 1 of 37	
Title: INSTRUCTIONS FOR USE – IFU, ACAPELLA CHOICE BLUE VIBRATORY PEP DEVICE (10031622-001)							

DOCUMENT VERIFICATION

Prior to using this document, the user is responsible for verifying that the revision and effective date are current.

REVISION HISTORY

Owning Site: MASH		Additional Sites Using Document: MMSP, MDBL			
Rev.	Effective Date	Change Order No.	Changes Made to Document		
100	25-SEP-2024	CO-10166624	PWR: 5750 This artwork replaces; 10010097-004 (AS10010097-004 Revision 100). Updates include the following changes; • Artwork Specification Label content updates align with FM-DVSOP8005-14: • SM Legacy/Portex branding changed to ICU Medical/Portex branding. • Manufacturer updated to ICU Medical, Nathan Lane. • "www.smiths-medical.com" changed to "www.icumed.com". • EU Rep address changed to Dutch address • acapella® changed to acapella™ in all languages • SM Legacy trademark statement updated to ICU trademark statement. Additional changes to artwork: • Distribution Statement updated in all languages • Copyright/Trademark date updated to 2024 Glossary: • Extended English removed • Single Patient Multiple Use symbol and text added in all languages • Quantity symbol and text changed to Packaging Unit and text in all languages		

NOTES

- A. Part Number: 10031622-001
- B. For detailed IFU specifications refer to the purchase specification.
- C. Artwork Author: Ross Griggs

English

acapella™ Choice Blue Vibratory PEP Device.
Instructions for Use

These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Precautions, before using this product. Failure to properly follow warnings, precautions and instructions could result in death or serious injury to the patient and/or clinician.
NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.

1. DESCRIPTION:

The **acapella™** Choice Blue Vibratory PEP Device (**acapella™** Choice Blue) is a **single patient** multiple use device that provides Positive Expiratory Pressure (PEP) Therapy for patients who have Cystic Fibrosis, COPD, asthma, and lung diseases with secretory problems, and patients with atelectasis. All patients must be capable of following instructions for Positive Expiratory Pressure Therapy or be under adult supervision. Review the diagram and the product to become familiar with all of the product features. Following use of acapella, the expectoration of phlegm from the airways may be collected for analysis and or diagnosis per healthcare practitioner procedures. The **acapella™** Choice Blue Vibratory PEP Device with one-way inspiratory valve and 22 mm male fittings consists of: (See Figure 1)

A. Detachable mouthpiece

B. Expiratory resistance/frequency adjustment dial

C. Cover

D. Detachable rocker assembly – DO NOT disassemble the rocker arm from its platform base

2. INDICATIONS:

The ICU Medical **acapella™** is intended for use as a Positive Expiratory Pressure (PEP) device for adults and children (5 years and above). It may also be used simultaneously with nebulized aerosol drug delivery.

3. CONTRAINDICATIONS:

Although no absolute contraindications to the use of PEP Therapy have been reported, the following should be carefully evaluated before a decision is made to initiate therapy:

• Inability to tolerate increased work of breathing

• Hemodynamic instability

• Intracranial pressure (ICP) > 20 mm Hg

• Acute sinusitis

• Recent facial, oral or skull surgery or trauma

• Epistaxis

• Esophageal surgery

• Active hemoptysis

• Untreated pneumothorax

• Nausea

• Known or suspected tympanic membrane rupture or other middle ear pathology

4. WARNINGS:

4.1 Use of this device at excessive pressures may have adverse effects. Expiratory pressures above 20 cmH₂O in patients sensitive to increased transpulmonary pressure may develop one or more of the adverse side effects listed below.

4.2 Expert clinical judgment should be exercised in the selection of the appropriate setting for each individual patient. Failure to match the appropriate resistance setting on the +/- dial indicator with the patients expiratory flow may result in failure

to achieve therapeutic objectives of vibratory PEP therapy or one or more adverse side effects listed below.

Adverse reactions may include:

• Increased work of breathing that may lead to hypoventilation and hypercarbia

• Increased cranial pressure

• Cardiovascular compromise

- Myocardial ischemia

- Decreased venous return

• Air swallowing with increased likelihood of vomiting and aspiration

• Claustrophobia

• Pulmonary barotraumas

4.3 The rocker assembly is not designed to be disassembled - do not attempt to separate the rocker arm from its platform base. Separating these components may cause the acapella™ Choice Blue to function improperly, which may deprive the patient of therapy and cause injury to the patient.

5. PRECAUTIONS:

5.1 Bleach is not recommended for use on the **acapella™** Choice Blue. It may deteriorate the nickel plated mechanism located in the interior of the device.

5.2 DO NOT MICROWAVE. The metal and magnet might ignite.

5.3 Visually inspect the device to ensure that the unit is free of contamination and foreign objects.

5.4 Verify all connections are secure.

6. SUGGESTED INSTRUCTIONS FOR USE:

6.1 Initial Settings:

6.1.1 If this is the first use of the **acapella™** Choice Blue device, ensure that the frequency adjustment dial is turned to the lowest frequency-resistance setting, indicated at "1". (See Figure 2)

6.1.2 Instruct the patient to relax while performing diaphragmatic breathing and inspiring a volume of air larger than normal tidal volume (but not to total lung capacity).

6.1.3 Direct the patient to exhale to Functional Residual Capacity (FRC) actively, but not forcefully, through the device.

6.1.4 Adjust the dial to increase the resistance.

6.1.5 The patient should be able to exhale for 3-4 seconds while the device vibrates. Clockwise adjustment increases the resistance of the vibrating orifice; which will allow the patient to exhale at a lower flow rate. Counter-clockwise adjustment decreases the resistance.

6.1.6 Once the proper range has been identified, the patient may be instructed to exhale harder or softer, or dial adjustments may be made to optimize the response the user "feels" from the vibratory pressure.

6.2 Procedure for the User:

6.2.1 Ensure the adjustment dial is set to the correct range as identified by your clinician.

6.2.2 Sit with elbows resting comfortably on table.

6.2.3 Place mouthpiece lightly in mouth.

• Be sure to maintain a tight seal on the mouthpiece during exhalation.

• Your clinician may recommend the use of a nose clip, if necessary.

6.2.4 Breathe from the diaphragm, as directed by your clinician, taking in a larger than normal breath, but not filling your lungs to capacity.

6.2.5 Hold your breath for 2-3 seconds.

6.2.6 Exhale actively, but not forcefully, through the device. Exhalation should last approximately 3 to 4 times longer than inhalation.

6.2.7 Perform 10-20 PEP breaths as recommended by clinician. Remove mouthpiece and perform 2-3 "huff" coughs to raise secretions as needed. Your clinician may direct you on proper cough technique.

6.2.8 Repeat steps 6.2.2 to 6.2.7 as prescribed.

6.3 Set-up nebulizer:

6.3.1 Review the diagram contained with the devices.

6.3.2 Connect **acapella™** Choice Blue and nebulizer as seen in Figure 3.
Note: Only use the set up reflected in Figure 3, only using the nebulizer mouthpiece, when connecting a nebulizer to **acapella™** Choice Blue.

6.3.3 Follow nebulizer manufacture's set-up instructions for each device.

6.3.4 Follow nebulizer manufacture's cleaning instructions contained with each device.

6.3.5 Inspect device(s) on a routine basis to ensure proper use and function.

6.3.6 If damaged, do not use.

6.3.7 Verify all connections are secure.

7. SPECIFICATIONS:

acapella™ Choice Blue

	Flow Rate	
Parameters	15 L/min	25 L/min
Pressure (cmH ₂ O) across all settings	2.97 - 6.19	6.55 - 11.76
Frequency (Hz) across all settings	6.98 - 13.94	8.97 - 15.43
Amplitude (cmH ₂ O) across all settings	11.23 - 18.73	21.22 - 29.39
Max Pressure (cmH ₂ O) across all settings	12.45 - 20.66	23.10 - 31.35

acapella™ Choice Blue with a Nebulizer

	Flow Rate	
Parameters	15 L/min	25 L/min
Pressure (cmH ₂ O) across all settings	4.13 - 11.37	9.04 - 23.74
Frequency (Hz) across all settings	8.00 - 21.94	11.29 - 23.72
Amplitude (cmH ₂ O) across all settings	3.62 - 20.90	7.36 - 28.09
Max Pressure (cmH ₂ O) across all settings	11.40 - 24.46	17.74 - 36.00

Note: Testing was performed at a fixed flow rate. Actual values are dependent upon individual effort.

The following specifications were established via a performance test using an eight-stage cascade impactor at a flow rate of 28 L/min equipped with USP induction port throat. Three (3) device samples were tested per each drug. Aerosol was sampled directly from the outlet. The specifications are listed below based upon a 95% confidence level.

For information on use, cleaning, and disposal of nebulizer see the original manufacturer's instructions for use.

Nebulizer Specification

	Standalone Nebulizer	Nebulizer with acapella™ Choice Blue	Standalone Nebulizer	Nebulizer with acapella™ Choice Blue	Standalone Nebulizer	Nebulizer with acapella™ Choice Blue
Salter Labs® 8900 Series Nebulizer (operated at 6 L/min)						
	Albuterol Sulfate		Cromolyn Sodium		Ipratropium Bromide	
MMAD (µm)	2.63	2.37	3.40	2.87	3.37	3.50
GSD	3.26	2.60	4.42	3.97	3.99	4.81
Total Dose (µg)	625	529	2830	2307	132	152
Total Respirable Dose (0.5-5)	305	264	1198	1047	60	55
Coarse Particle (>4.7 µm)	272	211	1353	1034	65	87
Fine Particle (<4.7 µm)	353	317	1478	1272	67	65
Ultra-Fine Particle (<1.0 µm)	106	127	518	517	22	24
Carefusion AirLife™ Misty Max 10 (operated at 8L/min)						
	Albuterol Sulfate		Cromolyn Sodium		Ipratropium Bromide	
MMAD (µm)	2.73	2.47	2.97	2.67	2.23	2.10
GSD	3.34	3.60	3.63	3.53	3.62	2.77
Total Dose (µg)	865	792	4803	5039	151	149
Total Respirable Dose (0.5-5)	496	410	2637	2406	94	92
Coarse Particle (>4.7 µm)	319	324	1908	2216	47	50
Fine Particle (<4.7 µm)	545	469	2895	2822	104	99
Ultra-Fine Particle (<1.0 µm)	165	170	863	981	32	28
Westmed UniHEART™ (operated at 4 L/min)						
	Albuterol Sulfate		Cromolyn Sodium		Ipratropium Bromide	
MMAD (µm)	1.63	1.37	1.40	1.10	1.87	1.83
GSD	2.52	2.74	2.96	3.06	3.62	3.13
Total Dose (µg)	538	472	1766	1817	97	87
Total Respirable Dose (0.5-5)	357	320	1019	1099	48	43
Coarse Particle (>4.7 µm)	102	70	351	348	39	32
Fine Particle (<4.7 µm)	437	402	1415	1469	58	55
Ultra-Fine Particle (<1.0 µm)	185	188	752	823	25	24

Note: Coarse particles (oropharyngeal deposition) and ultra-fine particles (exhaled) are not likely to deposit in the patient's airway and thus provide clinical benefit. The statistical significant values in the table above are based on the statistical significant analysis for the total respirable dose.

1

2

This document and all information contained herein is confidential, proprietary and may not be disclosed, disseminated or reproduced without written approval from Smiths Medical.

- INTERNAL -

Controlled Copy - Verify Revision & Effective Date are current before use

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69

70

 bringing technology to life		WORK INSTRUCTION FM-SOP-20002-04 Rev. 100		Doc. Number AS10031621-001		Revision 100	
WI Type: ARTWORK SPECIFICATION				Effective Date: 25-SEP-2024		Page 1 of 3	
Title: IFU, ACAPELLA VIBRATORY PEP THERAPY SYSTEM (10031621-001)							

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REVISION HISTORY

Owning Site: MASH		Additional Sites Using Document: MMSP, MDBL			
Rev.	Effective Date	Change Order No.	Changes Made to Document		
100	25-SEP-2024	CO-10166624	PWR: 5750 This artwork replaces 10018905-001 (AS10018905-001 Rev100) Updates include the following changes; • Artwork Specification Label content updates align with FM-DVSOP8005-14: • SM Legacy/Portex branding changed to ICU Medical/Portex branding. • Manufacturer updated to ICU Medical, Nathan Lane. • "www.smiths-medical.com" changed to "www.icumed.com". • EU Rep address changed to Dutch address • acapella® changed to acapella™ in all language • SM Legacy trademark statement updated to ICU trademark statement. Additional changes to artwork: • Title format for each language removed from front panel, moved to start of each language, adjusted to standard IFU Format (Language Identifier/Product Title/ Instructions For Use) • DEHP Statement updated per FM-DVSOP8005-03 • Distribution Statement updated in all languages • Copyright/Trademark date updated to 2024 • Single Patient Use statements removed Glossary changes: • Order of symbols corrected per FM-DVSOP8005-03 • Single Patient Multiple Use symbol and text added in all languages • Quantity symbol and text changed to Packaging Unit and text in all languages		

NOTES

- A. Part Number:** 10031621-001
- B.** For detailed IFU specifications refer to the purchase specification.
- C. Artwork Author:** Ross Griggs

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REVISION HISTORY

Owning Site: MASH		Additional Sites Using Document: MDBL	
Rev.	Effective Date	Change Order No.	Changes Made to Document
100	25-SEP-2024	CO-10166624	PWR: 5750 This artwork replaces 10001049-001 (DWG-10001049-001 - INSERT Revision 100). Updates include the following changes; • Artwork Specification Label content updates align with FM-DVSOP8005-14; • SM Legacy/Portex branding changed to ICU Medical/Portex branding. Additional changes to artwork: • Drawing Template converted to WI FM-SOP-20002-04 Rev. 100 format • Final sentence 'Smiths Medical's declarations regarding SVHCs in products can be found at www.smiths-medical.com/environment.' translation removed from each language (Confirmation email attached to CO) • English language description removed for each translation


NOTES

- A. Part Number:** 10031626-001
- B.** For detailed Insert specifications refer to the purchase specification.
- C. Artwork Author:** Ross Griggs

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 smiths medical bringing technology to life	WORK INSTRUCTION FM-SOP-20002-04 Rev. 100	Doc. Number AS10031626-001	Revision 100
WI Type: ARTWORK SPECIFICATION		Effective Date: 25-SEP-2024	Page 2 of 3
Title: INSERT – INSERT, REACH STATEMENT (10031626-001)			

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smiths medical bringing technology to life		WORK INSTRUCTION FM-SOP-20002-04 Rev. 100		Doc. Number AS10031621-001		Revision 100	
WI Type: ARTWORK SPECIFICATION				Effective Date: 25-SEP-2024		Page 1 of 3	
Title: IFU, ACAPELLA VIBRATORY PEP THERAPY SYSTEM (10031621-001)							

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REVISION HISTORY

Owning Site: MASH		Additional Sites Using Document: MMSP, MDBL			
Rev.	Effective Date	Change Order No.	Changes Made to Document		
100	25-SEP-2024	CO-10166624	PWR: 5750 This artwork replaces 10018905-001 (AS10018905-001 Rev100) Updates include the following changes; • Artwork Specification Label content updates align with FM-DVSOP8005-14: • SM Legacy/Portex branding changed to ICU Medical/Portex branding. • Manufacturer updated to ICU Medical, Nathan Lane. • "www.smiths-medical.com" changed to "www.icumed.com". • EU Rep address changed to Dutch address • acapella® changed to acapella™ in all language • SM Legacy trademark statement updated to ICU trademark statement. Additional changes to artwork: • Title format for each language removed from front panel, moved to start of each language, adjusted to standard IFU Format (Language Identifier/Product Title/ Instructions For Use) • DEHP Statement updated per FM-DVSOP8005-03 • Distribution Statement updated in all languages • Copyright/Trademark date updated to 2024 • Single Patient Use statements removed Glossary changes: • Order of symbols corrected per FM-DVSOP8005-03 • Single Patient Multiple Use symbol and text added in all languages • Quantity symbol and text changed to Packaging Unit and text in all languages		

NOTES

- A. **Part Number:** 10031621-001
- B. For detailed IFU specifications refer to the purchase specification.
- C. **Artwork Author:** Ross Griggs

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[illegible]

 bringing technology to life		WORK INSTRUCTION FM-SOP-20002-04 Rev. 100		Doc. Number AS10031621-001		Revision 100	
WI Type: ARTWORK SPECIFICATION				Effective Date: 25-SEP-2024		Page 1 of 3	
Title: IFU, ACAPELLA VIBRATORY PEP THERAPY SYSTEM (10031621-001)							

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REVISION HISTORY

Owning Site: MASH		Additional Sites Using Document: MMSP, MDBL			
Rev.	Effective Date	Change Order No.	Changes Made to Document		
100	25-SEP-2024	CO-10166624	PWR: 5750 This artwork replaces 10018905-001 (AS10018905-001 Rev100) Updates include the following changes; • Artwork Specification Label content updates align with FM-DVSOP8005-14: • SM Legacy/Portex branding changed to ICU Medical/Portex branding. • Manufacturer updated to ICU Medical, Nathan Lane. • "www.smiths-medical.com" changed to "www.icumed.com". • EU Rep address changed to Dutch address • acapella® changed to acapella™ in all language • SM Legacy trademark statement updated to ICU trademark statement. Additional changes to artwork: • Title format for each language removed from front panel, moved to start of each language, adjusted to standard IFU Format (Language Identifier/Product Title/ Instructions For Use) • DEHP Statement updated per FM-DVSOP8005-03 • Distribution Statement updated in all languages • Copyright/Trademark date updated to 2024 • Single Patient Use statements removed Glossary changes: • Order of symbols corrected per FM-DVSOP8005-03 • Single Patient Multiple Use symbol and text added in all languages • Quantity symbol and text changed to Packaging Unit and text in all languages		

NOTES

- A. Part Number:** 10031621-001
- B.** For detailed IFU specifications refer to the purchase specification.
- C. Artwork Author:** Ross Griggs

[illegible]

[illegible]