Intra-Aortic Balloon Catheter Kit

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tr İnteraortik Balon Kateter Kiti

This page left blank intentionally

en Intra-Aortic Balloon Catheter Kit

This Instruction for Use provides information for use of Arrow® Intra-Aortic Balloon (IAB) Catheters.

The following chart provides recommendations for choosing the IAB volume based on patient size. This information should be used as a guideline only as clinical studies have shown only limited correlation of patient height to aortic length. Clinical studies suggest that age (especially age > 70), ethnicity, body stature, patient sex and vascular status should be considered for IAB size selection. The clinician should determine appropriate IAB size based on good practice and hospital policies and procedures:

	30 cc	40 cc	50 cc
	147 cm - 162 cm 4'10" - 5'4"	162 cm - 182 cm 5'4" - 6'0"	> 182 cm > 6'0"
BSA 	< 1.8 m ²	> 1.8 m ²	> 1.8 m ²

Note: The Seldinger Technique is recommended for the percutaneous insertion of the IAB catheter. T.J. Conahan, A.J. Schwartz, R.T. Geer. "Percutaneous Catheter Introduction: The Seldinger Technique." J. Am. Med. Assoc. 237,1977.

The IAB Tray in conjunction with the Insertion Tray includes nearly everything required for percutaneous insertion. Standard surgical insertion technique can also be employed.

Safety and Efficacy Considerations:

Warning: Carefully read all instructions prior to use.

Precaution: Persons using these devices should be appropriately trained, skilled personnel. Training appropriateness is the determination of the prescribing physician.

The clinician must weigh benefits versus risks involved in employing this device based on current standards of good medical practice.

Intended Use of Device:

The Arrow® IAB is utilized for intra-aortic balloon counterpulsation therapy in the aorta, whereby balloon inflation, during diastole and deflation, during systole increases blood supply to the heart muscle and decreases work of the left ventricle.

Indications for Use:

The Arrow® Intra-Aortic Balloon (IAB) Catheter with the Intra-Aortic Balloon Pump as a control system is clinically indicated for use in any of the following conditions:

- Acute Coronary Syndrome
- Cardiac and Non-Cardiac Surgery
- Complications of Heart Failure

Contraindications:

Hemodynamically significant aortic regurgitation
Dissecting aortic aneurysm
Severe clotting disorders
Severe atherosclerotic peripheral vascular disease
End stage terminal disease

Potential Complications and Adverse Events:

- Vessel occlusion resulting in Limb Ischemia
- Aortic wall damage including aortic rupture or dissection
- Thrombotic complications, including:
 - Thrombus

- Emboli
 - Thrombocytopenia
 - Hemolysis
- Infection
 - Failure of the IAB membrane
 - IAB Rupture (Abrasion resulting in leak)
 - IAB entrapment due to clot in IAB making IAB removal difficult or requiring surgical removal
 - Improper positioning or sizing, resulting in;
 - Possible damage to aortic wall
 - Obstruction of arterial origins with related organ ischemia
 - Ineffective cardiac assist
 - Bleeding at insertion site or in other organs, resulting in a decrease in Hemoglobin and Hematocrit, requiring transfusion.

Warnings:

- Sterile, Single use:** Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.
- Thrombectomy may be indicated in the event of compromised distal circulation.
- If blockage or excessive resistance is met, use contralateral groin. Do not force. This may result in arterial tearing, dissection or balloon damage.
- An elevated blood pressure can expel dilator from sheath. Do not leave unattended.
- The provided sheath should not remain indwelling without internal support of a catheter or similar device to maintain patency.
- Do not use excessive force when inserting or removing balloon; arterial tearing, dissection, or balloon membrane damage may result.
- Do not withdraw guidewire against needle bevel to minimize the risk of damaging or severing the guidewire. Do not cut guidewire.
- Exercise extreme care in setup and flushing of central lumen arterial pressure monitoring line to minimize potential for introducing an embolus into the aortic arch. Stop counterpulsation prior to sampling or fast flushing. Never manually flush central lumen with a syringe. Should resistance be met, consider the line no longer patent and discontinue its use.
- Intra-aortic balloon membrane perforation may occur during IABP therapy. The occurrence and severity of IAB perforation is unpredictable and may be due to patient physiology, accidental contact with a sharp instrument or by contact with calcified plaque resulting in membrane surface abrasion and eventual perforation. Large perforations are rare. Small perforations can result in asymptomatic release of gas. Perforation can cause blood to appear in the balloon catheter and driveline tubing.
- If you suspect balloon perforation:
 - Immediately stop pumping. Consider tapering or discontinuing anticoagulation therapy.
 - Remove IAB from patient, using recommended removal technique.
- Leg ischemia can occur either during or after IABP therapy. It can be caused by thrombus formation, an intimal flap or layer separation. Persistent leg ischemia following balloon removal may require vascular surgical correction.

- Do not remove Arrow® IAB through hemostasis sheath introducer or hemostasis device. Once unwrapped (unfurled), balloon profile will not allow passage through the sheath and attempting removal in this manner may result in arterial tearing, dissection or balloon damage.
 - Always confirm proper balloon position with fluoroscopy and/or chest radiograph.
 - Augmentation that fails to fall within normal range within a few minutes of the start of counterpulsation may indicate:
 - Balloon positioned incorrectly.
 - Hemostasis sheath introducer positioned incorrectly.
 - Hemodynamically compromised patient condition.
 - Balloon membrane not fully opening (possible kink in catheter or membrane retained in sheath).
 - Never inject air into central lumen (female Luer fitting) to minimize the risk of an air embolus.
 - Balloon insertion into a tortuous or severely restricted vessel may prevent proper balloon unwrapping. If this should occur, attach supplied syringe to balloon Luer fitting. Aspirate and check for blood in driveline tubing. If no blood is seen, manually inflate and deflate balloon to its labeled volume.
 - Avoid using solvents (i.e. alcohol or acetone) in skin preparation as solvents may weaken or damage IAB components.
 - Improper balloon sizing may result in injury to patient or damage to balloon.
 - Do not use automatic dye injectors to instill dye into central lumen. This may result in damage to the central lumen necessitating IAB removal.
 - Failure to properly secure the IAB catheter may result in IAB migration. In order to maintain proper catheter positioning, follow the manufacturer's recommended methods and techniques.
 - The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.
- ### Precautions:
- Do not remove IAB from tray until it is to be inserted. The tray is designed to protect the IAB surface from damage or contamination, and maintain a tight wrap.
 - The one-way valve will maintain vacuum on the balloon and must remain in place until balloon is fully inserted.
 - To reduce the risk of kinking guidewire during insertion, advance guidewire through needle and into descending thoracic aorta before placing other devices, such as dilators, sheaths or IAB catheter.
 - Do not use excessive force during insertion of the guidewire or of the IAB over the guidewire. This could result in damage to the guidewire.
 - If the guidewire is compromised or damaged, use an alternate guidewire.
 - When inserting Arrow® IAB catheter into hemostasis sheath introducer, blood may occasionally seep along folds of the balloon. This is not a balloon leak. This is caused by arterial pressure forcing blood between folds of the IAB and wall of the sheath.

- If IABP therapy must be continued following removal of the IAB, a percutaneous insertion may be performed in contralateral femoral artery.
- A high arterial pressure can expel balloon; hold balloon in place until ready to remove it.
- The central lumen should be connected to a standard heparinized flush solution. This may reduce incidence of occlusion due to blood coagulation. The use of a transducer to monitor the central lumen is optional.
- Introduction of the IAB without use of a hemostasis sheath introducer is not recommended in patients with severe obesity, scarring of the groin or other contraindications to percutaneous insertion, where the distance between the skin and femoral artery exceeds 2 inches (5 cm).
- Use of the IAB beyond 30 days may increase the risk of IAB catheter kinks and balloon abrasion. The clinician should determine the risks vs. benefits of extended IAB use, as compared to IAB exchange after 30 days of use.
- Rx only.

A Suggested Procedure: Use sterile technique.

There are 3 options for insertion of the IAB: (1) sheathless, Note: The angle of insertion should not exceed 45 degrees. It is recommended to keep insertion angle as parallel to the leg as possible. This angle of insertion should be maintained throughout entire insertion; (2) using supplied hemostasis sheath introducer, or (3) if using an indwelling hemostasis sheath introducer is already in place in femoral artery, proceed to Step 4, SHEATHED INSERTION.

INSERTION SITE PREPARATION AND ACCESS:

Insertion may be performed in any location that is conducive to maintenance of a sterile field during insertion and appropriate post-insertion care.

1. Using standard preparation for Seldinger Technique, prep and drape groin that has the strongest pulse.
2. Anesthetize insertion site.
3. With a sterile scalpel blade, nick skin over common femoral artery.
4. With arterial needle, puncture anterior wall of common femoral artery.
5. Pass supplied guidewire through needle and advance guidewire to the aortic arch.
6. Remove needle, leaving guidewire in place.
7. Pass dilator over guidewire, through the skin and subcutaneous tissue, and into artery.
8. Remove dilator, leaving guidewire in place.
9. Clean blood from guidewire with a lint-free sponge.

Warning: If blockage or excessive resistance is met, use contralateral groin. Do not force. This may result in arterial tearing, dissection or balloon damage.

Precaution: To reduce the risk of kinking guidewire during insertion, advance guidewire through needle and into descending thoracic aorta before placing other devices, such as dilators, sheaths or IAB catheter.

Precaution: Do not use excessive force during insertion of the guidewire or of the IAB over the guidewire. This could result in damage to the guidewire.

Warning: An elevated blood pressure can expel dilator from sheath. Do not leave unattended.


Warning: Do not withdraw guidewire against needle bevel to minimize the risk of damaging or severing the guidewire. Do not cut guidewire.

Precaution: Do not remove intra-aortic balloon from tray until it is to be inserted. The tray is designed to protect the IAB surface from damage or contamination, and maintain a tight wrap.


IAB PREPARATION:

1. Connect one-way valve to male Luer on short driveline tubing attached to IAB. Insert supplied syringe into one-way valve. Slowly aspirate a full syringe of air.

Precaution: The one-way valve will maintain vacuum on the balloon and must remain in place until balloon is fully inserted. Remove syringe.


2. Remove catheter from tray, keeping it in line with balloon membrane.

3. Grasp catheter close to tray and pull it straight out of the holding sleeve.

Note: Keep catheter level with tray. Do not lift or bend during removal.


4. Remove stylet from central lumen (if applicable).
5. Flush central lumen with heparinized saline.
6. Wet balloon membrane with sterile saline prior to insertion to activate hydrophilic coating.

Note: The coating may become tacky immediately after wetting and becomes slippery within 8 to 12 seconds.

7. Proceed to Instructions for Use: a) for SHEATHLESS INSERTION or b) for SHEATHED INSERTION.

SHEATHLESS INSERTION:


1. Insert guidewire through central lumen opening in tip of IAB until it exits female Luer fitting on IAB.

2. Advance IAB over guidewire, through skin into artery, while controlling proximal end of guidewire. Grasp IAB close to skin line during insertion to minimize the risk of kinking the catheter. Advance IAB to proper position in descending thoracic aorta 1 - 2 cm distal to left subclavian artery. Place IAB as high above femoral bifurcation as possible.

Warning: If blockage or excessive resistance is met, use contralateral groin. Do not force. This may result in arterial tearing, dissection or balloon damage.

Precaution: If the guidewire is compromised or damaged, use an alternate guidewire.

3. Follow Steps 14 - 22 to complete insertion.

SHEATHED INSERTION:

1. Assemble sheath by inserting dilator into hemostasis sheath introducer. If using sheath with sidearm, flush sidearm with heparinized saline.


Although it is recommended to use sheath introducers supplied in the Insertion Kit, Arrow® IABs have been qualified for use with other manufacturer's introducers. Follow Arrow® IAB labeling for appropriate sheath selection and review sheath manufacturer's labeling and instructions for use.

2. Pass sheath and dilator over guidewire and advance into artery. Advance hemostasis sheath introducer until 1 inch (2.5 cm) remains outside insertion site. This will allow for easier handling and better control of bleeding.

Warning: If blockage or excessive resistance is met, use contralateral groin. Do not force. This may result in arterial tearing, dissection, or balloon damage.

Precaution: If the guidewire is compromised or damaged, use an alternate guidewire.

- After achieving vascular access, hold hemostasis sheath introducer in position and remove introducer dilator.

Note: When inserting IAB catheter through a sheath already in place in femoral artery, the clinician should always review sheath manufacturer's labeling and instructions for use.


If using a sheath with sidearm, aspirate 2 - 3 mL (cc) of blood from sidearm of sheath and then flush with heparinized saline solution. Using current hospital protocol, connect sidearm to a prepared standard arterial pressure monitoring assembly which delivers 3 mL (cc) of pressurized flush per hour.

Note: Aspirate 2 - 3 mL (cc) of blood prior to sampling blood or flushing sidearm to minimize the risk of embolus. If unable to aspirate blood, discontinue use of sidearm.


- To use an existing indwelling sheath, remove any catheters or devices inserted through the sheath per device manufacturer's Instructions for Use. Exchange guidewire with IAB supplied guidewire. Only supplied guidewire in IAB insertion kit should be used for IAB insertion. Some bleeding may occur through hemostasis valve when guidewire is used with another manufacturer's hemostasis sheath introducer or hemostasis device. Bleeding will subside once IAB catheter has been inserted into hemostasis valve and hemostasis sheath introducer or hemostasis device.

- Connect one-way valve to male Luer on short driveline tubing attached to IAB. Insert supplied syringe into one-way valve. Slowly aspirate a full syringe of air.

Precaution: The one-way valve will maintain vacuum on the balloon and must remain in place until balloon is fully inserted. Remove syringe.


- Remove catheter from tray, keeping it in line with balloon membrane.

- Grasp catheter close to tray and pull it straight out of the holding sleeve.

Note: Keep catheter level with tray. Do not lift or bend during removal.


- Remove stylet from central lumen (if applicable).

- Flush central lumen with heparinized saline.

- Remove peel-away hemostasis device as follows:

- Grasp wings of device leaving distal end of catheter pointing away.
- While pressing down along center of device with your thumbs, pull up on wings of device with your index fingers.


- Repeat this motion in other direction until hub is split in two.
- Peel the two halves of hemostasis device apart to remove.

- Wet balloon membrane with sterile saline prior to insertion to activate hydrophilic coating.

Note: The coating may become tacky immediately after wetting and becomes slippery within 8 to 12 seconds.

- Insert guidewire through central lumen opening in tip of IAB until it exits female Luer fitting on IAB.


Note: If using an indwelling hemostasis sheath introducer, backload supplied guidewire through central lumen of IAB catheter until J-tip protrudes slightly from tip of IAB. Draw guidewire back into IAB catheter approximately 1 inch (2.5 cm).


- Advance IAB into sheath while controlling proximal end of guidewire. Grasp IAB close to hub of introducer during insertion to minimize the risk of kinking the catheter.

Advance IAB to proper position in descending thoracic artery 1 - 2 cm distal to left subclavian artery. Place IAB as high above femoral bifurcation as possible.


Note: IAB catheter markings (if applicable) indicate relative position of the IAB. The double band indicates where IAB membrane has exited the 6" (15 cm) sheath. The IAB should always be inserted beyond this mark. For sheaths longer than 6" (15 cm), add the difference in sheath length to insertion depth. Below double band, markings indicate 1 cm distance, with numerical markings at every 5 cm distance. Catheter markings should be used as a relative indication of IAB position. Confirm IAB position with fluoroscopy and/or chest radiograph.


Warning: If blockage or excessive resistance is met, use contralateral groin. Do not force. This may result in arterial tearing, dissection, or balloon damage.

14. Maintaining IAB position, remove guidewire. Connect pressure tubing extension to balloon central lumen.


15. Aspirate 3 - 4 mL (cc) of blood and flush central lumen carefully with heparinized saline. Using current hospital protocol, connect pressure tubing extension to a prepared standard arterial pressure monitoring assembly which delivers 3 mL (cc) of pressurized flush per hour.
16. Slide and lock hemostasis cuff over hemostasis sheath introducer hub or hemostasis device hub. If using an alternate sheath, slide hemostasis cuff up to hub.


17. Remove one-way valve. Attach driveline tubing with female Luer fitting to male fitting. Twist Luer fittings to tighten.

Precaution: Twisting Luer fittings together reduces the risk of disconnection.

Warning: Always confirm proper balloon position with fluoroscopy and/or chest radiograph.

18. Attach opposite end of driveline tubing to control system using appropriate control system adapter. Follow control system manufacturer's instructions for detailed balloon purge and start-up procedure.
19. Suture both hemostasis cuff and proximal tie down to skin. Apply a sterile occlusive dressing to exposed catheter shaft per hospital protocol.

If using StatLock® Stabilization Device, refer to the Instructions for Use for manufacturer's recommendations. Both the hemostasis cuff and proximal tie down should be secured as described in Step 19, whether using sutures or an alternative securing device.

In the event that the balloon does not properly unwrap, disconnect driveline tubing from IAB catheter. Attach supplied 60 mL (cc) syringe into male end on driveline tubing attached to IAB catheter. Aspirate and check for blood in driveline tubing. If no blood is seen, manually inflate and deflate balloon to its labeled volume. Relieve any kinks or restrictions of sheath and IAB catheter at insertion site by pulling sheath back a half inch or rotating sheath a quarter turn.


Warning: Never inject air into central lumen (female Luer fitting) to minimize the risk of an air embolus.

Note: The Arrow® IAB in this package is equipped with a clear sleeve to allow repositioning of the IAB under clean conditions. Do not remove clear sleeve unless blood appears inside it.

When clinically necessary, the IAB may be repositioned by grasping IAB through clear sleeve and advancing unit to desired position.

Precaution: Do not attempt to reposition IAB with proximal tie down sutured to skin. Remove suture first. If IAB repositioning is necessary, use of a guidewire is recommended.

Warning: Confirm repositioning with fluoroscopy and/or chest radiograph.

FOR SHEATHLESS INSERTION ONLY:

20. Apply manual pressure to insertion site for 1 to 2 minutes to control bleeding. If bleeding persists, advance tip of pre-mounted hemostasis device into femoral artery up to taper, and then advance an additional 1 - 1.5 cm. The hemostasis device is designed for use whenever manual pressure is insufficient to achieve hemostasis. The shortest length of hemostasis device required to effectively control bleeding at insertion site should be inserted. Inserting entire hemostasis device could reduce the benefits of sheathless insertion since the hemostasis device will occupy space in the artery.

21. Attach cuff to hemostasis device.

22. If not using hemostasis device to control bleeding, either a) remove device or b) slide hemostasis device to insertion site skin line.

Intra-Aortic Balloon Removal Procedure:

Note: Successful removal of the Arrow® IAB depends upon adequate weaning procedures and assessment of the patient, usually over a 12 hour period prior to balloon removal.

Assess distal circulation several times before and after balloon removal. Periodic assessment of distal circulation is recommended, both during and after IAB counterpulsation.

If distal perfusion has been compromised during counterpulsation, or if percutaneous removal cannot be easily performed, surgical removal of percutaneous balloon may be indicated.

Discontinuance or reduction of anticoagulant therapy prior to balloon removal is highly recommended.

When removing balloon from patient, remove balloon and hemostasis sheath introducer or hemostasis device as a unit.

Warning: Do not remove the Arrow® IAB through a hemostasis sheath introducer or hemostasis device. Once unwrapped, balloon profile will not allow passage through the introducer and attempting removal in this manner may result in arterial tearing, dissection, or balloon damage.

BALLOON REMOVAL:

1. Carefully remove dressing from insertion site. Following standard hospital protocols, prep and drape insertion site for balloon removal. Remove all sutures and ties anchoring catheter to skin.

Precaution: A high arterial pressure can expel balloon; hold balloon in place until ready to remove it.

2. Turn pump drive OFF and disconnect balloon from control system.
3. Apply firm digital pressure to femoral artery distal to insertion site. This will help to minimize the potential for distal embolization when catheter is removed.
4. Grasp the Arrow® IAB and remove. If a hemostasis sheath introducer or hemostasis device has been used, withdraw the IAB until balloon comes in contact with sheath. A slight resistance will be encountered. Grasp the IAB and hemostasis sheath introducer or hemostasis device and remove as a unit.
5. Maintaining pressure distal to arterial puncture site, allow free flow of blood to flush puncture site for several seconds. Apply firm pressure proximal to arterial puncture site, and allow retrograde blood flow to flush puncture site for several seconds more.
6. After flushing arterial puncture site thoroughly, apply firm pressure to site for 30 minutes, or until hemostasis has been achieved.
7. Check distal pulses frequently and assess for signs of complications.
8. After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.

Supplemental Instructions for Use: Intra-Aortic Balloon Catheter with FiberOptix™ (Fiber Optic Arterial Pressure) Sensor

This Supplemental Section of Insertion Instructions for Use is specifically related to setup and use of the FiberOptix™ (Fiber Optic Arterial Pressure) Sensor.

Note: All other instructions pertaining to preparation and insertion of the IAB catheter should be followed.

Product codes that end in “LWS” indicate the IAB has a FiberOptix™ Sensor at the tip of the IAB.

Warnings:

- **Zero FiberOptix™ Sensor prior to IAB insertion. The FiberOptix™ Sensor cannot be zeroed after IAB insertion. Failure to properly zero FiberOptix™ Sensor may affect accuracy of FiberOptix™ Sensor arterial pressure values.**
- **The FiberOptix™ Sensor system is not suitable for use in the presence of flammable anesthetics.**
- **Use of FiberOptix™ Sensor with equipment other than Arrow® AutoCAT2 WAVE®/ AC3 Optimus™ IABP may lead to erroneous arterial pressure measurements or damage to FiberOptix™ Sensor and electronics. Only connect FiberOptix™ Sensor with Arrow® AutoCAT2 WAVE®/ AC3 Optimus™ IABP systems.**
- **Do not re-zero FiberOptix™ Sensor during use. This may affect accuracy of FiberOptix™ Sensor.**
- **Use of the FiberOptix™ Sensor at altitudes above 10,000 feet may result in inaccurate pressure values or loss of the FiberOptix™ Sensor signal. Have an alternative Arterial pressure source available for monitoring if this environmental condition is expected to occur.**

Precautions:

- **Use only the CAL key attached to the IAB. The CAL key is unique to each FOS sensor. Use of another CAL key may result in no FOS signal or inaccurate arterial pressure readings.**
- **Leave CAL key connected at all times. If CAL key is removed but same CAL key is reinserted, calibration and zero values are retained. If another CAL key is inserted, all previous calibration and zero values will be lost and FiberOptix™ Sensor reading accuracy may be affected.**
- **When connecting the Blue slider of the IAB into the FOS Slider into the (Silver) on the IABP, do not use excessive force during the connection as this may recess the FOS connector, resulting in an inability to establish the FOS AP signal. If resistance is met during the connection, pull the IAB slider back, ensuring it is straight within the slide connector. Then re-try the connection.**
- **CAUTION: Handle the FOS cable carefully as the fiber can be damaged by kinking, cutting or extreme bending of the cable. Do not touch the internal connections of the FOS sensor! Dirt, oils or other contaminants deposited on the connection can affect the light transmission. This may affect the ability to produce and display the AP waveform or the accuracy of hemodynamic values.**
- **DO NOT REMOVE the IAB from the tray until the FOS zero is complete (if the FOS sensor will be zeroed prior to insertion). This will maintain a tight wrap and allow easy passage through the sheath and reduce the risk of damage to the IAB or FOS fiber.**

- **PULL AND MAINTAIN vacuum on the IAB to reduce stress on the IAB catheter and FOS cable. This improves insertability by maintaining the smallest IAB wrap diameter for insertion into the sheath and reduces the risk of damaging the fiber in the IAB during removal from the tray.**
- **Do not use excessive force during IAB removal from the package and do not allow excessive bending of the IAB Catheter as it is removed from the tray. Keep the IAB as straight as possible and parallel to the tray during removal. This will reduce stress on the IAB and FOS cable.**

If the IAB is difficult to remove, verify the one-way valve is still on the Driveline tubing. It may be necessary to pull additional vacuum to tighten the IAB wrap.

Connection and Zeroing of FiberOptix™ Sensor to IABP System:

1. Prior to removing the IAB from the package, sterile user should pass the FiberOptix™ Sensor connector and CAL key to the pump Operator.

Note: Avoid manipulating the IAB during the zeroing process as this may affect the ability of the system to zero or affect the accuracy of the measured pressure values.

2. Connect CAL key to IABP

CAUTION: Use only the CAL key attached to the IAB. The CAL key is unique to each FOS sensor. Use of another CAL key may result in no FOS signal or inaccurate arterial pressure readings.

3. Slide the BLUE FOS connector into IABP receptacle until it clicks into place. A two beep tone will be heard when the sensor is fully connected and recognized by the IABP.

CAUTION: When connecting the Blue slider of the IAB into the FOS Slider into the (Silver) on the IABP, do not use excessive force during the connection as this may recess the FOS connector, resulting in an inability to establish the FOS AP signal. If resistance is met during the connection, pull the IAB slider back, ensuring it is straight within the slide connector. Then re-try the connection.

CAUTION: Handle the FOS cable carefully as the fiber can be damaged by kinking, cutting or extreme bending of the cable. Do not touch the internal connections of the FOS sensor! Dirt, oils or other contaminants deposited on the connection can affect the light transmission. This may affect the ability to produce and display the AP waveform or the accuracy of hemodynamic values.

4. After CAL key and FiberOptix™ Sensor are connected to IABP, while IAB tip is exposed to room air and before IAB is inserted, FiberOptix™ Sensor should be zeroed. This may occur automatically when sensor is connected and in room air. If pump is able to zero automatically, a message will appear on screen to indicate sensor is zeroed and a 4 beep tone will be heard. The FiberOptix™ Sensor icon will change from blue to green or the Zero icon will show zeroed.

Note: The IAB tip is exposed to air in the package. It is not necessary to remove the IAB from the tray to zero the FiberOptix™ Sensor.

Note: If the FOS sensor is NOT zeroed prior to use, the FOS Icon will show that the FOS sensor is not zeroed. Use an alternative AP source for treatment decisions. While this step is highly recommended for accuracy of the hemodynamic readings, it is not required for accurate WAVE® timing.

5. Pull vacuum on the IAB as described in IAB Preparation Section 1.

CAUTION: Pull and maintain vacuum on the IAB to reduce stress on the IAB catheter and FOS cable. This improves insertability by maintaining the smallest IAB wrap diameter for insertion into the sheath and reduces the risk of damaging the fiber in the IAB during removal from the tray.

6. Once FiberOptix™ Sensor is zeroed, IAB may be inserted.
7. Remove the IAB from the package keeping the IAB straight. Avoid bending the IAB as this may damage the IAB catheter or FOS fiber.

CAUTION: Do not use excessive force during IAB removal from the package and do not allow excessive bending of the IAB Catheter as it is removed from the tray. Keep the IAB as straight as possible and parallel to the tray during removal. This will reduce stress on the IAB and FOS cable.

If the IAB is difficult to remove, verify the one-way valve is still on the Driveline tubing. It may be necessary to pull additional vacuum to tighten the IAB wrap.

8. Follow remaining IAB instructions to position and secure IAB to patient and initiate IABP counterpulsation.

Note: All FOS IAB catheters have a fully functional central lumen and maintaining patency of the lumen is highly recommended to reduce the risk of clot formation at the IAB tip and to have a back-up AP source if needed.

Note: If the FOS AP signal is not available, connect the central lumen to a transducer and Zero the transducer. Connect the transducer to the IABP. AutoPilot mode will automatically select this waveform when it is available.

Removing Intra-Aortic Balloon with FiberOptix™ Sensor:

1. Disconnect CAL key and FiberOptix™ Sensor Connector from console with IABP volume connector.
2. Follow remaining instructions for IAB removal in main part of Instructions for Use.

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Corporate Office

Phone +1 610 225 6800, 550 E. Swedesford Road, Suite 400, Wayne, PA 19087, USA

Regional Offices

United States: Phone +1 919 544 8000, Toll Free 866 246 6990, cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

Latin America: Phone +1 919 433 4999, la.cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

International: Phone +353 (0)9 06 46 08 00, orders.intl@teleflex.com, Teleflex Medical Europe Ltd., IDA Business and Technology Park, Dublin Road, Athlone, Co Westmeath, Ireland

Australia 1300 360 226

Austria +43 (0)1 402 47 72

Belgium +32 (0)2 333 24 60

Canada +1 (0)800 387 9699

China (Shanghai) +86 (0)21 6163 0965

China (Beijing) +86 (0)10 6418 5699

Czech Republic +420 (0)495 759 111

France +33 (0)5 62 18 79 40

Germany +49 (0)7151 406 0

Greece +30 210 67 77 717

India +91 (0)44 2836 5040

Italy +39 0362 58911

Japan +81 (0)3 6632 3600

Korea +82 2 536 7550

Mexico +52 55 5002 3500

Netherlands +31 (0)88 00 215 00

New Zealand 0800 601 100

Portugal +351 22 541 90 85

Singapore (SEA non-direct sales countries)
+65 6439 3000

Slovak Republic +421 (0)3377 254 28

South Africa +27 (0)11 807 4887

Spain +34 918 300 451

Switzerland +41 (0)31 818 40 90

United Kingdom +44 (0)1494 53 27 61

For more information, please visit teleflex.com

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 · USA

Customer Service Tel:
1-919-544-8000
Email: cs@teleflex.com

Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

StatLock is a trademark or registered trademark of C.R. Bard, Inc.

Teleflex, the Teleflex logo, Arrow, the Arrow logo, AutoCAT2 WAVE, AC3 Optimus, FiberOptix, and WAVE are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries.

The products in this instruction booklet may not be available in all countries. Please contact your local representative.

© 2020 Teleflex Incorporated. All rights reserved.

PN B108-4800-130E Rev. 00 (08/20)