

VADI

Swivel flex tube

REF G-322012

Item Description : Reusable Catheter Mount

**Contents : 1x elbow connector 1x silicone tube
 1x straight connector**

Usage:

1. Connect silicone tube to the elbow connector.
2. Connect another end of silicone tube to the straight connector.
3. Connect the whole set catheter mount the endotracheal tube and the wye piece of breathing circuit.
4. Perform the leak test of whole set of breathing circuit according to the individual operating manual of ventilator.

Cleaning And Sterilization:

1. Dismantle the elbow connector, silicone tube, and straight& wye connectors.
2. Clean the parts with neutral detergents.
3. Dry them out by heated air or cool air .
4. Sterilization:
 - (1)Autoclave: 132°C / 4 minutes
 - (2)Pasteurization

VADI Medical Technology Co., Ltd. Yangmei

No. 198, Lane 298, Huandong Road, Zhongshan Village,
Yangmei District, Taoyuan City 32665, Taiwan

TEL: +886-3-475-9911

E-mail: sales@vadimedical.com.tw

FAX: +886-3-475-9922

Website: www.vadimedical.com.tw



Please see the seal.



Please see the seal.



Please see the seal.



Made In Taiwan
VA-1-G-322012.CFS_01