



Revision date: 17/11/2021

Revision number: 00/ Nov 2021

Bowie Dick Test Pack

Instructions For Use

Intended use:

Bowie dick test pack consists of a series of steam penetration (air removal) barriers in the center of which is a chemical indicator sheet. The color change of indicator sheet determines the effectiveness of pre vacuum sterilizer in removing air.

Directions for use:

1. The test pack is placed directly into an otherwise empty steam sterilizer chamber, with no retaining devices required.
2. During processing, the cycle must remove or displace the air from barrier material and replace it with steam throughout the pack.
3. A uniform change from blue to black indicates adequate steam penetration.

Shelf life:

2 years

Legal manufacturer:

Global SCP Medical Systems Pte. Ltd., Singapore