flexicare
TOTAL QUALITY - TOTAL CARE

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EN- STANDARD DUAL LIMB VENTILATOR BREATHING SYSTEMS

Flexicare Standard Dual Limb Ventilator Breathing Systems are intended for use on ventilated Adult, Paediatric and Infant/Neonatal patients by suitably trained/qualified personnel to connect a patient's airway to a Ventilator.The user should be familiar and trained in the selection of the appropriate Ventilator Breathing System, the use and the care of the product, and should be aware of any clinical contraindications.

Precautions
- Do not use if packaging is open, or packaging or contents are damaged, or after the expiry date.
- Before use, check all connectors are securely in place.
- Carry out a functional test to ensure there is no occlusion or leakage in the system.

Directions For Use
These instructions should be followed before each use in association with hospital policy and local guidelines:
1. Remove the patient end red safety cap before connecting to the patient airway device.
2. Use a push and twist action to connect components, and a twist and pull action to disconnect.
3. For Breathing Systems with one or two Self-sealing Water Traps, monitor the level of condensate and empty container as required.
4. For use with a Respiratory Humidifier Heater for Active Humidification, please refer to the manufacturer's Instructions for Use and ensure all necessary components are available.

Warnings / Cautions
- Do not reuse this device. Reuse may result in hazards such as, but not limited to; cross contamination, damage to device and/or suboptimal performance.
- The Breathing System is designed for use for up to a maximum period of 7 days, or less if dictated by local policies and guidelines.
- The Breathing System should be replaced immediately if it becomes visibly soiled, occluded, contains a high level of condensate (rain out) or develops a leak.
- Do not use a Heat and Moisture Exchanger (HME), Heat Moisture Exchanger with Filter (HMEF) or Pleated Paper Filter between a patient and a humidification device.
- This product may not be decontaminated in any way. This includes rinsing, washing or decontamination using gas, heat, steam, or boiling water.
- For Breathing Systems with one or two Self-sealing Water Traps, ALWAYS ensure Water Trap(s) are positioned at the lowest point at the limb(s) and are below the level of the patient's head.
- ALWAYS treat condensate as infectious material.
- DO NOT cover or surround the Breathing System with bedding or other materials.
- TAKE CARE not to stretch, occlude or crush the Breathing System limbs.

Any serious incident that occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State of the user and/or patient.

Protective Filters
- Ventilators may be further protected from contamination by a validated* bacterial / viral filter at the machine end of the system, fitted to either the inspiratory, expiratory or both machine ports.
- An appropriately validated* bacterial / viral filter may also be used at the patient end of the system.
- A Heat and Moisture Exchanger (HME) or Heat Moisture Exchanger with Filter (HMEF) may only be used on systems without humidification.
- The bacterial / viral filter must be changed as dictated by the manufacturer's Instructions for Use and/or hospital policy. Change filter immediately if ventilation is compromised or significant pressure drop or rise is observed.
- All filters should be used in accordance with the manufacturer's instructions.
- When using nebulised drugs, resistance to flow should be closely monitored and the filter replaced after treatment as a precaution against blockage.
- DO NOT use a filter between the patient and the nebuliser.

*A validated bacterial/viral filter meets BS EN ISO 23328 Breathing system filters for anaesthetic and respiratory use. Parts 1 & 2.

Record Keeping & Expiry Date
- Each Ventilator Breathing System is supplied with an adhesive label for record keeping and traceability. Complete the date, time, change date, and lot number at the time of opening and attach the label to the Breathing System at the machine end.
- Do not use after the end of the expiry date printed on the packaging or after the change date noted on the adhesive label.

Storage
Store in a cool, dry place out of direct sunlight.

Disposal
After use, dispose of in accordance with local policies and / or as advised by your local healthcare provider.

Technical Data
Please contact Flexicare® for any additional technical data.

Visit www.flexicare.com/downloads for symbols guide and information on the risks associated with the re-use of single use devices.

All information in this leaflet is correct at the time of printing.

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Flexicare Medical Limited
Flexicare Deutschland GmbH