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CE 0123

Table of Contents

1. Symbols ····· 1
1.1 Control Buttons ····· 1
1.2 Device Symbols ····· 1
2. Warning, Caution and Important Tip ····· 3
3. Intended Use ····· 3
4. Contraindications ····· 4
5. Clinical benefit ····· 5
6. Specifications ····· 5
7. Available Therapies ····· 8
8. Glossary ····· 9
9. Model ····· 11
10. Package Contents ····· 12
11. System Features ····· 13
12. First Time Setup ····· 14
12.1 Placing the Device ····· 14
12.2 Installing the Air Filter and Filter Cap ····· 15
12.3 Connecting Power Supply ····· 15
12.4 Connecting to Power Cord Locker ····· 16
12.5 Assembling the Tube and Mask ····· 17
12.6 Using Oxygen with the Device ····· 17
12.7 Inserting the SD Card (Only for the device that equipped with SD card) ····· 18
12.8 Using the Optional Kits ····· 19
12.9 Starting Treatment ····· 19
13. Routine Use ····· 20
13.1 Connecting the Tube ····· 20
13.2 Adjusting the Tube ····· 20
13.3 Turning on the Airflow ····· 20
13.4 Heating the Water ····· 20
13.5 Using the Ramp Button ····· 20
13.6 Turning the Device Off ····· 20
14. Heated Humidifier ····· 21
14.1 Filling the Water Chamber ····· 21
14.1.1 Removing the Water Chamber ····· 21
14.1.2 Filling the Water Chamber ····· 21
14.1.3 Inserting the Water Chamber ····· 22
14.2 Emptying the Water Chamber ····· 23
14.3 Setting the Humidity Level ····· 23
15. Navigating the Patient Menu ····· 24
15.1 Steps to Navigate the Patient Menu ····· 24
15.1.1 Accessing the Main Interface ····· 24
15.1.2 Bringing up the Initial Setup Interface ····· 25
15.1.3 Accessing the Setup Interface ····· 25
15.1.4 Selecting Options ····· 25
15.1.5 Adjusting Options ····· 26
15.1.6 Confirming Adjustments ····· 26
15.1.7 Turning Pages ····· 26
15.1.8 Exiting the Patient Menu ····· 27
15.2 Options in the Patient Menu and Corresponding Descriptions ····· 28
16. Alarm ····· 29
16.1 Alarm Grades and Description ····· 29
16.2 Visual Alarm ····· 29
16.3 Auditory Alarm ····· 29

16.4 Alarm Silence ····· 30
16.5 Alarm Messages and Description ····· 30
16.6 Reposition of Alarming ····· 33
16.7 Alarm Journal ····· 33
17. Cleaning and Disinfection ····· 33
17.1 Cleaning the Mask and Headgear ····· 34
17.2 Cleaning the Optional Kits ····· 34
17.3 Cleaning the Water Chamber ····· 34
17.4 Cleaning the Tube ····· 35
17.5 Replacing the Air Filter ····· 35
17.6 Cleaning the SpO_2 Probe ····· 35
17.7 Disinfection ····· 35
17.7.1 Disinfecting the mask and headband ····· 36
17.7.2 Disinfection of the SpO_2 Probe ····· 36
18. Traveling with the Device ····· 37
19. Transferring the Device to another Patient ····· 37
20. Reordering ····· 37
21. Technical Support ····· 38
22. Disposal ····· 38
23. Troubleshooting ····· 39
23.1 Common Problems in Patients and Corresponding Solutions ····· 39
23.2 Common Problems in the Device and Corresponding Solutions ····· 41
24. EMC Requirements ····· 42
25. Limited Warranty ····· 45

1. Symbols

1.1 Control Buttons

 Ramp Button

 Mute Button

 Knob

1.2 Device Symbols

 Follow Instructions for Use

 Consult Instructions for Use

 Type BF Applied Part (mask and SpO_2 probe)

 Class II (Double Insulated)

 For indoor use only

 AC Power

 DC Power

IP22 ≥ 12.5 mm Diameter, Dripping (15º tilted)

 There is high voltage, beware of electric shock

 Hot Surface

 Serial Number

 Manufacturer

 Authorized Representative in the European Community

 Do not use if package is damaged and consult instructions for use

 Maximum Water Level

 Disassembly is prohibited

CE 0123	CE Marking
1	Single Patient Multiple Use
LOT	Batch Code
	Non-Ionizing Radiation
SD	SD Card
	Air Inlet
	Air Outlet
	Caution
	Marking of Electrical and Electronic Equipment
BMC	Logo of BMC Medical Co., Ltd.
MD	Medical Device
UDI	Unique Device Identifier
#	Model Number
CN	Made in China
70°C -25°C	Temperature Limit
93% 15%	Humidity Limitation
1060hpa 760hpa	Atmospheric Pressure Limitation
	Importer Information
SpO2	Alarm Inhibit (Low SpO_2 alarm is not available)

2. Warning, Caution and Important Tip

WARNING!

Indicate the possibility of injury to the user or operator.

CAUTION!

Indicate the possibility of damage to the device.

IMPORTANT TIP!

Indicate the possibility that such operation may affect the effectiveness or ease of use of the device.

Warnings, Cautions, and Important Tips appear throughout this manual as they apply.

3. Intended Use

BPAP System (G2S BPAP Series) is a Bi-level PAP (Bi-level Positive Airway Pressure) device, which is intended to provide non-invasive ventilation for patients with obstructive Sleep Apnea (OSA) and Respiratory Insufficiency. These devices are intended for adult patient by prescription in the home or hospital/institutional environment.

The device is to be used only on the instruction of a licensed health care professional. Your home care provider will make the correct pressure settings according to your health care professional's prescription.

WARNINGS!

• This device is intended for adult use only.

• This device is not intended for life support.

• The instructions in this manual are not intended to supersede established medical protocols.

• Do not bring the device or accessories into a Magnetic Resonance (MR) environment as it may cause unacceptable risks to the patient or damage to the device or MR medical devices. The device and accessories have not been evaluated for safety in an MR environment.

• Do not use the device or accessories in an environment with electromagnetic equipment such as CT scanners, Diathermy, RFID and electromagnetic security systems (metal detectors) as it may cause unacceptable risks to the patient or damage to the device. Some electromagnetic sources may not be apparent, if you notice any unexplained changes in the performance of this device, if it makes unusual or harsh sounds, disconnect the power cord and stop using it. Contact your home care provider.

• Don't near active HF SURGICAL EQUIPMENT and the RF shielded room of an ME SYSTEM for magnetic resonance imaging, where the intensity of EM DISTURBANCES is high.

• Several accessories are available to make your treatment with this device as convenient and comfortable as possible. To ensure that you receive the safe and effective therapy prescribed for you, use only the manufacturer accessories.

• The optional SpO_2 Kit used with the device together is intended to be used to obtain readings of SpO_2 value as part of routine checkups.
• Any serious incidents in relation to this device should be reported to the manufacturer and the competent authority in your country.

CAUTIONS!

• This device is restricted to sale by or on the order of a physician.
• The device is intended for use by operators trained or experienced in similar equipment.
• The patient is an intended operator.
• Cleaning and disinfection can be performed by the patient.

IMPORTANT TIP!

• Read and understand the entire user manual before operating this system. If you have any questions concerning the use of this system, contact your home care provider or health care professional.

4. Contraindications

If you have any of the following conditions, tell your doctor before using this device:
• Insufficient respiratory drive to endure brief interruptions in non-invasive ventilation therapy
• Acute sinusitis or otitis media
• Epistaxis causing a risk of pulmonary aspiration
• Conditions predisposing to a risk of aspiration of gastric contents
• Impaired ability to clear secretions
• Hypotension or significant intravascular volume depletion
• Pneumothorax or pneumomediastinum
• Recent cranial trauma, cerebrospinal fluid leak or surgery
• Apparent non-cooperation or extreme stress
• Dehydration

The following side effects may occur during treatment:
- Dryness of the mouth, nose and throat
- Abdominal bloating
- Ear or sinus discomfort
- Eye irritation
- Skin irritation due to the use of a mask
- Chest discomfort

CAUTION!

• Contact your health care professional if symptoms of sleep apnea recur. Contact your health care professional if you have any questions concerning your therapy.

IMPORTANT TIPS!

• An irregular sleep schedule, alcohol consumption, obesity, sleeping pills, or sedatives may aggravate your symptoms.

• Please use the mask which meets ISO 17510:2015 and ISO 18562.

5. Clinical benefit

(1) Relieve symptoms in patients with Respiratory Insufficiency.

(2) Relieve symptoms in patients with sleep apnea and hypopnea.

(3) The clinical benefit of humidification is a reduction of noninvasive ventilation related side effects.

6. Specifications

Device Size

Dimensions: 274 mm × 184 mm × 115 mm

Weight: 1.9 kg

Water capacity: To maximum fill line 360 mL

Product Use, Transport and Storage

	Operation	Transport and Storage
Temperature:	5°C to 35°C (41°F to 95°F)	-25°C to 70°C (-13°F to 158°F)
Humidity:	15% to 93% Non-condensing	15% to 93% Non-condensing
Atmospheric Pressure:	760 to 1060 hPa	760 to 1060 hPa

Heated Humidifier

Humidifier Settings: Off, 1 to 5 (95°F to 154.4°F / 35°C to 68°C)

Humidifier Output: No less than 10 mg H_2O/L

Environmental Conditions: Maximum airflow, 35°C, 15% relative humidity

Maximum Delivered Gas Temperature: ≤ 43°C

Mode of Operation

Continuous

Work Mode

CPAP, S, AutoS, AutoCPAP, S/T, T

SD Card

The SD card can record usage data and fault information.

AC Power Consumption

100 V - 240 V ~, 50 Hz / 60 Hz, Max 2 A

Main Device offer to Communications Port

5 V ⎓ 2.0 A

Type of Protection against Electric Shock

Class II Equipment

Degree of Protection against Electric Shock

Type BF Applied Part

Degree of Protection against Ingress of Water

IP22

Pressure Range

IPAP: 4.0 ～ 20.0 hPa (only applies to G2S B20A, G2S B20S, G2S B20T);
4.0 ～ 25.0 hPa (only applies to G2S B25A, G2S B25S, G2S B25T, G2S B25VT, G2S LAB);
4.0 ～ 30.0 hPa (only applies to G2S B30T, G2S B30VT, G2S B30AT);
in 0.5 hPa increments.

EPAP: 4.0 ～ 20.0 hPa (only applies to G2S B20A, G2S B20S, G2S B20T);
4.0 ～ 25.0 hPa (only applies to G2S B25A, G2S B25S, G2S B25T, G2S B25VT, G2S B30T, G2S B30VT, G2S B30AT, G2S LAB);
in 0.5 hPa increments.

CPAP and AutoCPAP mode: 4.0 ～ 20.0 hPa

Under single fault conditions, ≤ 30 hPa for CPAP mode, ≤ 40 hPa for the rest modes.

Pressure Display Accuracy

± (0.8 hPa＋4%)

Static Pressure Stability

± 0.5 hPa

Ramp

The ramp time ranges from 0 to 60 minutes.

Sound Pressure Level

< 30 dB(A), when the device is working at the pressure of 10 hPa, uncertainty: 2 dB(A).

Sound Power Level

< 38 dB(A), when the device is working at the pressure of 10 hPa, uncertainty: 2 dB(A).

Maximum Flow

For models with a maximum pressure of 20 hPa, the flow table is as follows.

Test Pressure (hPa)	4	8	12	16	20
Measured Pressure at the Patient Connection Port (hPa)	3	7	11	15	19
Average Flow at the Patient Connection Port (L/min)	90	150	150	150	145
When the working pressure is set to the values listed in the table, the average flow rate at the patient end should be greater than 80% of the corresponding flow value in the table.					

For models with a maximum pressure of 25 hPa, the flow table is as follows.

Test Pressure (hPa)	4	10	15	20	25
Measured Pressure at the Patient Connection Port (hPa)	3	9	14	19	24
Average Flow at the Patient Connection Port (L/min)	90	150	150	145	110
When the working pressure is set to the values listed in the table, the average flow rate at the patient end should be greater than 80% of the corresponding flow value in the table.					

For models with a maximum pressure of 30 hPa, the flow table is as follows.

Test Pressure (hPa)	4	11	17	24	30
Measured Pressure at the Patient Connection Port (hPa)	3	10	16	23	29
Average Flow at the Patient Connection Port (L/min)	90	150	150	120	85
When the working pressure is set to the values listed in the table, the average flow rate at the patient end should be greater than 80% of the corresponding flow value in the table.					

SpO_2

Range: 0 ～ 100%

The margin of error for SpO_2 between 70% and 100% is ±3%. No strict accuracy requirements for SpO_2 below 70%.

Pulse Rate

Range: 40 ～ 240 BPM

Margin of Error: ±1%

Wavelengths

Red: 663 nanometers

Infrared: 890 nanometers

Maximal Optical Output Power

Less than 1.5 mW maximum average.

Automatic leak compensation

The pressure error ≤ ± 1 hPa when the air leakage is no more than 70 L/min.

Tube

Length: 6 ft. (1.8 m)

The Form and the Dimensions of the Patient Connection Port

The 22 mm conical air outlet complies with ISO 5356-1.

Power of heating components

38 W

7. Available Therapies

The device delivers the following therapies:

CPAP – Delivers Continuous Positive Airway Pressure; CPAP maintains a constant level of pressure throughout the breathing cycle.

AutoCPAP – Delivers CPAP therapy and provides an air pressure no less than the prescribed one based on the patient's needs.

AutoS – A bi-level mode which responds to both your inhalation and exhalation. The differential pressure of IPAP and EPAP are presetted by a home care provider. While working in auto mode, the device will automatically adjust the IPAP and EPAP if it detects a sleep apnea.

S – A bi-level mode which responds to both your inhalation and exhalation by increasing pressure when you start to inhale and decreasing pressure when you start to exhale. There is no automatic delivery of breathing gas if you do not inhale. IPAP (Inspiratory Positive Airway Pressure) and EPAP (Expiratory Positive Airway Pressure) are preset by a home care provider.

S/T – A bi-level mode which responds to both your inhalation and exhalation by increasing pressure when you start to inhale and decreasing pressure when you start to exhale. If you do not start inhaling within a set time, the device will automatically start the process of inhalation. When the device starts the process of inhalation, it controls the time of inhalation and automatically decreases the pressure for exhalation within a set time.

T – A bi-level mode in which the device automatically starts inhalation and exhalation, and automatically controls the time of inhalation and that of exhalation according to the preset parameter.

8. Glossary

Apnea

A condition marked by the cessation of spontaneous breathing.

AutoCPAP

Adjust CPAP pressure automatically to improve patient comfort based on monitoring of apnea and snoring events.

Auto Off

When this feature is enabled, the device automatically discontinues therapy whenever the mask is removed.

Auto On

When this feature is enabled, the device automatically initiates therapy when you breathe into the mask.

CPAP

Continuous Positive Airway Pressure.

EPAP

Expiratory Positive Airway Pressure.

IPAP

Inspiratory Positive Airway Pressure.

iCode

A feature designed to give access to compliance and therapy management information. "iCode" consists of six separate codes displayed in the Patient Menu, each code being a sequence of numbers. "iCode QR" and "iCode QR+" display two-dimensional codes.

LPM

Liters Per Minute.

OSA

Obstructive Sleep Apnea.

Patient Menu

The display mode in which you can change patient-adjustable device settings, such as the starting pressure of the Ramp feature.

Ramp

A feature that increases patient comfort at the beginning of treatment. It begins from a low pressure and then gradually increases to the prescribed setting pressure so that the patient can fall asleep more comfortably.

Rise Time

The time it takes for the device to change from EPAP to IPAP. You can adjust this time for your comfort.

Res Rate

Respiratory Rate. Number of breaths per minute.

Reslex

A therapy feature that is enabled by your home care provider to provide pressure relief during exhalation.

Standby State

The state of the device when power is applied but the airflow is turned off.

min

Means the time unit "minute".

h

Means the time unit "hour".

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Denotes date.

9. Model

Model	Product Contents		Work Mode	Maximum Work Pressure (hPa)
	Main Device	Optional Accessories		
G2S B20S	Main device (2.4-inch TFT)	Tube, Mask, Cellular Module, SpO_2 Kit, SpO_2&GPRS Kit	CPAP, S	20
G2S B20A	Main device (2.4-inch TFT)		CPAP, S, AutoS	20
G2S B20T	Main device (2.4-inch TFT)		CPAP, S, S/T, T	20
G2S B25S	Main device (3.5-inch TFT)		CPAP, S	25
G2S B25A	Main device (3.5-inch TFT)		CPAP, S, AutoS	25
G2S B25T	Main device (3.5-inch TFT)		CPAP, S, S/T, T	25
G2S B25VT	Main device (3.5-inch TFT)		CPAP, S, S/T, T	25
G2S LAB	Main device (3.5-inch TFT)		CPAP, S, AutoS, AutoCPAP, S/T, T	25
G2S B30T	Main device (3.5-inch TFT)		CPAP, S, S/T, T	30
G2S B30VT	Main device (3.5-inch TFT)		CPAP, S, S/T, T	30
G2S B30AT	Main device (3.5-inch TFT)		CPAP, S, AutoS, AutoCPAP, S/T, T	30

10. Package Contents

After unpacking the system, make sure you have everything shown here (Different models of the product may contain different components):

No.	Articles	Qty.	Notes
1	Device	1	
2	Air Filter	2	
3	Power Adapter	1	
4	Power Cord	1	
5	Mask	1	Optional
6	Cellular Module	1	Optional
7	SpO_2 Kit	1	Optional
8	SpO_2&GPRS Kit	1	Optional
9	Tube	1	Optional
10	SD Card	1	Optional
11	Carrying Case	1	Optional
12	Accompanying Documents	1	
13	Power Cord Locker	1	

All parts and accessories are not made of natural rubber latex.

The expected service life of the device is five years from first date of use, if the device is used, maintained, cleaned and disinfected in strict accordance with the User Manual.
The expected service life of the tube and mask refer to the user manual of the mask and L1 tube. The service cycle of the water chamber is 6 months, and the water chamber can withstand cleaning up to 180 times and disinfection 30 times. The WiFi kit and Cellular Module expected service life is one year from first date of use.
The mask and SpO_2 Probe are applied part.

WARNINGS!

• This device should only be used with the mask and accessories manufactured or recommended by the manufacturer. The use of inappropriate masks and accessories may affect the performance of the device and impair the effectiveness of treatment.
• Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this device could result in increased electromagnetic emissions or reduced electromagnetic immunity of this equipment and result in improper operation.
• Use of this device adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this device and the other equipment should be observed to verify that they are operating normally.
• When the insulation layer of the SpO_2 probe cable is damaged, do not connect the probe to the patient.
• Do not stack the long tubing or SpO_2 Kit lead near the patient's neck, as it may wrap around the patient's head or neck during sleep.
• Do not attach any equipment to the device unless recommended by the manufacturer or your health care provider.

• Please contact the manufacturer for an SD card if needed.
• The manufacturer cannot guarantee the normal function of the device, nor the safety and effectiveness of the device when it exceeding the expected service life.

IMPORTANT TIPS!

• If any of the above parts are missing, contact your home care provider.
• Contact your home care provider for additional information on the available accessories of this device. When using optional accessories, always follow the instructions enclosed with the accessories.
• The tube with a diameter of 15 mm or 22 mm is available.

11. System Features


Fig. 11-1

Name	Function
Humidifier Indicator	Indicate the humidity level. There are five levels in total. The number of indicator lights that light up is directly proportional to the humidity level. If the indicator lights are off, it means the humidifier is turned off
Mute Button	Press this button to mute the alert. However, if the problem causing the alert is not solved, the alert will sound again two minutes later
Knob	Start treatment and adjust device settings
Ramp Button	Enable the Ramp feature
Display Screen	Display operation menus, information, monitoring data, etc.
Power Indicator	Indicate the power supply status


Fig. 11-2

Name	Function
SD Card Slot	Insert the SD card into this slot
Air Outlet	Deliver pressurized air; connects to the tube
Communications Port	Connected to external equipment (only for connection of KS-CM01 SpO_2 Kit, WL-200 Cellular Module and SG-200 SpO_2&GPRS Kit to the device)
DC Inlet	An inlet for the DC power supply
Filter Cap (Air Inlet)	Place the cap on the air filter, which is used to filter dust and pollen in the air entering the device

CAUTION!

• The pictures in this manual are only for reference, if they are different from the material objects, the latter shall prevail.

12. First Time Setup

12.1 Placing the Device

Place the device on a firm, flat surface.

WARNINGS!

• If the device has been dropped or mishandled, if the enclosure is broken, or if water enters the enclosure, disconnect the power cord and discontinue use. Contact your home care provider immediately.

• If the room temperature is above 95°F (35°C), the airflow generated by the device may exceed 109.4°F (43°C). The room temperature must be kept below 95˚F (35°C) while the patient is using the device.

CAUTIONS!

• If the device has been exposed to very hot or very cold temperatures, allow it to acclimate to room temperature (approx. 2 hours) before beginning setup.

• Make sure the device is away from any heating or cooling equipment (e.g., forced air vents, radiators, air conditioners).

• The device is not suitable for use in high humidity environments. Make sure that no water enters the device.

• Make sure that bedding, curtains, or other objects (such as pests) are not blocking or entering the filter or vents of the device.

• Keep pets or children away from the device and avoid small objects being inhaled or swallowed.

• To avoid explosion, this device must not be used in the presence of flammable gases (e.g., anesthetics).

• Tobacco smoke may cause tar to build up in the device, which could lead to the malfunctioning of the device.

• Air must flow freely around the device to allow it to function properly.

• Do not move or tilt the device when there is water in the water chamber.

12.2 Installing the Air Filter and Filter Cap

(1) Attach the air filter to the filter cap, as shown in Fig. 12-1.


Fig. 12-1

(2) Install the filter cap containing the air filter to the device, as shown in Fig. 12-2.


Fig. 12-2

CAUTIONS!

- The air filter must be in place when the device is operating.
- When installing the air filter and filter cap, the device must be unplugged.

12.3 Connecting Power Supply

(1) Insert the plug of the power adapter into the DC Inlet on the back of the device;

(2) Connect the power cord to the power adapter;

(3) Plug the other end of the power cord into the power outlet.


Fig. 12-3

Note: The length of the power cord and power adapter is 1.5 m and 1.8 m respectively without the function of preventing electromagnetic interference.

WARNINGS!

• The device is powered on for use when the power cord and power adapter are connected. Use the **Knob** to turn the blower On / Off.

• Using the device at an AC voltage outside the specified range (see Section 6 "AC Power Consumption") may damage the device or cause device failure.

• Connect to the proper power source for proper operation of the device.

CAUTION!

• Inspect the power cord often for any signs of damage. Replace a damaged cord immediately.

IMPORTANT TIPS!

• After interruption and restoration of the power supply, the device will restore its pre-interruption working status automatically.

• To remove AC power, disconnect the power cord from the power outlet.

12.4 Connecting to Power Cord Locker

(1) Assemble the power cord locker to the device directed by the positioning groove.


Fig. 12-4

(2) Open the power cord locker, plug the power cord to the power supply, and press the locker downward to fix the power cord into the power port.


Fig. 12-5

The function of the locker is to prevent the power cord from falling off from the power port.

12.5 Assembling the Tube and Mask

(1) Connect one end of the tube to the air outlet of the main device, as shown in Fig. 12-6.


Fig. 12-6

(2) Connect the other end of the tube to the mask according to the user manual for the mask. Wear the mask.

WARNINGS!

• If multiple persons are going to use the device (e.g., rental devices), a low-resistance, main flow bacteria filter should be installed in-line between the device and the tube. Pressures must be verified by your home care provider when using spare or optional accessories.

• If you are using a mask with a built-in exhalation port, connect the mask's connector to the tube.

• If you are using a mask with a separate exhalation port, connect the exhalation port to the tube. Position the exhalation port so that the released air is blows away from your face. Connect the mask's connector to the exhalation port.

• If you are using a full-face mask (a mask covering both your mouth and nose), the mask must be equipped with a safety (entrainment) valve.

• To minimize the risk of CO_2 rebreathing, the patient should observe the following instructions:

- It is recommended to use the accompanying tube and mask provided by the manufacturer.
- Do not wear the mask for more than a few minutes while the device is not operating.
- Use only masks with vent holes. Do not block or try to seal the vent holes in the exhalation port.

CAUTION!

• When using the tube with a diameter of 15mm, set the tube of the Patient Menu to 15 mm (see Section 15.2 "Options of the Patient Menu and Corresponding Descriptions").

12.6 Using Oxygen with the Device

Oxygen may be added at the mask connection. Please observe the instructions listed below when using oxygen with the device.

WARNINGS!

• Connect the oxygen tube to the oxygen inlet of the mask.

• The oxygen supply must comply with the local regulations for medical oxygen.

• Turn on the device before turning on the oxygen. Turn off the oxygen supply before turning off the device. Explanation of Warning: When the device is turned off, but the oxygen flow remains, oxygen can accumulate inside the device's enclosure and pose a fire hazard. Turning off the oxygen supply before turning off the device will prevent oxygen accumulation in the device and reduce the risk of fire. This warning applies to most devices.

• Oxygen supports combustion. Keep the device and the oxygen source away from heat, open flames, any oily substances, or other sources of ignition. DO NOT smoke in the area near the device or the source.

• Sources of oxygen should be more than 1 m away from the device.

• When using oxygen with this system, a Pressure Valve must be placed in-line with the patient circuit between the device and the oxygen source. The pressure valve helps prevent the backflow of oxygen from the patient circuit into the device when the unit is off. Failure to use the pressure valve could result in a fire hazard.

• Do not connect the device to an unregulated or high pressure oxygen source. The pressure of oxygen source does not exceed the working pressure of the device.

12.7 Inserting the SD Card (Only for the device that equipped with SD card)

Insert the SD card into the SD Card Slot, as shown in Fig. 12-7.


Fig. 12-7

If the SD card is inserted correctly, a symbol indicating correct insertion will appear on the main screen of the device, as shown in Fig. 12-8.


Fig. 12-8

If the SD card is inserted incorrectly, a symbol indicating incorrect insertion will appear on the main screen of the device, as shown in Fig. 12-9.


Fig. 12-9

CAUTIONS!

• If no SD card is inserted, neither of the symbols will appear on the main screen of the device.

• To avoid data loss or any damage to the SD card, the SD card can only be removed after the main device stops delivering air.

12.8 Using the Optional Kits

Optional Kits contain KS-CM01 SpO_2 Kit, WL-200 Cellular Module and SG-200 SpO_2&GPRS Kit. For more details, please refer to the corresponding user manual.

CAUTIONS!

• The device can only use the SG-200 SpO_2&GPRS Kit, KS-CM01 SpO_2 Kit and WL-200 Cellular Module which configured for the device.

• When the power failure, no special operation is required after restore power supply. The device will continue to work according to the setting before.

12.9 Starting Treatment

Connect the device to a power outlet, press **the Knob** , and the device will start delivering air.

WARNINGS!

• Be sure to follow your home care provider's instructions on adjusting the settings! To order any accessories not included with this device, contact your home care provider.

• DO NOT connect any ancillary equipment to this device unless recommended by the manufacturer or your home care provider. If you suffer from chest discomfort, shortness of breath, stomach bloating, or severe headache when using the device, contract your home care provider or qualified medical personnel immediately.

13. Routine Use

13.1 Connecting the Tube

Connect the power cord, power adapter, and tube properly in accordance with the instructions in the First Time Setup. Connect the mask and headgear according to the user manual of the mask.

CAUTION!

• Before each use, examine the tube for any damage or foreign object. If necessary, clean the tube to remove the foreign object. Replace any damaged tube. Make sure that the mask does not leak.

13.2 Adjusting the Tube

Lie down on your bed, and adjust the tube so it is free to move if you turn over during sleep. Adjust the mask and headgear until you have a comfortable fit and until there are no airflow leakage around the mask.

13.3 Turning on the Airflow

Press **the Knob** to turn on the airflow. The screen will display treatment pressure and other information.

13.4 Heating the Water

Pay attention to the humidifier indicator lights when using the humidifier. The indicator lights indicate the On / Off state of the humidifier. It is off when all indicator lights go out.

CAUTION!

• Observe the water level in the water chamber before using the humidifier. Make sure there is sufficient water in the water chamber, and avoid heating the humidifier with an empty water chamber.

13.5 Using the Ramp Button

Every time **the Ramp Button** is pressed, the pressure will drop to the initial pressure, and then gradually rise to the prescribed treatment pressure according to the preset ramp time, so as to make it easy for the patient to fall asleep. The screen displays a real-time countdown of the remaining ramp time in minutes.

CAUTIONS!

• You can press **the Ramp Button** as often as you wish during sleep.

• The ramp feature is not prescribed for all users.

13.6 Turning the Device Off

Take off the mask and headgear, press and hold **the Knob** for two seconds, and the device will stop delivering air. Disconnect the power cord from the power outlet to turn off the device.

CAUTIONS!

• Do not position the device where it is difficult to operate the disconnection device.

• To isolate the device from the supply mains, disconnect the plug.

14. Heated Humidifier

The humidifier can be obtained from your home care provider. The humidifier can reduce nasal dryness and irritation by adding moisture (and heat if applicable) to the airflow.

14.1 Filling the Water Chamber

14.1.1 Removing the Water Chamber

Grab the water chamber, and pull it out of the device, as shown in Fig. 14-1.


Fig. 14-1

WARNING!

• Turn the device off and the heating plate and water to cool for approximately 15 minutes.

14.1.2 Filling the Water Chamber

Open the cap, as shown in Fig. 14-2, and fill the water chamber with approximately 360 ml of water, as shown in Fig. 14-3. Make sure that the water does not exceed the maximum water level line.


Fig. 14-2


Fig. 14-3

WARNING!

• Change water before every use and do not surpass the maximum water level line. Otherwise, the liquid may enter the tube.

CAUTIONS!

• Empty the water chamber when the heated humidifier is not in use.

• Distilled water is recommended.

14.1.3 Inserting the Water Chamber

Close the cap after it is filled with water, as shown in Fig. 14-4, and return it to the device, as shown in Fig. 14-5.


Fig. 14-4


Fig. 14-5

WARNING!

• For safety, the device must be placed on a flat surface below the height of the patient's head when he is lying on a bed, so that the condensation flows back to the water chamber rather than remaining in the tubing which can cause droplet spraying.

CAUTIONS!

• Avoid moving or tilting the device when the water chamber has water in it.

• Take precautions to protect furniture from water damage.

14.2 Emptying the Water Chamber

(1) **Removing the water chamber** according to instructions in 14.1.1.

(2) **Emptying the water chamber:** Open the cap, as shown in Fig. 14-6, and pour any remaining water out of the water chamber.


Fig. 14-6

CAUTION!

• Empty and air-dry the water chamber when the device is not in use.

(3) **Inserting the Water Chamber** according to instructions in 14.1.3.

14.3 Setting the Humidity Level

After the device is powered on, turn **the Knob** to turn the heated humidifier on or off and to adjust the humidity level according to instructions in the screen of the device.

There are five humidity levels available, and the number of blue indicator lights that light up is directly proportional to the humidity level. If none of the indicator lights light up, it means that the heated humidifier is turned off.

The temperature of the water in the water chamber maintains a constant set level. Three indicator lights light up when the humidity is adjusted to Level 3, as shown in Fig. 14-7.


Fig. 14-7

WARNING!

• Do not touch the heating plate of the device when it is in operation, otherwise you may get burned. Stop heating when the heated humidifier is not in use.

CAUTIONS!

• Generally speaking, the humidity level inside the mask is low when the water temperature is low.

• The greater the difference between the temperature inside the air tubing and the room temperature, the more likely condensation will occur in the tubing.
• If there is only a small amout of condensed water droplets in the tubing in the morning after treatment, the humidity level is appropriate; if there is a large amount of condensed water droplets inside the tubing and / or the mask, the humidity level is too high and should be set lower; Nasal dryness means that the humidity level is too low and should be set higher.

15. Navigating the Patient Menu

15.1 Steps to Navigate the Patient Menu

15.1.1 Accessing the Main Interface

Connect the power cord and power adapter properly. The screen displays the Main Interface shown in Fig. 15-1, or the Main Interface shown in Fig. 15-2, or the Main Interface shown in Fig. 15-3, or the Main Interface shown in Fig. 15-4.

Fig. 15-1

Note: The above interface only applies to G2S B20S and G2S B20A.

Fig. 15-2

Note: The above interface only applies to G2S B20T.

Fig. 15-3

Note: The above interface only applies to G2S B25A and G2S B25S.


Fig. 15-4

Note: The above interface only applies to G2S B25T, G2S B25VT, G2S B30T, G2S B30VT, G2S B30AT and G2S LAB.

15.1.2 Bringing up the Initial Setup Interface

From the Main Interface shown in Fig. 15-1 or Fig. 15-2 or Fig. 15-3 or Fig. 15-4, press and hold **the Ramp Button** for three seconds. The screen displays the Initial Setup Interface of the Patient Menu, as shown in Fig. 15-5.


Fig. 15-5

The first icon on the left side of the screen indicates the Main Interface, the second icon indicates the Initial Setup Interface, and the third icon indicates the iCode Interface. As you turn **the Knob**, the cursor will switch among the three icons, and the interface displayed on the screen changes accordingly.

15.1.3 Accessing the Setup Interface

When the cursor is on the icon, the screen displays the Setup Interface. Access the Setup Interface by pressing **the Knob**. The first option on the Setup Interface is then displayed in blue, as shown in Fig. 15-6.


Fig. 15-6

15.1.4 Selecting Options

As you turn **the Knob** clockwise, the cursor moves downwards from one option to another. As you turn it counterclockwise, the cursor moves upwards. When the cursor is on a

certain option, press **the Knob**, and the option is then displayed in yellow, meaning that the option is now ajustable, as shown by the **Humidifier** option in Fig. 15-7.


Fig. 15-7

15.1.5 Adjusting Options

Adjust the option by turning **the Knob**. As shown in Fig. 15-7, the **Humidifier** option is selected. As you turn **the Knob** clockwise, the number increases, indicating a higher humidity level. As you turn **the Knob** counterclockwise, the number decreases, indicating a lower humidity level. The **Humidifier** option is still displayed in yellow, as shown in Fig. 15-8.


Fig. 15-8

15.1.6 Confirming Adjustments

Press **the Knob** to confirm your adjustment for a particular option. The option is then displayed in blue, as shown in Fig. 15-9.


Fig. 15-9

15.1.7 Turning Pages

When the cursor is on **Mask Type**, the last option shown in Fig. 15-9, the remaining options will appear on a new page if you continue to turn **the Knob** clockwise, as shown in Fig. 15-10.


Fig. 15-10

Note: are page turning symbols.

15.1.8 Exiting the Patient Menu

(1) Returning to the Initial Setup Interface

Move the cursor to the **Back** option by turning **the Knob**, as shown in Fig. 15-11.


Fig. 15-11

Press **the Knob** to make the cursor jump to the second icon on the left side of the screen. The screen displays the Initial Setup Interface, as shown in Fig. 15-12.


Fig. 15-12

(2) Returning to the Main Interface

Turn **the Knob** to move the cursor to the **Home** option, as shown in Fig. 15-13.


Fig. 15-13

Press **the Knob** to exit the Patient Menu. The screen will display the Main Interface shown in Fig. 15-1 or Fig. 15-2 or Fig. 15-3 or Fig. 15-4.

15.2 Options in the Patient Menu and Corresponding Descriptions

Option	Range	Description
Humidifier	Off / 1 ～ 5	There are five humidity levels available. As the number increases, the humidity rises accordingly. "Off" means the humidifier is turned off
Tube	15 mm / 22 mm	Set diameter of the tube by adjusting this option
Reslex	Off / 1 ～ 3	This feature enables the device to automatically reduce the treatment pressure when the patient exhales, so as to make him/her more comfortable. The higher the number, the more pressure the device reduces. "Off" means this feature is disabled
Ramp Time	0 ～ Max Ramp	In order to increase comfort and help the patient fall asleep easily, the pressure can be increased gradually, when the Ramp feature is enabled. The ramp time during which the initial pressure rises to the preset treatment pressure can be adjusted. As you turn **the Knob** to the nearest point, the number increases or decreases by five minutes. The screen displays a real-time countdown of the remaining ramp time in minutes
Delay	On / Off	When the humidifier is on, this feature allows the airflow to continue for about 15 minutes at a low pressure (about 2 hPa) after you press **the Knob** to discontinue the treatment. In this process, the vapor left in the humidifier will be blown away to avoid any damage to the device. When this feature is set to "Off," which means it is disabled, the device will stop delivering air instantly after you press **the Knob**
Time	00:00 — 23:59	Set time by adjusting this option
Brightness	High / Low	Set the brightness of the screen by adjusting this option
Mask Type	Full Face / Nasal / Pillow / Other	There are three mask types available, Full Face (full-face mask), Nasal (nasal mask), and Pillow (nasal pillow mask). The patient can choose other suitable masks as well. When selecting masks other than the above three types of masks, the patient can set the mask type as Other
iCode	iCode / iCode QR / iCode QR+	iCode provides access to the patient's compliance data during a recent time period. The iCode mode displays data in number sequence, and the iCode QR / iCode QR + mode displays data in two-dimensional codes
Use Time	0 ～ 50000 h	Displays how long has the device been used by the patient. The use time can be erased
About	Model / SN / Firmware version / ID / PIN	Show the relevant information about the device, this is only for users to view, can not be modified. Model: the device model; SN: Serial Number of the device; Firmware version: Software version of the device; ID: Contains information such as gallery and language; PIN: Personal identification code

16. Alarm

This chapter describes device the alarms of the device and proper responses the operators should make to different alarms.

After the device is running, disconnect its power supply by unplugging the power cord, an audible alert will sound like "beep beep beep, beep-beep, beep beep beep, beep-beep", which means that the alarming system of the device works normally.

WARNING!

• To avoid the failure of the alarm, do not set the alarm limit to the limit value.

16.1 Alarm Grades and Description

The alarm grades and description of this equipment are shown as follows:

Grade	Sign of Grading	Description
High	!!!	Requires operator to make instant response
Intermediate	!!	Requires operator to make instant on-time response
Low	!	Requires operator to be more cautious about the change of the state of equipment

16.2 Visual Alarm

The visual alarm grades are indicated by the background of the alarming information on the top of the screen and the color of the LED light under the Mute Button, which is shown as follows:

Grade	Visual	Description
High	Red	Red light flickers—high-grade alarming
Intermediate	Yellow	Yellow light flickers—intermediate alarming
Low	Yellow	Yellow light indicates in a fixed manner—low-grade alarming

16.3 Auditory Alarm

The alarming sounds at different grades will occur and are described as follows:

Grade	Auditory	Description
High	●●● ●● ●●● ●●	beep beep beep beep-beep beep beep beep beep-beep
Intermediate	● ● ●	beep beep beep
Low	●	beep

In accordance with the requirements of relevant standards, the volume of the audible alarm signal meets the requirements. The sound pressure range of the measured auditory alarm signal is described as follows:

Alarm Condition	Measured sound pressure level (dB)	A-weighted sound pressure level averaged over the measurement surface (dB)	Remarks
High priority	52.2	38.5	Maximum volume
Median priority	51.8	39.6	Maximum volume
Low priority	51.8	37.2	Maximum volume

16.4 Alarm Silence

When the device sounds an alarm, press the alarming Mute Button  and it will become silent for 100 to 120 seconds. Then the alarm will sound again immediately at the end of the silence. If the Mute Button is re-pressed during the silence period, the silence period will be re-timed.

16.5 Alarm Messages and Description

Alarm Message	Alarm Priority	Alarm Type	Description
Power Failure!!!	High Priority	Technology Alarm	An audible alert will sound within 6 s if the device is accidentally disconnected from power supply when it is delivering air. Alarming duration time is no less than 30 s. **Note:** (1) The alert will not sound if power failure occurs when the device is in standby state. (2) No alert message will appear on the screen during a power failure.
Device fault!!!	High Priority	Technology Alarm	An audible alert will sound if no airflow comes out of the machine; the screen will display "**Device fault!!!**".
Tube disconnected!!! (only applies to G2S B20T, G2S B25T, G2S B25VT, G2S B30T, G2S B30VT, G2S B30AT and G2S LAB)	High Priority	Function Alarm	When the airflow is on, an audible alert will sound if the tube accidentally detached, the screen will display "**Tube disconnected!!!**".

High Pressure!!!	High Priority	Function Alarm	When the airflow is on, an audible alert will sound if the airway pressure exceeds the alarm limit; the screen will display "**High Pressure!!!**". **Note:** The thresholds for different models: Off, 5 ～ 21 hPa applies to G2S B20T, in 0.5 hPa increments, the default setting is "**20 hPa**". Off, 5 ～ 26 hPa applies to G2S B25T, G2S B25VT and G2S LAB, in 0.5 hPa increments, the default setting is "**25 hPa**". Off, 5 ～ 31 hPa applies to G2S B30T, G2S B30VT and G2S B30AT, in 0.5 hPa increments, the default setting is "**30 hPa**".
Low Pressure!!	Middle Priority	Function Alarm	When the airflow is on, an audible alert will sound if the airway pressure is below the alarm limit; the screen will display "**Low Pressure!!**". **Note:** The limens for different models: Off, 3 ～ 19 hPa applies to G2S B20T, in 0.5 hPa increments, the default setting is "**4 hPa**". Off, 3 ～ 24 hPa applies to G2S B25T, G2S B25VT and G2S LAB, in 0.5 hPa increments, the default setting is "**4 hPa**". Off, 3 ～ 29 hPa applies to G2S B30T, G2S B30VT and G2S B30AT, in 0.5 hPa increments, the default setting is "**4 hPa**".
Low RR!!! (only applies to G2S B20T, G2S B25T, G2S B25VT, G2S B30T, G2S B30VT, G2S B30AT and G2S LAB)	High Priority	Function Alarm	When the airflow is on, an audible alert will sound if the respiratory rate is below the alarm limit; the screen will display "**Low RR!!!**". Setting range: Off, 4 ～ 40 BPM, in 1 BPM increments, the default setting is "**6 BPM**". **Note:** This function is available under the work mode of S/T or T.
Leak!!	Middle Priority	Function Alarm	When the airflow is on, an audible alert will sound if the air leak rate exceeds 150 L/min; the screen will display "**Leak!!**". The alarming duration time is no less than 30 s.

Mask Blocked!! (only applies to G2S B20T, G2S B25T, G2S B25VT, G2S B30T, G2S B30VT, G2S B30AT and G2S LAB)	Middle Priority	Function Alarm	When the airflow is on, an audible alert will sound if the vents of the mask are blocked; the screen will display "**Mask Blocked!!**".
Low MV!! (only applies to G2S B20T, G2S B25T, G2S B25VT, G2S B30T, G2S B30VT, G2S B30AT and G2S LAB)	Middle Priority	Function Alarm	When the airflow is on, an audible alert will sound if the minute volume is below the alarm limit; the screen will display "**Low MV!!**". Setting range: Off, 1 ～ 30 L/min, in 1 L/min increments, the default setting is "**1 L/min**".
Low Input Voltage!!	Middle Priority	Technology Alarm	If the voltage supplied by power adaptor is lower than 22 V, an audible alert will sound and the screen will display "**Low Input Voltage!!**".
High RR!! (only applies to G2S B20T, G2S B25T, G2S B25VT, G2S B30T, G2S B30VT, G2S B30AT and G2S LAB)	Middle Priority	Function Alarm	When the airflow is on, an audible alert will sound if the respiratory rate exceeds the alarm limit; the screen will display "**High RR!!**". Setting range: Off, the setting value of Low RR ～ 80 BPM, in 1 BPM increments, the default setting is "**40 BPM**". **Note:** This function is avaliable under the work mode of S/T or T.
Humidifier Failure!!	Middle Priority	Function Alarm	When humidifier is applied, an audible alert will sound when the humidifier fails to work in 10 minutes; the screen will display "**Humidifier Failure!!**".
Please Change Filter!	Low Priority	Technology Alarm	When the Filter Alert feature is enabled, an audible alert will sound if the preset replacement time is reached but the air filter is not replaced; the screen will display "**Please Change Filter!**". The default setting is "**Off**".
SD Card Full!	Low Priority	Technology Alarm	The screen will display "**SD Card Full!**" if the SD card has reached its maximum capacity.
Reinsert SD card!	Low Priority	Technology Alarm	The screen will display "**Reinsert SD card!**" if the SD card fails to work.
SpO_2 faults!	Low Priority	Technology Alarm	Disconnect the SpO_2 Probe from the SpO_2 Probe equipment and place it in series with a circuit with which each SpO_2 Probe wire can be opened or shorted to any other SpO_2 Probe wire. The screen will display **"SpO_2 faults!"**.

Note: the delay time of alarming system of this device is no more than 1 s.

16.6 Reposition of Alarming

After the alarm faults are cleared, the residual alarm messages still remain (alarming messages are shown on the top of the screen without any visual and auditory alarm). Turn the button leftwards or rightwards to reduce the residual alarming messages.

16.7 Alarm Journal

The alarm journal is designed to record the latest 6 alarm messages. Retained inside the device, the alarm journal will not be lost after a power interruption and the latest alarm message will replace the previous one, retaining 6 messages.

WARNINGS!

• Before using the device, the oeprator should check whether the current alarm preset value is applicable to each patient. Such preset value can only be changed by a medical professional and cannot be modified by the patients at home.

• In the case of a power failure or power cut for no more than 30 seconds, the last set alarm value will be restored at the next oepration.

17. Cleaning and Disinfection

WARNINGS!

• Regular cleaning of the device and its accessories is very important for the prevention of respiratory infections.

• To avoid electric shock, always unplug the device before cleaning.

• Use mild soap that is nontoxic to humans.

• Follow the manufacturer's instructions on cleaning the mask and tube and on determining the frequency of cleaning.

• Before cleaning, check that the device is disconnected from the power supply, whether the power cord is unplugged, and whether the water chamber of the humidifier has cooled down. Make sure that the heating plate has cooled down to room temperature, so that you do not get burned.

• Do not open or modify the device. There are no user serviceable parts inside. Repairs and service should only be performed by an authorized service agent.

• The device shall not be serviced or maintained while a patient is using it.

CAUTIONS!

• Overheating of the materials could lead to early fatigue of the materials.

• Do not use solutions containing chlorinated lime, chlorine, or aromatic to clean the device and its accessories. Liquid soap containing moisturizing agents or antimicrobials should not be used either. These solutions may harden cleaned materials or reduce their lifespan.

• Do not clean with water above 80°C (176°F) or dry the device and its accessories in an environment with temperature above 80°C (176°F). High temperatures could reduce product life.

• Do not immerse the device in any fluids.

17.1 Cleaning the Mask and Headgear

For details, refer to the cleaning instructions in the user manual of the mask.

17.2 Cleaning the Optional Kits

For details, refer to the cleaning instructions in the user manual of the corresponding user manual.

CAUTIONS!

• It is recommended to clean the SpO_2&GPRS Kit / SpO_2 Kit / Cellular Module once a week.

17.3 Cleaning the Water Chamber

Preparation tools: soft bristle brush, drinking quality water, mild liquid detergent

Remove the Water Chamber according to instructions in 14.1.1.

(1) **Opening the Water Chamber:** Open the cap of the water chamber, as shown in Fig. 17-1.


Fig. 17-1

(2) **Cleaning the Water Chamber:** Clean the water chamber with a soft bristled brush for one minute after soaking the water chamber in detergent solution (see table below) for 5 minutes. Pay particular attention to all crevices and cavities. Then rinse with running water for 5 minutes. Wipe it dry with a soft cloth or air dry out of direct sunlight.

Detergent	Water temperature
Alconox (diluted at 1%)	Warm water (approx. 113 to 140°F or 45 to 60°C)

WARNINGS!

• Emptying and cleaning the water chamber daily will help prevent mold and bacteria growth.

• Allow the water in the chamber to cool down to room temperature before removing it from the device.

CAUTIONS!

• Clean the water chamber only after the water in it cools. Make sure that no water enters the device.

• After cleaning, rinse the water chamber throughly in clean water to make sure that no soap

residue is left; then wipe it dry with a lint-free cloth, so as to prevent calcareous accumulations.

• Check the water chamber for any leak or damage. Replace the water chamber if there is any damage.

• It is recommended to do daily cleaning of the water chamber.

17.4 Cleaning the Tube

For details, refer to the cleaning instructions in the user manual of the L1 tube.

CAUTIONS!

• It is recommended to clean the tube once a week.

17.5 Replacing the Air Filter

(1) Open the air filter cap to remove the air filter.

(2) Put a new air filter in the filter area, and then place the filter cap back properly.

CAUTIONS!

• To avoid material damage, do not place the spare air filter in direct sunlight, humid environments, or temperatures below the freezing point. The air filter should be replaced every 6 months (It may be replaced more frequently based on actual sanitary conditions).

• Operating the device with a dirty air filter may stop it from working properly and may cause damage to the device.

• Replacing the air filter and filter cap, device must be unplugged.

17.6 Cleaning the SpO_2 Probe

For details, refer to the cleaning instructions in the user manual of the SpO_2 Probe.

17.7 Disinfection

Generally speaking, if you have strictly followed the above cleaning instructions, you do not have to disinfect the device and / or humidifier. If the device is contaminated or used in clinical trials, you may purchase disinfectants from a home medical equipment company to disinfect the device.

Disinfection of the Water Chamber:

Prior to disinfection, clean the water chamber according to Section 17.3 "Cleaning the Water Chamber". In the procedures below, only one disinfection process needs to be performed each time.

Preparation before disinfection

Preparation tools: soft bristle brush, drinking quality water, applicable disinfectant

Chemical Disinfection:

(1) Clean the water chamber according to steps in cleaning instructions.

(2) Put the Johnson's CIDEX OPA 0.55% into a plastic box, which can completely immerse the water chamber;

(3) Immerse the water chamber in the CIDEX OPA 0.55% for 12 min at 20°C.
(4) Rinse the water chamber with 8L running water for 3 times to remove the residual disinfectant.
(5) Wipe it dry with a soft cloth or air dry out of direct sunlight.

UHT Disinfection (90°C water)

(1) Open the cap of the water chamber
(2) Put the water chamber of the device in the sink of the Automatic Washer Disinfector.
(3) Immerse the water chamber in the water chamber.
(4) Heat the water to 90°C and immerse the water chamber for 5 min.
(5) Wipe it dry with a soft cloth or air dry out of direct sunlight.
Note: Use an EN ISO 15883-compliant AWD.

17.7.1 Disinfecting the mask and headband

For details, refer to the disinfection instructions in the user manual of the mask.

17.7.2 Disinfection of the SpO_2 Probe

For details, refer to the cleaning instructions in the user manual of the SpO_2 Probe.

WARNINGS!

• After disinfection, rinse any disinfected component in clean water thoroughly, especially components in close contact with the patient such as the mask, headgear, and tube, so as to prevent disinfectant residuals from damaging the skin or respiratory tract or causing allergies.
• Disinfection of this device and its components other than recommended is not permitted, otherwise the manufacturer will not be able to verify the safety or performance of the device.
• In order to prevent cross-infection of patients or contamination of equipment, BSF (Breathing System Filter) that meet ISO 23328-1:2003 and ISO 23328-2:2002 standards and have medical device registration certificates can be used.
(1) BSF need to be replaced before different patients using this device.
(2) When using the BSF, please install and operate it according to the instructions of the BSF, and pay attention to adjust the output pressure setting of the device according to resistance of the BSF to ensure delivering normal treatment pressure.
(3) Atomization or humidification will increase the resistance of the BSF. The operator must often monitor the resistance increase and blockage of the BSF to ensure delivering normal treatment pressure.
• If you use ozone or other cleaning and disinfection methods that are not recommended by BMC, BMC will not be able to verify the safety or performance of the equipment.

CAUTIONS!

• Disinfectants tend to damage the materials and reduce the life of components. Try to select the appropriate disinfectant, and follow the disinfectant manufacturer's instructions and recommendations.

• After disinfection, check the disinfected component for any signs of damage. Replace any damaged component immediately.

18. Traveling with the Device

(1) Use the BMC carrying case to carry the device and accessories along with you. Do not put them in your checked baggage.

(2) This device operates on a power supply of 100 V - 240 V and 50 Hz / 60 Hz, and is suitable for use in any country in the world. No special adjustment is required, but you will need to find out the types of the power sockets at your destination. If necessary, bring a power socket adaptor which can be purchased at electronics stores.

(3) Remember to bring a spare air filter and emergency documentation (filled and signed by your home care provider) about this device. If you plan to travel by air, remember to bring the multilingual emergency documentation about respiratory therapy, in case that the border and customs officers in your destination country inspect the device. With emergency documentation, you can prove to them that it is a medical device.

(4) Security Stations: For convenience at security stations, there is a note on the bottom of the device stating that it is medical equipment. It may be helpful to bring this manual with you to help security personnel understand the device.

CAUTIONS!

• Empty the water chamber of the humidifier before packing the device for your trip; to prevent any remaining water from entering the device.

• Using the device at an incorrect altitude setting may result in higher airflow pressures than the specified setting. Always verify the altitude setting when traveling or relocating.

• If the device is used when the atmospheric pressure is outside the specified range (See Section 6), the accuracy of the leak alert will be affected.

19. Transferring the Device to another Patient

If the device is transferred to another patient, components in close contact with the previous user, including the mask, headgear, tube, and air filter, should be replaced to prevent cross-infection.

20. Reordering

Contact your home care provider to order accessories or replacement filters.

The device does not require routine service.

WARNINGS!

• If you notice any unexplained changes in the performance of the device, if it is making unusual or harsh sounds, if it has been dropped or mishandled, if the enclosure is broken, or if water has entered the enclosure, please stop using the device and. Contact your home care provider.

• If the device fails to work properly, contact your home care provider immediately. Never attempt to open the enclosure of the device. Repairs and adjustments must be performed by the manufacturer-authorized service personnel only. Unauthorized service could cause injury, invalidate the warranty, or result in costly damage.

• If necessary, contact your local authorized dealer or BMC Medical Co., Ltd. for technical support and documents.

21. Technical Support

Please contact BMC directly if you need the circuit diagram of the device or the list of components for certain purposes such as maintenance or connection to other equipment. BMC will provide the technical documents in whole or in part according to your needs.

22. Disposal

This product complies with the requirements of Directives 2011/65/EU and (EU) 2015/863 (RoHs) on the restriction of hazardous substances.

According to the European Directive 2012/19/EU (WEEE) on electrical and electronic equipment, products cannot be discarded indiscriminately and must be disposed of in accordance with the laws and regulations of the country in which disposal occurs.

The crossed-out wheeled bin symbol  indicates that the product bearing this symbol may not be disposed of together with general household waste, but instead requires separate disposal. This requirement for separate disposal is based on the European Directive 2012/19/EU for electrical and electronic equipment. You can hand in the product at a municipal collection point, for example. This reduces the impact on natural resources and prevents contamination of the environment through the release of hazardous substances.

For further information regarding product disposal, please contact your specialist distributor.

23. Troubleshooting

The table below lists common problems you may encounter with the device and possible solutions to resolve them. If none of the corrective actions solve the problem, please contact your home care provider.

WARNING!

• Due to sensor performance degradation or loose equipment performance degradation, resulting in abnormal use or other problems, please stop using and contact the manufacturer.

23.1 Common Problems in Patients and Corresponding Solutions

Problem	Possible Cause	Solution(s)
Dry, cold, runny, and blocked nose; having a cold	The nose reacts to the airflow and cold. Due to fast airflow, the air becomes cold, resulting in the irritation of nasal mucosa and subsequent dryness and swelling	Increase the humidity setting of the humidifier. Contact your physician, and continue treatment unless the physician suggests the opposite
Dry mouth and throat	It may be because the patient sleeps with the mouth open, and the pressurized air flows out through the mouth, causing dryness of the nasal passage and throat	Use a chin strap to prevent the mouth from opening during sleep, or use a full-face mask. Contact your physician for details
Eye irritation	The mask may not be the correct size or type, or the mask may be incorrectly positioned resulting in an air leak	Narrow the distance between the forehead support of the mask and the forehead. Note that adjusting the mask too tight may leave marks on the patient's face. Add additional filling to the mask so it does not leak. Contact your home care provider for an appropriate mask. Add additional filling to the mask if necessary
	Mask cushion (the soft part of the mask) hardens	Replace the mask or mask cushion
Facial reddening	The mask is too tight	Loosen the headgear
	The distance between the forehead support of the mask and the forehead is not correct	Try a different distance. The angle and size of the forehead support differ according to the type of masks
	Wrong mask size	Contract your home care provider for a correct-size mask

Problem	Possible Cause	Solution(s)
Facial reddening	The patient is allergic to the materials of the mask	Contact your physician and home care provider. Use a latex-free mask. Place a lining between the skin and mask
Water in mask	When the humidifier is used, the humidified air tends to condense in the cold tube and mask if the room temperature is low	Turn the humidity setting down, or raise the room temperature. Place the tube under the quilt, or use the tube cover. Hang the tube loosely, and make sure that the lowest part of the tube should be lower than the patient's head
Nasal, sinus, or ear pain	Sinus or middle ear inflammation	Contact your physician immediately
Discomfort due to inability to adapt to the treatment pressure	The patient will feel uncomfortable when the treatment pressure is higher than 13 hPa. However, the treatment pressure is determined by the patient's conditions, and the device will not be able to treat sleep apnea if the treatment pressure is set too low	It takes a maximum of four weeks for the patient to adapt to pressurized air. Relax and breathe through the nose. If the problem still exists, contact your physician
Obstructive sleep apnea symptoms recur	It may be because the patient sleeps with the mouth open, and the pressurized air flows out through the mouth, causing a blockage in the respiratory tract	Use a chin strap to prevent the mouth from opening during sleep, or use a full-face mask. Contact your physician for details
The device is too noisy	The tube is not connected properly	Reconnect the tube properly
Air delivered from the device is abnormally hot	The air inlet of the device may be partially blocked, leading to insufficient airflow into the device	Replace the air filter (see 17.5 Replacing the Air Filter), and clean the air inlet Place the device in an area where air flows freely, and make sure that the device is at least 20 centimeters away from the wall, curtain, or other things

23.2 Common Problems in the Device and Corresponding Solutions

Problem	Possible Cause	Solution (s)
The device does not work when it is turned on	The Auto On / Off feature is enabled	Take a few deep breaths with the mask on, and the device will start automatically
	Power is not connected properly	Ensure that the power cord, power adapter, and the device are connected properly
	There is no voltage	Check whether a power outage occurs by turning on a light or other means. If you are sure the fuse in the device is broken, contact your home care provider for repair
	Cannot find any cause	Contact your home care provider
The device is working, but the pressure inside the mask differs from the set treatment pressure	The tube is not connected properly	Reconnect the tube properly
	There may be holes in the mask or pressure sensing tube	Contact your home care provider
	It is a faulty device	Contact your home care provider
The device produces very low pressures	The air inlet of the device may be blocked	Replace the air filter (see 17.5 Replacing the Air Filter), and clean the air inlet. Make sure the air inlet is unblocked
	The treatment pressure has been changed accidentally	Contact your physician
	When the Ramp feature is enabled, it takes some time for the initial pressure to rise to the treatment pressure. This is normal	If necessary, disable the Ramp feature, or set the ramp time shorter
After the device is turned on, the screen displays intermittently, or displays nothing at all	The operating system of the device needs to be readjusted or restarted	Unplug the power cord of the device, and re-plug it 20 seconds later
The device is in standby, and will not start	The operating system of the device needs to be readjusted or restarted	Unplug the power cord of the device, and re-plug it 20 seconds later

24. EMC Requirements

The cables must be provided by BMC. The information of each cable is as follows:
Power Cord: 1800 mm ± 30 mm, unshielded.

Guidance and manufacturer's declaration - electromagnetic emissions	
The device is intended for use in the electromagnetic environment specified below. The user of the device should ensure that it is used in such an environment.	
Emissions Test	**Compliance**
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies

WARNINGS!

• This device may be interfered with by other equipment, even if the equipment complies with CISPR EMISSION requirements.

• Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the G2S BPAP series, including cables specified by the manufacturer. Otherwise, the performance degradation of this equipment could be resulted.

• During the operation of the device, due to electrostatic interference, the following phenomena may occur: (1) Temporary loss of function or performance degradation, such as abnormal screen display, etc. The device will return to normal after being restarted; (2) Automatic restart of the device. These phenomena will not affect the normal use of the device, nor will they cause permanent performance degradation or loss of function of the device.

• This device is not intended for use with MRI equipment and is not allowed to be used in MR environments.

Guidance and manufacturer's declaration - electromagnetic immunity		
The device is intended for use in the electromagnetic environment specified below. The user of the device should make sure that it is used in such an environment.		
Immunity Test	**IEC 60601 Test Level**	**Compliance Level**
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air
Electrical fast transient / burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input / output lines	±2 kV for power supply lines ±1 kV for input / output lines
Surge IEC 61000-4-5	±1 kV line (s) to line (s) ±2 kV common mode	±1 kV line (s) to line (s) ±2 kV common mode
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U_T; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_T; 1 cycle 70% U_T; 25 / 30 cycle At 0° 0% U_T; 250 / 300 cycle	0% U_T; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_T; 1 cycle 70% U_T; 25 / 30 cycle At 0° 0% U_T; 250 / 300 cycle
Power frequency (50 Hz / 60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m
Note: U_T is the AC mains voltage prior to application of the test level.		

Guidance and manufacturer's declaration - electromagnetic immunity		
The device is intended for use in the electromagnetic environment specified below. The user of the device should make sure that it is used in such an environment.		
Immunity Test	**IEC 60601 Test Level**	**Compliance Level**
Conducted RF IEC 61000-4-6	3 V 0.15 MHz ～ 80 MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz	3 V 0.15 MHz ～ 80 MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	10 V/m 80 MHz to 2.7 GHz
Note 1: At 80 MHz and 800 MHz, the higher frequency range applied. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		
[a] Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device. [b] Over the frequency range 150 kHz to 80 MHz, the field strengths should be less than 10 V/m.		

Recommended separation distances between RF wireless communications equipment

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between RF wireless communications equipment and the device as recommended below, according to the maximum output power of the communications equipment.

Frequency (MHz)	Maximum Power (W)	Distance	IEC 60601 Test Level	Compliance Level
385	1.8	0.3	27	27
450	2	0.3	28	28
710	0.2	0.3	9	9
745				
780				
810	2	0.3	28	28
870				
930				
1720	2	0.3	28	28
1845				
1970				
2450	2	0.3	28	28
5240	0.2	0.3	9	9
5500				
5785				

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Test frequency	Modulation	Immunity Test Level (A/m)
30 kHz [a]	Pulse modulation [b] CW	8
134.2 kHz	Pulse modulation [b] 2.1 kHz	65 [c]
13.56 MHz	50 kHz	7.5 [c]

[a] This test is applicable only to ME EQUIPMENT and ME SYSTEMS intended for use in the HOME HEALTHCARE ENVIRONMENT.
[b] This carrier shall be modulated using a 50% duty cycle square wave signal.
[c] r.m.s., before modulation is applied.

25. Limited Warranty

BMC Medical Co., Ltd. warrants that the device shall be free from defects of workmanship and materials and will perform in accordance with the product specifications for a period of one (1) year for main device and three (3) months for all accessories from the date of sale by BMC Medical Co., Ltd. to the dealer. If the product fails to perform in accordance with the product specifications, BMC Medical Co., Ltd. will repair or replace, at its option, the defective material or part. BMC Medical Co., Ltd. will pay customary freight charges from BMC Medical Co., Ltd. to the dealer location only. This warranty does not cover damage caused by accident, misuse, abuse, alteration and other defects not related to material or workmanship.

BMC MEDICAL CO., LTD. DISCLAIMS ALL LIABILITY FOR ECONOMIC LOSS, LOSS OF PROFITS, OVERHEAD OR CONSEQUENTIAL DAMAGES WHICH MAY BE CLAIMED TO ARISE FROM ANY SALE OR USE OF THIS PRODUCT. SOME STATES DO NOT ALLOW THE EXCLUSION OR LIMITATION OF INCIDENTAL OR CONSEQUENTIAL DAMAGES, SO THE ABOVE LIMITATION OR EXCLUSION MAY NOT APPLY TO YOU.

To exercise the rights under this warranty, contact the local authorized dealers or:

MANUFACTURER:
BMC Medical Co., Ltd.
Room 110 Tower A Fengyu Building, No. 115 Fucheng Road, Haidian, 100036 Beijing, PEOPLE'S REPUBLIC OF China
Tel: +86-10-51663880
Fax: +86-10-51663880 Ext. 810
URL: en.bmc-medical.com
E-mail: intl@bmc-medical.com

EU AUTHORISED REPRESENTATIVE:
Shanghai International Holding Corp. GmbH (Europe)
Eiffestraβe 80, 20537 Hamburg, Germany
Tel: 0049-40-2513175
Fax: 0049-40-255726

V3.01
Issue date: February 7, 2023