Revision date: 19/12/2021

Revision number: 04/ Dec 2021

Disposable Bedsheet (Single & Multi Layers)

Instructions For Use

Product code

- MDC729421321050005 (1- Layer)
- MDC729421321050001 (3- Layer)
- MDC729421321050003 (5- Layer)
- MDC729421321050002 (7- Layer)

Product description

- Single and multiple layers that fits all hospital beds and is Waterproof & Disposable
- Eliminates waterproof mattress pad usage
- The product is latex, lead & phthalate free
- Consists of proprietary nonwoven blend of polyester, viscose and bamboo
- This product is packed and sealed individually

Category

Patient bedding

Intended use

For incontinence, bed-ridden patients or any person where a disposable bed sheet can be used

Direction of use



1. Put on the fitted sheet to the mattress. The head of the bedsheet is indicated with InnoQ Tag



2. To peel the top layer, pull down from the top corner and peel along the seams to reveal a new sheet underneath



3. Dispose the used layer responsibly

Material compositions

- 55% Polyester fibers
- 40% Rayon viscose fibers
- 5% Sustainable Bamboo Fibers with a Polyester spun bond skirt and latex free elastic

Features

- 32% softer than traditional bed sheets
- 100% waterproof
- Bacterial will not transmit between layers
- Microporous laminate allows air to pass through
- Latex free

Testing

Independent laboratory tested for flammability, CPSC conformity, absence of heavy metals, absence of phthalates and formaldehyde

Shelf life

5 years

Storage condition

- Dry storage to up 70°C
- Do not refrigerate

Legal manufacturer

Global SCP Medical Systems Pte. Ltd., Singapore