File No.: CE-007W-01 02
Date of Issue: 2021.01.14

Disposable Anesthesia Breathing Circuit

Instruction for Use

Intended Use

This product is used for anesthesiology, ICU and emergency room use to transfer medical air.

Type

Smoothbore Circuit, Corrugated Circuit, Collapsible Circuit, Coaxial Circuit, Catheter Mount, Duo-Limb circuit, Smoothbore Circuit (Non-sterile), Corrugated Circuit (Non-sterile), Collapsible Circuit (Non-sterile), Coaxial Circuit (Non-sterile), Catheter Mount (Non-sterile), Duo-Limb circuit (Non-sterile)

Precautions

1. Check if the product is within the period of validity before use.
2. Do not use the product if the packing is damaged.
3. Check if all connectors are secure after opening the packing.
4. Check the connectors and circuits which connect to anesthesia machine or breathing machine, ensure there is no leakage.
5. This product is for single use only, dispose it after use.

Directions for Use

1. Open the package by tearing off from the bottom edge V-shape opening, DO NOT use scissors or cutter to open the package to avoid product damage
2. Before connecting the product to patient, securely attach all connections by using a twisting action to ensure a leak free and tight fit
3. Following the machine instruction to test the product, make sure the product is compatible with the machine used and confirm that the product has no leakage and occlusion
4. Make sure the ports of the product such as gas sampling port, pressure port and temperature port are closed when not connected to the monitoring line
5. Patient should be monitored constantly whenever the product is in use

Warnings

Do not reuse, the product is **FOR SINGLE USE ONLY**，reuse of this product can result in cross-infection.

Storage

Storage relative humidity below 80%, clean, dark, non-corrosive gases, ventilated clean room.

Contraindications

There are no known contraindications.

Shelf life

3 years.

Symbol Definition

				
DO NOT REUSE	USE BY	BATCH CODE	DATE OF MANUFACTURE	STERILIZED USING ETHYLENE OXIDE
				
MANUFACTURER	CAUTION	AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY	CONSULT INSTRUCTIONS FOR USE	CATALOGUE NUMBER
				
KEEP DRY	DO NOT USE IF PACKAGE IS DAMAGED	KEEP AWAY FROM SUNLIGHT	CE Mark	

Imported and Distributed by:

PT. IDS Medical Systems Indonesia
IDS Medical Systems (Singapore) Pte. Ltd.
IDS Medical Systems (Hong Kong) Co. Ltd.
IDS Medical Systems (Hong Kong) Co. Ltd. (Taiwan Branch)
IDS Medical Systems (M) Sdn Bhd
IDS Medical Systems (Thailand) Company Ltd.
IDS Medical Systems Philippines, Inc.
IDS Medical Systems (Vietnam) Limited



Zhejiang Haisheng Medical Device Co., Ltd.



Shanghai International Holding Corp. GmbH (Europe)
Eiffestrasse 80,20537 Hamburg Germany