HUDSON RCI®

AQUAPAK®

REF 400400, 400405, 403700, 403705, 403709, 404000, 404045 — STERILE A

REF 403133, 403128, 031-28F — STERILE R

REF 403728, 404128, 404428, 044-28F, 037-28F, 404433, 404135, 404435, 404133, 403770, 403735, 403733, 403739 — STERILE R / STERILE A

Rx ONLY

2 STERILIZE 0°C PHT DEHP

CE 0120 EC REP Teleflex Medical IDA Business and Technology Park Dublin Rd. Athlone, Co. Westmeath Ireland

Teleflex MEDICAL
3015 Carrington Mill Blvd
Morrisville, NC 27560 USA
USA: (866) 246 6990
International: +1 919 544 8000

© Copyright 2007, 2010, 2014, 2016, 2018 Teleflex
L02485 R00

en

Prefilled Sterile Nebulizer

Indication for use: Adds sterile water or saline solution in aerosol form to a patient's breathing gases.

Instructions for use:

NOTE: Use only Hudson RCI AQUAPAK® sterile nebulizer reservoirs.

1. Thread the Nebulizer Adaptor onto the reservoir and tighten.
2. Grasp the return tube behind the thumb flanges on the puncture pin and remove the pin sheath.
3. Insert the pin into the puncture port on the reservoir and push until the thumb flanges rest against the reservoir.
4. Thread the nebulizer onto the flow-meter.
5. Attach the patient tube to the nebulizer outlet.
6. Turn the entrainment collar on the nebulizer to the desired oxygen concentration.
7. Set the flow-meter to the desired setting and verify the flow of gas through the device.
8. The nebulizer is ready for use when liquid flows through the return tube.

Warnings:

The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Single use: Do not reuse, reprocess or re-sterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.

Contents of reservoir for inhalation therapy only, not for parenteral administration. Contains no antimicrobial agents.

Do not use if package is damaged.

PHT DEHP **Contains DEHP**

AQUAPAK and Hudson RCI are registered trademarks of Teleflex Incorporated or its affiliates.

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