VADI

Connector

REF G-311019-2

MD

Reusable Connector

Intended use : use to connect breathing tube for delivering medical gas from ventilator to patient.

Material: Polysulfone

Procedure:

1.Connect the connector to the breathing circuit.

2. Ensure all connections are secure.

3. Place the test lung on the outflow port of the circuit and test the function of the ventilator before connect to the patient.

4. Then the circuit can be connected to the patient.

Caution:

1. Check all connects are tight before use.

2. Please follow waste disposal guideline locally

3. This product is used under physician's order.

4. Please follow the infection guideline to use the device.

5. Please follow the attached to do the clean, disinfection and sterilization information.

6. Allow autoclave maximum 20 times.

VADI Medical Technology Co., Ltd. Yangmei
No. 198, Lane 298, Huandong Road, Zhongshan Village,
Yangmei District, Taoyuan City 32665, Taiwan
E-mail: sales@vadimedical.com.tw TEL: +886-3-475-9911
Website: www.vadimedical.com.tw FAX: +886-3-475-9922

Please see the seal.

Please see the seal.

LOT Please see the seal.

Made In Taiwan

VA-1-G-311019-2.CFS_01