| Language | Page |
|---|---|
| English | 1 |
| Hrvatski | 2 |
| Srpski | 3 |
| Slovenščina | 4 |
| Русский | 5 |

# MEDRAD® Salient
# MEDRAD® Vistron Plus
**Contrast Injection System**

## Connector Tubing
## INSTRUCTIONS FOR USE

| Catalog No. | Description |
|---|---|
| ZY5154 | LPCT with Dual Check Valves |

Shelf Life: 5 years from date of sterilisation
Primary Packaging: Tyvek/Nylon+PE blister pack
Secondary Packaging: Cardboard shipper box (50 per)
Sterilization Method: Ethylene Oxide (EtO)
Internal volume: 5.4ml
Pressure: 300 psi
Latex Free
Plasticizer: ATBC

**Imaxeon Pty Ltd**
Unit 1, 38-46 South Street
Rydalmere NSW 2116
Australia

EC REP **Bayer Medical Care B.V.**
Avenue Céramique 27
6221 KV Maastricht
The Netherlands

CE 2797

© 2018-2020 Bayer. This material may not be reproduced, displayed, modified or distributed without the express prior written consent of Bayer.
Bayer, the Bayer Cross, Imaxeon, MEDRAD FluiDots, MEDRAD and FluiDots are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries.
Bayer reserves the right to modify the specifications and features described herein or to discontinue any product or service identified in this publication at any time without prior notice or obligation. Please contact your authorized Bayer representative for the most current information.


| Symbol | Description |
|---|---|
| | Consult instructions for use (ISO 15223-1, 5.4.3) |
| | Non-Pyrogenic Fluid Path (ISO 15223-1, 5.6.2 & 5.6.3) |
| | Do not use if package is opened or damaged (ISO 15223-1, 5.2.8) |
| | Single Use Only (ISO 15223-1, 5.4.2) |
| | Date of Manufacture/Sterilisation (ISO 15223-1, 5.1.3) |
| | Use By (ISO 15223-1, 5.1.4) |
| LOT | Lot Number (ISO 15223-1, 5.1.5) |
| REF | Catalog Number (ISO 15223-1, 5.1.6) |
| STERILE EO | Sterilised with Ethylene Oxide (ISO 15223-1, 5.2.3) |
| LATEX | Not made with natural rubber latex (ISO 15223-1, 5.4.5, Annex B, B.2) |
| CE 2797 | Indicates the device conforms to the requirements of the European Medical Device Directive 93/42/EEC |
| **Warning** | **Warning** - Advises you of circumstances that could result in injury or death to the patient or the operator |
| **Caution** | **Caution** - Advises you of circumstances that could result in damage to the device |
| | Quantity (IEC TR 60878, 2794) |
| | Manufacturer (ISO 15223-1, 5.1.1) |
| | Fragile, handle with care (ISO 15223-1, 5.3.1) |
| | Keep dry (ISO 15223-1, 5.3.4) |
| | This side up (ISO 7000-0623, ISO 780- No. 3) |
| | Temperature range Primary label: Storage Secondary label: Transport (ISO 15223-1, 5.3.7) |
| | Humidity range Primary label: Storage Secondary label: Transport (ISO 15223-1, 5.3.8) |
| | Do Not Resterilize (ISO 15223-1, 5.2.6) |
| | Attention: Refer to warnings and cautions on Instructions for Use packaged in each carton. (ISO 15223-1, 5.4.4) |

**Introduction:**
Read the information contained in this section. Understanding the information will assist you in operating the device in a safe manner.
**Important Safety Notice:** This device is intended to be used by individuals with adequate training and experience in diagnostic imaging studies.
**Intended Use:** The contents of this package are intended to be used in the delivery of contrast media or saline. They are indicated for single-use on one patient only.
**Contraindications:** These devices are not intended for multiple patient use, drug infusion, chemotherapy, or any other use for which the device is not indicated.
**Restricted Sale:** Rx Only

⚠ **Warnings**
- **Air embolization can cause death or serious injury to the patient.** Do not connect a patient to the injector until all trapped air has been cleared from the syringe and fluid path.
- **Biological contamination can result from reusing disposable items or failure to follow aseptic technique.** Properly discard disposable items after use, or if there is any possibility that contamination may have occurred.
- **Patient or operator injury may result if package is opened or damaged, or if damaged components are used.** Visually inspect contents and package before each use.
- **Patient or operator injury may result from contrast media leaks or tubing ruptures**. Ensure that the fluid path is open; do not exceed pressures identified on the front of the package. Use of greater pressures or occlusions in the fluid path may result in leaks or ruptures.


⚠ **Cautions**
- Component damage or leaks may occur if not installed properly. Ensure all connections are secure; do not overtighten. This will help minimize leaks, disconnection, and component damage.
- Refer to the Injector Operation Manual for further instructions.

**Connector Tube One-Way Check Valves**
The tubing incorporates one-way check valves in both short ends of the tubing. This ensures fluid flows only from the connected syringe to the distal end of the tubing.

**Connector Tube Installation**
1. Ensure all air is purged from the syringe.
2. Remove the connector tube from the package. Remove the dust covers on the luer fittings. Attach the connector tube to the syringe, 1/4 turn to 1/2 turn maximum. Do not install with excessive force.
3. Ensure that the connector luer fitting is secured to the tip of the syringe and verify that the tubing is not kinked or obstructed.

Attach the straight portion of the Y-connector to the contrast (Syringe A) and the extension to the saline (Syringe B).
4. Ensure all air is purged.
5. Rotate the injector head downward.
6. Connect to the patient and press ARM.
7. The injector will ask for CHECK FOR AIR.

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