Rx ONLY CE 2797

PORTEX

en EzPAP® Positive Airway Pressure Therapy System
fr Système de thérapie par pression expiratoire positive EzPAP®
de EzPAP® Atemwegs-Überdrucktherapiesystem
it Sistema terapeutico di ventilazione a pressione positiva EzPAP®
es Sistema de terapia EzPAP® de presión positiva en vías respiratorias
pl Sistema de Terapia de Pressão Positiva Das Vias Respiratórias EzPAP®
nl EzPAP® positief luchtwegdruk-therapiesysteem
sv EzPAP®-system för PAP-behandling (Positive Airway Pressure)
da EzPAP® system til PAP-behandling
no EzPAP® Positive Airway Pressure Therapy System
fi Positiivisen painehoidon EzPAP®-järjestelmä

Manufacturer:
Smiths Medical ASD, Inc.
6000 Nathan Lane North
Minneapolis, MN 55442 USA
Tel: 1 800 258 5361 (US/CA)
Tel: +1 614 210 7300

www.smiths-medical.com

EC REP European Representative:
Smiths Medical Czech Republic a. s.
Olomoucká 306, Hranice 1 - Město,
753 01 Hranice, Czech Republic
Tel: +44 (0)1233 722100

2019-07
10016174-002

smiths medical

en English

1. DESCRIPTION:
The EzPAP® Positive Airway Pressure Therapy System is a non-sterile single patient use, hand held respiratory therapy device. Positive Airway Pressure Therapy System has the ability to provide supplemental oxygen when used with compressed oxygen. The device is also available with various mask sizes.

The EzPAP® system consists of:
1. Positive Airway Pressure Device
 - Ambient Air Inlet
 - Gas Inlet Port
 - Pressure Monitoring Port With Cap
2. Patient Interface -Mouthpiece (or Mask in some models)
3. 7 foot (2.13m) O2 tubing to connect gas inlet port (B) to wall flow met

2. INDICATIONS FOR USE:
The EzPAP® Positive Airway Pressure Therapy System is indicated for use the treatment and prevention of atelectasis. The EzPAP® Positive Airway Pressure Therapy System is for patients who are capable of following directions for Positive Airway Pressure (PAP) therapy.

3. CONTRAINDICATIONS:
Although no absolute contraindications to the use of PAP Therapy have been reported, the following should be carefully evaluated before a decision is made to initiate therapy:
- Inability to tolerate increased work of breathing
- Intracranial pressure (ICP) > 20mm Hg
- Recent facial, oral or skull surgery or trauma
- Esophageal surgery
- Untreated pneumothorax
- Known or suspect hepatic membrane rupture or other middle ear pathology
- Hemodynamic instability
- Acute sinusitis
- Epistaxis
- Active hemoptysis
- Nause

4. WARNINGS:
4.1 Use of this device at excessive pressures may have adverse effects. Expert clinical judgement should be exercised in the selection of settings for each individual patient. If using compressed oxygen as the gas supply, do not exceed the clinically prescribed oxygen flowrate.
4.2 DO NOT OCCLUDE AMBIENT AIR INLET. DO NOT EXCEED 15 LPM GAS INPUT OR 60 PSI COMPRESSED GAS INPUT.
Adverse reactions may include:
Increased work of breathing that may lead to hypoventilation and hypercarbia; increased
cranial pressure; cardiovascular compromise; myocardial ischemia; decreased venous return; air swallowing with increased likelihood of vomiting and aspiration; claustrophobia; skin break down and discomfort from mask; and pulmonary barotrauma.
4.3 This product contains DEHP – data derived from animal and epidemiological studies outline the following risks, but are not limited to: Adverse effects on fetal development, kidney and liver toxicity; negative impact on fertility. Special attention on the application of phthalate products should be given for the following patient populations: pediatrics, pregnant and nursing women. For further information on phthalates please refer to the Smiths Medical website: www.smiths-medical.com/phthalates

5. PRECAUTIONS:
5.1 DO NOT AUTOCLAVE.
5.2 DO NOT IMMERSE COMPLETELY IN ANY LIQUID.
5.3 DO NOT DISASSEMBLE DEVICE/FILTER FOR MOUTHPIECE, MASK OR TUBING.

6. INSTRUCTIONS FOR USE:
6.1 Set-up: Review the diagram and the product to become familiar with all the product features.

6.1.1 Air source:
6.1.1.1 Connect one end of tubing (supplied with the device) to a 0-15 LPM flow meter that is connected to 50-60 psi medical air (oxygen if prescribed).
6.1.1.2 Connect the patient end of the tubing coming from the flow meter to the device's gas inlet port (B on diagram).

6.1.2 Pressure monitoring (recommended):
Connect the device's pressure port (C) to a pressure gauge (with tubing). This is recommended for initial set-up to become aware of pressures transmitted to the patient. If pressure monitoring is not being done, seal port with attached pressure monitoring port cap. Pressure monitoring port should not be exposed to ambient air while in use. Attach the mouthpiece or facemask to rounded end of the device. Serrated end (A) is to be exposed to ambient air at all times - DO NOT OCCLUDE. Inspect device on a routine basis to ensure proper use and function. If damaged, do not use. Verify all connections are secure.

6.2 Initial Flow Meter Setting:
6.2.1 Instruct the patient to relax while performing diaphragmatic breathing.
6.2.2 Set flow rate to 5 LPM on the wall air or oxygen flow meter.
6.2.3 Place mouthpiece in patient's mouth and instruct patient to breathe easily against the pressure from the device. If using a mask, apply mask tightly but comfortably over nose and mouth.
6.2.4 While monitoring airway pressure, slowly adjust flow meter until desired expiratory airway pressure is reached.
6.2.5 Slow inhalations and exhalations by the patient will maintain adequate airway pressure during the breathing cycle.
6.2.6 Maintain breathing through the EzPAP® Positive Airway Pressure Therapy System to accomplish clinical goals.
6.2.7 Several repetitions may be needed to ensure that individual patient needs are being met.
NOTE: For each flow meter setting:
 A. The greater the inspiratory flow through the device - the higher the airway pressure.
 B. The greater the inspiratory flow through the device - the lower the airway pressure.

6.3 EzPAP® Procedure:
6.3.1 Ensure the flow meter is set as identified in previous steps.
6.3.2 Have patient sit with elbows resting comfortably on table.
6.3.3 Place mouthpiece lightly into mouth.
 • Be sure patient maintains a tight seal on the mouthpiece during therapy.
 • Recommend the use of a nose clip, if necessary.
 • If using a mask, apply mask tightly, but comfortably, over nose and mouth.
6.3.4 Instruct the patient to breathe diaphragmatically.
6.3.5 Perform breathing through the device as recommended.
EzPAP® Positive Airway Pressure Therapy System is intended to be used on an intermittent basis (not for prolonged periods of time).
6.3.6 Repeat steps 6.3.2 to 6.3.5 as prescribed.

7. FOR USE WITH NEBULIZER:
Suggested Instructions For Use With Small Volume Nebulizer:
CAUTION: It is important that device instructions for the EzPAP® Positive Airway Pressure Therapy System, as well as the intended nebulizer, are understood prior to use. The selected nebulizer should be used according to the manufacturer's instruction for use in a positive pressure system.
Set-up:
Review the diagrams contained with the devices.
A possible nebulizer and EzPAP® Positive Airway Pressure Therapy System set up is reflected to the right. Follow set-up instructions for each device. Follow cleaning instructions contained with each device. Inspect device on routine basis to ensure proper use and function. If damaged, do not use. Verify all connections are secure.

8. Cleaning Instructions:
Clean as needed. Do not immerse in any liquid. Wash the parts by swabbing with warm soapy water. Rinse and shake dry. Air dry the device by resting the unit on its side with the mouthpiece removed until all liquid droplets are evaporated before continuing use.

9. DISPOSAL:
Dispose of **EzPAP®** Positive Airway Pressure Therapy System in a safe manner according to Federal/State/Local regulations and guidelines for disposal of contaminated medical waste.

10. Temperature Specification: EzPAP® system w/disposable manometer
Operation Temperature: -18°C/0°F to 50°C/122°F
Storage Temperature: -40°C/-40°F to 60°C/140°F

Ordering Information:
Each system includes EzPAP, one 7 ft. oxygen tube, one pressure port cap and one of the following:

EzPAP® System w/Mouthpiece (10 per Case)	23-0747
EzPAP® System w/disposable manometer (10 per Case)	23-0757
EzPAP® System w/Pediatric Mask (1 per Case)	23-1747
EzPAP® System w/Medium Mask (1 per Case)	23-2747
EzPAP® System w/Large Mask (1 per Case)	23-3747
EzPAP® Kit (1 per Case)	23-6000

The Smiths Medical and Portex design marks and EzPAP are trademarks of Smiths Medical.
The symbol ™ indicates the trademarks registered in the U.S. Patent & Trademark Office & certain other countries.

©2013–2019 Smiths Medical. All rights reserved.

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Znaki graficzne Smiths Medical i Portex oraz EzPAP są znakami towarowymi przedsiębiorstw Smiths Medical. Symbol ™ wskazuje, że znak towarowy jest zarejestrowany w amerykańskim urzędzie ds. patentów i znaków towarowych. U.S. Patent and Trademark Office) oraz w niektórych innych krajach.