Operation Manual

MEDRAD® Intego

PET Infusion System

Report any serious incident that has occurred in relation to this device to Bayer (radiology.bayer.com/contact) and to your local European competent authority (or, where applicable, to the appropriate regulatory authority of the country in which the incident has occurred).

1 Important Safety Notice to Users of the MEDRAD® Intego PET Infusion System 1

2 Introduction 3

2.1 Certifications 3
2.2 Intended Use 3
2.3 Indications for Use 3
2.4 Restricted Sale 3
2.5 Disclaimers 3
2.6 Required Training 3
2.7 Definitions 4
2.8 Symbols and Icons 4
2.8.1 Warnings 6
2.8.2 Display Icons 7
2.8.3 Status Indicators 9
2.8.4 Packaging 10

3 System Basics 13

3.1 System Overview 13
3.2 Safety Features 13
3.3 Shielded Chamber Components 15
3.4 System Components 16
3.5 Source Administration Set (SAS) Components 17
3.6 Patient Administration Set (PAS) Components 18
3.7 Vial Shield Components 19
3.8 Calibration Source Holder 19
3.9 Managing Power States 20
3.9.1 Powering On the System 20
3.9.2 Powering Off the System 21
3.9.3 Placing the System Into Standby 21
3.9.4 Taking the System out of Standby 21
3.10 Securing the Intego PET Infusion System 21
3.10.1 Locking the Device 22
3.10.2 Unlocking the Device 22
3.10.3 Locking the Shielded Chamber 23
3.11 Navigating the User Interface 23
3.11.1 Screen Overview 23
3.11.2 System Preparation Screen 24
3.11.3 Schedule Screen 24
3.11.4 Dosing Screen 25
3.11.5 Configuration Screen 25
3.12 Moving the Intego PET Infusion System 26
3.12.1 Moving the System Using the Drive Override 27

4 System Configuration 29

4.1 System Settings 29
4.2 Dosing Settings 31
4.3 Maintenance Settings 34
4.4 Security Settings 35
4.5 Communication Settings 36

5 Daily Setup 39

5.1 Placing the System into Clinical Mode 39
5.2 Using the Schedule 41
5.2.1 Manually Entering Schedule Information 42
5.2.2 Importing the Schedule Information via the USB Port 42
5.2.3 Importing the Schedule from the Modality Worklist 43
5.2.4 Editing the Schedule 44
5.2.5 Exporting the Infusion History 45
5.3 Performing Daily QC 45
5.4 Entering RP Assay and Saline Information 48
5.5 Installing the SAS 50
5.6 Installing the Multi-Dose Vial 57
5.7 Installing the PAS 58
5.8 Priming the SAS 59

6 Patient Infusion 61

6.1 Adding a Patient / Editing Patient Information 61
6.2 Selecting a Patient from the Schedule 62
6.3 Priming the PAS 63
6.4 Selecting the Flow Rate 64
6.5 Performing a Saline Test Injection 65
6.6 Entering the Requested Dose Activity 66
6.6.1 Personalized Dose Entry 66
6.6.2 Manual Dose Entry 66
6.6.3 Requested Activity and the Activity Bar 67
6.7 Preparing an RP Dose 68
6.8 Infusing or Discarding the Dose 69
6.8.1 Start and Monitor Infusion 70
6.8.2 Infusion Completion 71
6.8.3 Discard Completion 72
6.9 Sending Secondary Capture Images to PACS 72
6.10 Monitoring the Vial 73
6.11 Replacing the Saline Bag 74

7 Vial Shield and SAS Removal 77

7.1 Vial Shield and SAS Removal 77

8 Training Mode 79

8.1 Training Mode 79

A Cleaning and Maintenance 83

A.1 Cleaning Guidelines 84
A.2 Recommended Maintenance Schedule 84
A.3 Dose Calibrator Linearity Check 85
A.4 Dose Calibrator Calibration 90
A.5 Geometry Check 93
A.6 Calibrating the Pinch Valves 93
A.7 Replacing the Printer Paper 95

A.8 Battery Maintenance and System Storage....96
A.9 Disposal of Equipment....97

B Specifications....99

B.1 Mechanical....99
B.2 Power Cable Specifications....99
B.3 MEDRAD® Intego PET Infusion System to IT Network Connections....100
B.4 Radiation Shielding Profile....100
B.4.1 Vial Shield Radiation Profile....100
B.4.2 Radiation Exposure Reduction Specification and Testing Results....100
B.5 Environmental....101
B.6 Electrical....102
B.6.1 Protection Against the Ingress of Fluids....102
B.6.2 Mode of Operation....102
B.7 Fluid Delivery....103
B.8 Mobility....105
B.9 MEDRAD® Intego Wireless Connection Feature....106
B.9.1 Required Characteristics....106
B.9.2 Required Configuration....107
B.9.3 Intended Information Flow....107
B.9.4 Network Communication Technical Specifications....108
B.9.5 FCC RF Radiation Exposure Statement....108

C Compliance to IEC 60601-1-2 / 2nd, 3rd, and 4th Editions....109

D Troubleshooting Tips....113

D.1 Recovering from Priming Issues....113
D.2 Recovering from Activation of Air Detector....113
D.2.1 During Priming....113
D.2.2 During a Saline Test Inject....114
D.2.3 During an Infusion....114
D.3 Recovering a Dose Due to a System Failure....114
D.3.1 With a Dose of RP in the Dose Calibrator....114
D.3.2 With RP in the Vial Shield....115
D.4 PAS Occlusion Recovery....116
D.5 System Messages....117
D.5.1 Critical Error Messages....117
D.5.2 Recoverable Errors / Messages....118

E Vials and Vial Shields....123

F Components and Catalog Numbers....133

G Country-Specific Information....135

G.1 Australia....135
G.2 New Zealand....135
G.3 Brazil....135
G.4 South Korea....135
G.5 Thailand....135

G.6 Mexico 135
G.7 Serbia 135
G.8 Singapore 136
G.9 Ukraine 136
G.10 Taiwan 136
G.11 Oman 136

1 Important Safety Notice to Users of the MEDRAD® Intego PET Infusion System

This manual and the equipment it describes are for use only by qualified medical professionals with training and experience in Nuclear Medicine procedures. It is intended as a guide for using both the MEDRAD® Intego (Intego) PET Infusion System and dedicated Intego PET Infusion System disposables.

The Intego PET Infusion System contains a gross air detection feature, which is intended to assist qualified medical professional users/operators during the set up procedure to ensure that all air is out of the system.

The safe and effective use of the Intego PET Infusion System to a large degree depends upon factors solely under the control of the medical professionals using the system. There is no substitute for a properly trained and vigilant device user. It is important that the operating instructions and user warnings supplied with the Intego PET Infusion System be read, understood, and followed.

Before starting any PET infusion procedure, the device user should be trained in the particular procedures to be performed and should be familiar with the medical literature related to procedures and the potential complications and risks verses the benefits of utilizing radiopharmaceutical fluid infusion procedures.

This manual is intended as an extension of the user interface of the Intego PET Infusion System to provide procedural and technical information. Additional Intego training information is available in the following formats:

- On-site in-service sessions
- Service manual
- Package inserts (IFU)

Contact Bayer if any of these resources are needed.

Operating specifications and feature availability may vary by country. Check with your local product representative and country-specific operating instructions.

NOTE: Contact Bayer for Installation information.

2 Introduction

2.1 Certifications

The MEDRAD® Intego PET Infusion System is equipped to operate at 100-240 VAC, 50/60Hz, 300 VA. The Intego PET Infusion System is designed to comply with IEC 60601-1 (2nd and 3rd Edition Amendment 1) and IEC 60601-1-2 (2nd, 3rd, and 4th Edition) standards, including national differences. Special precautions regarding Electro-magnetic Compatibility (EMC) are required for installation and use of this system. Detailed EMC information can be found in "Compliance to IEC 60601-1-2 / 2nd, 3rd, and 4th Editions".

2.2 Intended Use

The Intego PET Infusion System is intended to deliver accurate doses of ^{18}F-Fluorodeoxyglucose (^{18}F-FDG) or ^{18}F-Sodium Fluoride (^{18}F-NaF) radiopharmaceuticals and commonly used flushing solutions to patients during molecular imaging (nuclear medicine) diagnostic procedures. The Intego PET Infusion System is also intended to provide effective radiation shielding to medical personnel from Fluorine-^{18} (^{18}F) radiation exposure during nuclear medicine diagnostic procedures.

2.3 Indications for Use

The Intego PET Infusion System is indicated for the administration of ^{18}F-FDG, ^{18}F-NaF, and commonly used flushing solutions to patients during molecular imaging (nuclear medicine) procedures.

NOTE: The Intego PET Infusion System is intended for use with ^{18}F-FDG or ^{18}F-NaF. Intego PET Infusion System disposables intended for use with ^{18}F-FDG may be used to deliver ^{18}F-FDG or ^{18}F-NaF. Please contact Bayer for more information.

2.4 Restricted Sale

The United States Food and Drug Administration (FDA) restricts sale of this system to physicians or those with written authorization from a physician.

2.5 Disclaimers

External wiring and modifications disclaimers: Bayer disclaims liability for any modifications or interfaces with other equipment that are not in conformity with the specifications and information contained within this manual.

Anyone who connects additional equipment to the device or configures a medical system is responsible that the system complies with the relevant requirements of IEC 60601-1. An accessory or equipment connected to the device must be certified to either IEC 60601-1 (Operator or Patient Environment Use) or, outside the patient environment, the level of safety must be equivalent to equipment complying with their respective IEC or ISO safety standards, e.g. IEC 62368-1 or IEC 60950-1 (Operator Environment Use Only), and must comply with the relevant requirements according to IEC 60601-1. Consult Bayer for any modifications to the equipment.

This manual applies to the Intego PET Infusion System, Catalog Number INT SYS 200. Read all the information contained in this manual. Understanding this information assists in operating the Intego PET Infusion System in a safe manner.

Screen images in this manual are for illustration purposes only. Actual screens may vary.

2.6 Required Training

This device is intended to be used by individuals with training and experience in nuclear imaging studies.

2.7 Definitions

The following are definitions of the terms WARNING, CAUTION, and NOTE found throughout this document:

⚠WARNING	Indicates that the information is a warning. Warnings advise of circumstances that could result in injury or death to the patient or operator. Read and understand the warnings before operating the Intego PET Infusion System.
⚠CAUTION	Indicates that the information is a caution. Cautions advise of circumstances that could result in damage to the system or improper functioning of the system. Read and understand the cautions before operating the Intego PET Infusion System.
NOTICE	Indicates that the information is a notice. Notices advise you of circumstances that could result in damage to the device. Read and understand the notices before operating the Intego PET Infusion System.
Note	Indicates that the information that follows is additional important information, a tip that helps the clinician to recover from an error, or points to related information within the manual.

2.8 Symbols and Icons

The symbols and icons discussed in the sections below describe the requirements to which the Intego PET Infusion System conforms, how warnings are displayed in the manual, and the icons used on the equipment and equipment packaging.

CLASSIFIED c UL us E353299	Medical - General Medical Equipment As To Electrical Shock, Fire, and Mechanical Hazards Only In accordance with ANSI/AAMI ES60601-1 (2005) + AMD 1 (2012) CAN/CSA-C22.2 No. 60601-1 (2014)
CE 0086	Indicates that this system conforms to requirements of the European Medical Device Directive 93/42/EEC.
IPX1	IPX1 Code that specifies the degree of protection provided by the enclosure against vertically falling water drops (IEC 60529).
	Indicates separate collection for Electrical and Electronic Equipment per Directive 2002/96/EC. Refer to the following website for additional information: www.weee.bayer.com
CLASS 1	Indicates the Intego PET Infusion System is a Class 1 electrical device as determined by IEC 60601-1 standards.
RP	Radiopharmaceutical (^{18}F-FDG or ^{18}F-NaF).
	Stop button (IEC TR 60878, 5110A)
	On/Shutdown/Standby button (IEC TR 60878, 5009)

Identifies a type BF applied part complying with IEC 60601-1 standards (IEC 60417, 5333)

Fuse rating. See Intego address label for specific fuse rating information. (IEC TR 60878, 5016)

Consult instructions for use. (ISO 15223-1, 5.4.3)


See accompanying documentation. This symbol indicates the user shall refer to the instructions-for-use to ensure safe operation. (ISO 7010, M002)


This symbol informs the user of the correct manner of use.


This symbol informs the user of the incorrect manner of use.


Label showing the correct and incorrect routing of the Source Administration Sets (SAS) through the Saline Pump.


Label showing the correct and incorrect insertion of the Saline Spike into the Saline Container.


Label showing the correct and incorrect routing of the tubing through the RP Pump.


Label showing the correct and incorrect direction to attach the Needle Cartridge into the Needle Cartridge Holder.


Label showing the correct and incorrect insertion of the tubing into the Air Detector.

Label showing the correct and incorrect orientation of the T-Connector and the tubing in each of the Pinch Valves.

Label showing the correct and incorrect method of inserting the tubing into a Pinch Valve.

Label showing the correct method for installing a Pinch Valve Cover onto a Pinch Valve.

Indicates Drive Speed Switch is in fast speed range position.

Indicates Drive Speed Switch is in slow speed range position.

2.8.1 Warnings

Warning: Indicates a pinch or crush hazard. (ISO 7000, W024)

Crush hazard. This symbol indicates there is the potential for crush injury.

Warning: Indicates hazardous voltages. (ISO 7010, W012)


Radiation exposure hazard. Both symbols indicate that opening the Shielded Chamber when there is RP in the Intego PET Infusion System could expose the operator or patient to radiation. (ISO 7010, W003)

DO NOT place greater than 1000 ml Saline Container on the Saline Hook.

Sitting on the unit is prohibited. (ISO 7010, P018)


Maximum weight of the injector system and accessories during normal use.
(ISO 7000, 1321B; ISO 15223-1, 5.4.3)

2.8.2 Display Icons

Radiation activity icon:

- Not highlighted in yellow- Radioactivity may or may not be present within the Intego PET Infusion System
- Highlighted in yellow - Radioactivity is present within the Dose Calibrator and assay information has been entered into the Intego PET Infusion System

Fluid delivery battery status icon:

- Not highlighted in yellow – Battery not present or battery is completely depleted.
- Highlighted in yellow – Unit is operating on battery power

Plug icon:

- Battery is present but the unit is operating on AC power/charging

Fluid icon:

- If present, fluid is currently being pumped within the Intego PET Infusion System

Wireless status icon:

- Not highlighted in yellow - Wireless is not connected
- Highlighted in yellow - Wireless is connected

Calendar icon. Touch this icon to enter the desired date.

Day selector icon. Toggles between today and yesterday.

Configuration icon. Located on the configuration button; used to open the Configuration screen.

Reset icon. Located on the reset button. Resets the RP field to the configured default value or the current scheduled dose and resets the patient infusion process back to entering patient information.

AM/PM icon. Press this icon to set AM or PM when entering time data.

Approximately icon. Used to identify estimated values.

Attention icon. Used to identify items that require clinician attention.

Partial infusion delivered.


Complete infusion delivered.


Add appointment icon. Located on the ADD APPOINTMENT button. Used to add an appointment to the schedule.

Import worklist icon. Located on the IMPORT WORKLIST button. Used to import the worklist.

Import dose plan icon. Located on the IMPORT DOSE PLAN button. Used to import a dose plan from a USB memory device.

Export history icon. Located on the EXPORT HISTORY button. Used to export infusion history to a USB memory device.

Clear schedule icon. Located on the CLEAR SCHEDULE button. Used to delete the entire schedule.

HIDE/UNHIDE appointment button. Used to hide or unhide an appointment from view on the Dosing screens.

Shows how many appointments are hidden.

Transmission in progress icon. Used to indicate secondary captures for the patient are in transmission to PACS.

Transmission successful icon. Used to indicate secondary captures for the patient have been successfully transmitted to PACS.


Transmission error icon. Used to indicate an error in the transmission of secondary captures to PACS for the patient.
The system will retry information transmission three times. There is a 30 second delay between tries. If the system cannot upload the data after three attempts, this icon will display.

Transmission not possible icon. Used to indicate transmission of secondary captures to PACS for the patient is currently not possible.
Transferring to PACS will only be attempted if Patient Name, Patient ID, and Accession Number are entered. If one of these is not entered, this icon will display.


2.8.3 Status Indicators

The following Status Indicators are used on all Intego PET Infusion Systems:

System Power and Dose Calibrator Status (Green) (IEC TR 60878, 5406)

- On - Intego PET Infusion System is On and Dose Calibrator is ready for use
- Blinking - Intego PET Infusion System is On and Dose Calibrator is warming up
- Off - Intego PET Infusion System is shut down

Fluid Delivery Battery Backup Status (Amber)

- On - Intego PET Infusion System is using battery backup
- Blinking - Only 5 minutes or less remain on battery before the Intego PET Infusion System will completely shut down. Connect the Intego PET Infusion System to AC power.
- Off - Intego PET Infusion System is not using battery backup

Drive System Status (Blue)

- On - Drive System is available for use
- Blinking - Low Drive System battery - Approximately 3 minutes or less remain before the Drive System becomes unavailable for use. Connect the Intego PET Infusion System to AC power.
- Off - Drive System is not available for use

NOTE: The Shielded Chamber Lid must be closed and latched in order to use the Drive System.

System Battery Charging Status (Violet)

- On - Batteries are charging
- Off - Batteries are not charging

2.8.4 Packaging

Catalog Number (ISO 15223-1, 5.1.6)

Serial Number (ISO 15223-1, 5.1.7)

Attention: Refer to warnings and cautions on Instructions for Use packaged in each carton. (ISO 15223-1, 5.4.4)

Single use only (ISO 15223-1, 5.4.2)

For use with one vial of media only

Do not use if package is opened or damaged (ISO 15223-1, 5.2.8)

Lot Number (ISO 15223-1, 5.1.5)

Date of Manufacture (ISO 15223-1, 5.1.3)

Non-pyrogenic Fluid Path (ISO 15223-1: 5.6.2, 5.6.3)

Sterilized using irradiation (ISO 15223-1, 5.2.4)

Device contains a sterile fluid path that has been sterilized using ethylene oxide (ISO 15223-1, 5.2.9)

Use by (ISO 15223-1, 5.1.4)

Atmospheric pressure range (ISO 15223-1, 5.3.9)

Humidity range (ISO 15223-1, 5.3.8)

Temperature range (ISO 15223-1, 5.3.7)

Device contains phthalates (BS EN 15986:2011)

Do not stack (IEC TR 60878, 2402)

This side up (ISO 7000, 0623)

Keep dry (ISO 15223-1, 5.3.4)

Fragile, handle with care (ISO 15223-1, 5.3.1)

Manufacturer (ISO 15223-1, 5.1.1)

Authorized representative in the European community (ISO 15223-1, 5.1.2)

3 System Basics

3.1 System Overview

The Intego PET Infusion System delivers ^{18}F-FDG or ^{18}F-NaF to patients during a PET or PET/CT diagnostic procedure. In addition, the system provides effective radiation shielding to medical personnel from ^{18}F radiation exposure during nuclear medicine diagnostic procedures.

The Intego PET Infusion System meets the following clinical needs:

- For a typical 555 MBq (15 mCi) infusion per patient, it limits ^{18}F radiation exposure for medical personnel to less than 60 µSv (6 mRem) finger dose and 3 µSv (0.3 mRem) whole body dose.
- Flexibility to program the required dose either by activity only or by activity per patient weight.
- Ability to deliver ^{18}F-RP within ±10% of the prescribed dose and within ±2% of the measured dose, excluding Dose Calibrator calibration factor. (If using Dose Calibrator accuracy of ±3%, then able to deliver ^{18}F-RP within ±4.7% of the measured dose.)
- Capability to retain and print infusion history and dispensing records.

The Intego PET Infusion System is a self-contained, self-powered mobile cart. A Multi-Dose Vial of RP (up to 27.75 GBq (750 mCi)) is stored within a Shielded Chamber within the body of the system. A SAS is installed within the Shielded Chamber at the same time a new Multi-Dose Vial of the RP indicated for use is installed. Just prior to an infusion, the system measures a dose of RP along with a saline flush in the Dose Calibrator. Once the correct radiation level is achieved, the RP dose and saline are infused into the patient via a PAS. Rechargeable batteries provide sufficient power to keep the Dose Calibrator warm so that the system can be unplugged and moved to a new location without having to completely power Off.

The Intego PET Infusion System consists of the following:

- RP Pump
- Saline Pump
- Dose Calibrator
- Air Detector
- System Shielding
- Vial Shielding

NOTE: Radiation shielding performance is achieved by using the Vial Shield designed by Bayer.

3.2 Safety Features

The Intego PET Infusion System incorporates the following safety features to help protect patients and operators while the system is in use.


NOTE: These features are intended to augment the safety program of a site.

- The RP is stored in a Shielded Chamber within the Intego PET Infusion System. A patient's dose of the RP indicated for use is created automatically while the dose remains within the Shielded Chamber, greatly reducing the operator's exposure to radiation. Any waste material left over after an infusion remains within the Shielded Chamber until time for disposal. By employing a SAS, which is also contained within the Shielded Chamber, it is not necessary to replace the tubing set directly connected to the Multi-Dose Vial for each patient.
- Prior to infusing a dose of the RP indicated for use into the patient, the Intego PET Infusion System measures the dose activity to ensure that the correct dosage will be infused.
- The Intego PET Infusion System contains a Waste Container within the Shielded Chamber. If it is necessary to discard a dose of RP, the dose is transferred into this Waste Container to help prevent radiation exposure to the patient or operator.

NOTE: Estimates of the amount of activity in the waste container represent a worst-case value. The amount of activity in the waste container may be lower than the estimate shown on the display.


- The Shielded Chamber Lid is secured by a heavy-duty latching system, which reduces the likelihood of unintentionally opening the Shielded Chamber while radioactivity is contained within. In addition, the Shielded Chamber Lid can be locked to prevent unauthorized access to the Shielded Chamber.
- The Intego PET Infusion System Display can be locked to prevent unauthorized access to the operating system.
- The Intego PET Infusion System has a saline test inject feature that may be used to check vein patency prior to an infusion.
- The Intego PET Infusion System features an air detection system that automatically disarms the Intego PET Infusion System if air is detected in the SAS.
- The Intego PET Infusion System PAS features a one-way check valve that prevents backflow of fluids into the SAS.
- The Intego PET Infusion System has its own powered Drive System that is controlled by the operator. When the system is preparing a patient dose or infusing fluid to the patient, the Drive mechanism is automatically disabled to prevent unintended movement of the Infusion System.

3.3 Shielded Chamber Components

#	Name	Description
1	SAS Track	Recessed path used to hold the SAS and to prevent it from being damaged as the Shielded Chamber Lid is opened or closed.
2	Vial Shield Compartment	The chamber where the Vial Shield is placed.
3	Needle Insertion Device	A tungsten shielded lid with a fold-down handle that holds the SAS Needle Cartridge.
4	RP Pump	Precision pump used for RP dose preparation.
5	SAS Confluence Holder	Secures the SAS Confluence in the correct location and orientation.
6	T-Connector Holder	Secures the SAS T-Connector in the correct location and orientation.
7	Dose Calibrator	
8	Waste Pinch Valve	Controls fluid flow to the Waste Container.
9	PAS Pinch Valve	Controls fluid flow to the PAS.
10	Waste Storage	The area where the SAS Waste Container is placed.
11	Air Detector	Used to detect air in the fluid as it passes from the SAS into the PAS.
12	Air Detector Holder	Secures the SAS within the Air Detector.
13	Swabbable Valve Holder	Secures the Swabbable Valve in the correct location for joining the PAS and SAS together.
14	Saline Tube Holder	Holds the Saline Tube in the PAS Compartment.


3.4 System Components

#	Name	Description
1	Display	Color touch screen display used to control the Intego PET Infusion System, view delivery data, and report radioactivity within the system.
2	USB Port	Provides the operator the capability of importing and exporting data.
3	Shielded Chamber	Lead shielded compartment that houses the SAS, RP, and various other system components.
4	Shielded Chamber Lid	A sliding lid that allows access to the components and RP contained in the Shielded Chamber.
5	Shielded Chamber Lid Latch	Used to secure the Shielded Chamber Lid from unintentional opening.
6	PAS Compartment	Not Shown - A separate compartment of the Shielded Chamber. Provides shielding for the PAS (located under the PAS Compartment Cover)
7	PAS Compartment Cover	Hinged cover to access the PAS Compartment.
8	On/Shutdown/Standby Button	Turns the Intego PET Infusion System On and Off. Also used to place the Intego PET Infusion System into and out of Standby. Also used to lock and unlock the system.
9	Stop Button	Used to stop current operation of the system. (E.g. When infusing a patient, pressing the red stop button will pause the infusion. It will not discontinue the procedure. When transporting the system, pressing the red stop button will re-set the motor drive system).
10	Power Switch	Rocker switch used to turn AC power On or Off.
11	Saline Pump	Peristaltic pump used to circulate saline in the SAS.
12	Saline Hook	Used for hanging the Saline Container.
13	Printer	Used to print patient-specific infusion records.

#	Name	Description
14	Cable Storage	Posts used for wrapping the power cable when the Intego PET Infusion System is being moved or when not in use.
15	Caster Brakes	Used to secure the Intego PET Infusion System from movement when parked.
16	Drive Controller	Thumb wheel to control forward/reverse variable speed.
17	Drive Engage Switch	Used to enable the Drive System
18	Drive Speed Switch	Used to select the speed range of forward movement.
19	Drive Override	Mechanism to override the Drive System.
20	Status Indicators	Provides battery, Dose Calibrator, and Drive System status.

3.5 Source Administration Set (SAS) Components



#	Name	Description
1	Saline Spike	Inserted into the Saline Container.
2	Saline Tube	Section of tube that is placed along the SAS Track and into the Saline Pump.
3	RP Tube	Section of tube from the Needle Cartridge to the Confluence that is placed in the RP Pump.
4	Needle Cartridge	Contains the needs used to draw the RP into the SAS.
5	SAS Coil	Inserted into the Dose Calibrator. It holds the RP dose being measured by the Dose Calibrator.
6	Pre-Coil Tube	Section of tube from the Confluence to the SAS Coil.
7	Post-Coil Tube	Section of tube from the SAS Coil to the T-connector.
8	Waste Container	Collects any discarded fluid.
9	T-Connector	Union of Waste Tube, Patient Tube, and Post-Coil Tube.
10	Waste Tube	Section of tube from the T-connector to the Waste Container.
11	Patient Tube	Section of tube from the T-connector to the Swabbable Valve.
12	Swabbable Valve	Enables aseptic connection between the SAS and PAS.
13	Saline Disconnect Luer	Disconnects the Saline Tube from the Saline Container.
14	Confluence	Union of RP Tube, Saline Tube, and Pre-Coil Tube.

#	Name	Description
15	Pinch Valve Cover	Aids in proper installation of tubing through pinch valve.


3.6 Patient Administration Set (PAS) Components

#	Name	Description
1	Removable Prime Tube	Captures excess saline during PAS priming.
2	Integrated Filter	The filter may or may not be integrated into the PAS. PAS without integrated filter (user must manually attach filter): 
3	One-Way Check Valve	Connected to the patient catheter or needle. Prevents fluid backflow.
4	PAS Tubing	
5	Luer Connector	Connects the PAS to the SAS Swabbable Valve.


3.7 Vial Shield Components

#	Name	Description
1	Vial Shield	Multi-dose vial radiation shielding vessel designed to work specifically with the Intego PET Infusion System..
2	Vial Cap	Removable cap that provides access to the RP vial.
3	Access Cap	Removable cap that provides access to the RP vial septum.
4	Carrying Handle	Removable handle for transporting the Vial Shield and removing the Access Cap

NOTE: For more information regarding Vial Shields, refer to "Vials and Vial Shields".

3.8 Calibration Source Holder

Name	Description
Calibration Source Holder	Used to insert calibration sources into the Intego PET Infusion System Dose Calibrator

3.9 Managing Power States

3.9.1 Powering On the System

⚠WARNING

To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth.

Fire hazard. To avoid an electrical fire, ensure the correct type of fuse is used for replacement. The fuse must be replaced by qualified personnel only. Refer to Section B.6 for fuse type.


⚠CAUTION

Do not use the Intego PET Infusion System immediately after it has been brought indoors from extreme outside temperature.

- Condensation may cause electrical damage to the system. Allow the Intego PET Infusion System to stabilize at room temperature before use.

1. Plug the Intego PET Infusion System power cable into an AC outlet and move the Power Switch to the On position.
2. Press and hold the ON/SHUTDOWN/STANDBY button until a beep sounds. The Introduction screen will appear.

3.9.2 Powering Off the System

1. Press the ON/SHUTDOWN/STANDBY button.
2. A dialog box appears on the screen with the message “Intego Power and Security Options.” Press the SHUTDOWN button.

NOTE: If the Display is not functional, the ON/SHUTDOWN/STANDBY button can be pressed and held for 15 seconds to power Off the Intego PET Infusion System.

3.9.3 Placing the System Into Standby

1. Press the ON/SHUTDOWN/STANDBY button.
2. A dialog box appears on the screen with the message “Intego Power and Security Options.” Press the STANDBY button.

NOTE: Placing the system into Standby eliminates the need to wait for the Dose Calibrator to warm up and allows the system to be driven. If the system is Shutdown, the user will need to wait up to 15 minutes to use the Dose Calibrator and the system cannot be driven.

NOTE: If the Intego PET Infusion System is in Standby and the power cable is not plugged into AC power within 30 minutes, the system will completely power Off.

NOTE: If the power cable is removed from an AC outlet, the Display will remain powered for 60 minutes. After 60 minutes, the system will automatically go into Standby.

3.9.4 Taking the System out of Standby

1. Press the ON/SHUTDOWN/STANDBY button to activate the Intego PET Infusion System.


3.10 Securing the Intego PET Infusion System

The Device Lock feature secures the display and disables the Drive System. A Standard lock can be used to secure the Shielded Chamber.

3.10.1 Locking the Device


NOTE: To enable the Device Lock feature, see "System Configuration".

1. Press the ON/SHUTDOWN/STANDBY button.
2. A dialog box appears on the screen with the message “Intego Power and Security Options.” Press the LOCK button. A dialog box appears on the screen with the message “Device Locked.

3.10.2 Unlocking the Device

1. If the Intego PET Infusion System is powered Off or in Standby, press the ON/SHUTDOWN/STANDBY button.
2. In the “Device Locked” dialog box, enter the four digit lock code and press the ENTER button.

NOTE: The factory set lock code is 7237. It can be used to override any user configured code.

3.10.3 Locking the Shielded Chamber


A standard lock can be inserted in the mechanism shown below, if locking the Shielded Chamber is desired.

3.11 Navigating the User Interface

This section provides an overview of the main screens the operator will use with the Intego PET Infusion System.


3.11.1 Screen Overview

#	Name	Description
1	Title Bar	Located at the top of the screen. It identifies the current screen: System Preparation, Schedule, Dosing, or Configuration.
2	Status Bar	Located below the title bar. Information about the current state during operation is displayed here.
3	Navigation Bar	Located at the bottom left of the screen. It is used to move among the System Preparation, Schedule, and Dosing screens. Press the SALINE button or the RP button to access the System Preparation screen.
4	Configuration Button	Located at the bottom center of the screen. It is used to access the Configuration screen.
5	Reset Button	Appears only on the Dosing screen. Located at the bottom right of the screen. It resets the RP field to the configured default value or the current scheduled dose and returns the user to the Patient Information panel.
6	Status Icons	Located at the top right of the screen. Refer to the "Symbols and Icons" section.
7	Battery Meter	Located at the top right of the screen. Indicates the amount of time that the system will operate on battery power.

3.11.2 System Preparation Screen


This screen enables the user to enter RP assay and saline information. The Multi-Dose Vial data, Saline Container remaining volume, activity in the SAS tubing, and the volume and activity contained in the Waste Container are monitored from the System Preparation screen. Daily QC and SAS priming are initiated from this screen.

3.11.3 Schedule Screen

This screen enables the user to add appointments to the schedule, import a dose plan via a USB memory device, export the patient infusion history via a USB memory device, and clear appointments from the schedule.

The user can wirelessly import the modality worklist from the Schedule screen if the worklist software feature has been licensed.

3.11.4 Dosing Screen

This screen enables the user to monitor the patient schedule, prepare and deliver the patient dose using Personalized Dosing or Manual Dosing, enter patient and operator data, prime the PAS, perform test injections, change the flow rate, and monitor the RP activity being delivered to the patient.


The user can monitor the transmission status of secondary captures to PACS from the Dosing screen if the PACS software features have been licensed.

3.11.5 Configuration Screen

This screen enables the user to set preferences, such as language, date and time, units, audio level, default values, lock code, and Personalized Dosing formulas. Daily QC source assay information is entered on the Configuration screen. Dose Calibrator quarterly linearity and yearly calibration, and pinch valve calibration are performed from the Configuration screen.

The user can manage communications for the worklist and PACS from the Configuration screen if software features have been licensed.

NOTE: The Navigation Bar does not appear on the Configuration screen. To navigate from the Configuration screen, it is necessary to press the OK button. The user will be returned to the screen being used prior to the Configuration screen.

3.12 Moving the Intego PET Infusion System

⚠WARNINGS

Injury to patient may occur if system is moved without disconnecting the PAS from the patient.

- Disconnect the PAS from the patient prior to moving the Intego PET Infusion System.

Crush hazard. Using the Drive Override disengages the Drive System brakes.

- Ensure the Intego PET Infusion System is not on an incline or left unattended when using the Drive Override and re-engage the Drive Override after moving the system. Engage the Caster Brakes.

Injury or damage to the system may occur if it is moved without the Shielded Chamber Lid secured.

- Do not move system unless the Shielded Chamber Lid is latched.

Injury or damage may occur if the Intego PET Infusion System is used on steep inclines.

- Do not use the system on a floor with an incline greater than 10 degrees.

Radiation exposure hazard.

- Do not move the Intego PET Infusion System unless the Shielded Chamber Lid is latched.

Trip hazard.

- Injury may occur. Do not operate Drive System if the power cable is plugged into an AC outlet.

The system is not to be serviced or maintained while in use with a patient.

⚠CAUTIONS

Do not plug the Intego PET Infusion System power cable into an extension cord or multi-outlet power strip.

- Only use the power cable supplied with the system.

Do not hang clothing or other pieces of equipment from the Intego PET Infusion System.

- Hanging clothes or equipment may impair the movement of the system.

Do not move any portion of the Intego PET Infusion System by force without disengaging the caster brakes or brake release.

Avoid making fast turns or sudden starts and stops.

1. Disconnect the PAS from the patient.
2. Unplug the power cable from the AC outlet and wrap the power cable on the Cable Storage.

 NOTE: If the Intego PET Infusion System is not placed in Standby and the power cable is unplugged from the AC outlet, the Display will remain powered for up to 60 minutes after which the system will go into Standby.

 NOTE: If the Intego PET Infusion System is in Standby and the power cable is not plugged into AC power within 30 minutes, the Intego PET Infusion System will completely Shutdown.

3. Disengage the Caster Brakes.

 NOTE: The Intego PET Infusion System cannot be moved using the Drive System if the Device Lock feature is enabled and the Display is locked.

NOTE: The Intego PET Infusion System cannot be moved using the Drive System unless the Shielded Chamber Lid is closed and the Shielded Chamber Lid Latch is engaged.

NOTE: The Drive System cannot be engaged when priming the SAS or PAS, performing a test injection, preparing an RP dose, or delivering an RP dose to the patient or Waste Container.

NOTE: The Intego PET Infusion System must be powered On or in Standby to use the Drive System. If the system is powered Off, the Drive System will not operate.

NOTE: The battery must be charged for the Drive System to function. The blue Status Indicator will be illuminated when the Drive System is available for use. If the blue Status Indicator is blinking, plug the Intego PET Infusion System power cable into an AC outlet to charge the battery.

4. Choose either slow or fast on the Drive Speed Switch.
5. Firmly grip both of the handles on the Intego PET Infusion System at the Operator Position. Squeeze the Drive Engage Switch.
6. To move forward, use the thumb of the right hand to slowly rotate the Drive Controller forward while keeping the Drive Engage Switch pressed. To move backward, use the thumb of the right hand to slowly rotate the Drive Controller backward while keeping the Drive Engage Switch pressed.

 NOTE: If the blue Status Indicator begins blinking while moving the Intego PET Infusion System, there are approximately 3 minutes remaining before the Drive System will stop functioning. When the blue Status Indicator begins blinking, move the system to an AC outlet and plug in the power cable.

7. Move the Intego PET Infusion System by carefully pushing the handles to the left and right to steer. The speed can be varied using the Drive Controller. The further the Drive Controller is rotated from the center position, the faster the system will move.

 NOTE: Use forward motion to cross thresholds while moving the Intego PET Infusion System.


8. Once the Intego PET Infusion System is in its new location, release the Drive Engage Switch or the Drive Controller to stop.
9. Apply the Caster Brakes.
10. Unwrap the power cable from the Cable Storage and plug into an AC outlet.

 NOTE: To ensure that the battery is fully charged and has the maximum capacity, keep the Intego PET Infusion System connected to AC power at all times except when the system is being moved.

3.12.1 Moving the System Using the Drive Override

NOTE: If the Drive System is not available, the Intego PET Infusion System can be moved by using the Drive Override.

1. Disconnect PAS from patient.
2. Unplug the power cable from the AC Outlet. Wrap the power cable on the Cable Storage
3. Pull the Drive Override out until it is fully extended (approximately 15 cm [6 inches]).

4. Disengage the Caster Brakes.
5. Firmly grip both handles on the Operator Position of the Intego PET Infusion System.
6. Move the Intego PET Infusion System using firm, steady force. Avoid making jerking movements or sudden starts and stops. Leave approximately 1 m (3 ft) per side to complete turns.
7. Once the Intego PET Infusion System is in its new location, apply the Caster Brakes.
8. Push the Drive Override in to the normal operating position.
9. Unwrap the power cable from the Cable Storage and plug into an AC outlet.


4 System Configuration

To configure the Intego PET Infusion System, press the CONFIGURATION button on the Navigation Bar.


NOTE: To save configuration changes, press the OK button. To return to the factory settings, press the DEFAULT button.

4.1 System Settings

On the left side of the Configuration screen, press the SYSTEM button to access the configuration options. Languages, Date/Time, Units, and Audio tabs will be displayed along the top.

1. Press the Languages tab to access the language configuration options.
 a. Configure the Intego PET Infusion System by pressing the button with the desired language.
2. Press the Date/Time tab to access the date and time configuration options.


 a. Configure the date format by pressing the Date Format field. Five date formats will be displayed. Press the button showing the desired date format.

b. Configure the time format by pressing the Time Format field. A 12-hour and a 24-hour format will be displayed. Press the button showing the desired time format.

c. Configure the current date and time by pressing the EDIT button. The current date is entered by pressing the Calendar icon. The current time is entered by pressing the Current Time field.

NOTE: The system does not automatically adjust for Daylight Savings Time. The system time must be adjusted manually.


3. Press the Units tab to access the weight and activity unit configuration options.

a. Configure the weight units by pressing either the POUNDS button or the KILOGRAMS button.

b. Configure the activity units by pressing either the CURIE button or the BECQUEREL button.

4. Press the Audio tab to access the audio configuration options.

a. Configure the audio level by pressing the HIGH button, NORMAL button, or LOW button.


4.2 Dosing Settings

⚠WARNING

Injury or death may result from injecting air.

- If an additional flush volume is used with extension tubing, the Intego PET Infusion System will not prime the additional tubing automatically. Air may be injected into the patient if the extension tubing is not primed before use. The additional tubing must be primed before connecting to the patient by selecting the “Prime PAS” function until all air is purged from the PAS.

On the left side of the Configuration screen, press the DOSING button to access the dosing options. RP, Saline, Case, and Printing tabs will be displayed along the top.

1. Press the RP tab to access the RP configuration options.

 a. Configure the default dose by pressing the Default Dose field and entering the desired value via the keypad.


 NOTE: If a vial schedule is entered using the Schedule screen, the Default Dose will be overridden with the schedule value. For more information on when the Default Dose is used and how a user may override it, refer to the "Patient Infusion" section.

 b. Configure the Intego PET Infusion System to use Personalized Dosing by default by pressing the Default Setting ON button.

 c. Configure the Intego PET Infusion System to use Manual Dosing by default by pressing the Default Setting OFF button.

 d. Up to five Personalized Dosing formulas can be entered by pressing the EDIT FORMULAS button.


 - Press the Multiplier field to enter the RP activity per unit weight and then press ENTER on the numeric key-pad.
 - Press the Minimum field to enter the minimum RP activity that can be prepared regardless of the patient weight and then press ENTER on the numeric keypad.
 - Press the Maximum field to enter the maximum RP activity that can be prepared regardless of the patient weight and then press ENTER on the numeric keypad.
 - Press the OK button.

NOTE: If more than one formula is created, a formula name must be entered when entering the data. Press a NEW FORMULA button and enter the formula name via the keypad.

NOTE: At least one formula must be entered to use the Personalized Dosing.

2. Press the Saline tab to access the saline configuration options.

a. Configure the default Saline Container volume by pressing the Default Volume field and entering the desired value via the keypad.


b. Configure the additional flush volume by pressing the Additional Flush Volume field and entering the desired value via the keypad.

NOTE: The additional flush volume is added to the fixed 35 ml saline flush used during a patient infusion.

c. Configure the test inject volume by pressing the Test Inject Volume field and entering the desired value via the keypad.

d. Configure the default flow rate by pressing either the 1.0mL/s or 0.5mL/s button.


3. Press the Case tab to access the Case configuration options.

a. Configure the Intego PET Infusion System by pressing the button with the desired default infusion site.

NOTE: If a default infusion site is not configured, the user will be given the option to select the infusion site when entering patient information on the Dosing screen. If a default infusion site is configured, it can be overridden when entering patient information on the Dosing screen.

4. Press the Printing tab to access the Printing configuration options.

a. Configure the Auto Print by pressing the ON or OFF button.

b. Press the Quantity field and enter the desired number of labels to automatically be printed if the Auto Print option is set to On.

NOTE: If the Auto Print option is not configured to be On to print labels, the user may still print manually.


4.3 Maintenance Settings

⚠WARNING

Shock hazard.

- The USB port is intended for connection of USB-powered devices, such as memory sticks, having no other electrical power connections. Connection of devices having other electrical power connections to this port may result in an electrical shock hazard.

1. On the left side of the Configuration screen, press the MAINTENANCE button to access the maintenance options. Daily QC, Linearity, Calibration, and Service tabs will be displayed along the top.
2. Press the Daily QC tab to access the Daily QC configuration options.

a. Press the Lot Number field and enter the lot number of the ^{137}Cs calibration source to be used for Daily QC.

b. Press the Activity field and enter the activity of the ^{137}Cs calibration source to be used for Daily QC.

c. Press the EDIT button to enter the assay date and time.


d. Press the Calendar icon next to the Reference Date and enter the assay date of the ^{137}Cs calibration source to be used for Daily QC.

NOTE: When reading the date from the ^{137}Cs calibration source assay information, be cognizant of the format. Some formats are dd/mm/yy and others are mm/dd/yy.

e. Press the Reference Time field and enter the assay time of the ^{137}Cs calibration source to be used for daily QC. If a time is not included with the assay information, enter 12:00.

NOTE: Pressing the Linearity and Calibration tabs will access these configuration screens. Descriptions on how to perform linearity testing, pinch valve, and yearly calibration are provided in the Cleaning and Maintenance Appendix.

3. Press the Service tab to access the service configuration options.


a. To configure a Service reminder, press the Calendar icon next to the Reminder Date field and enter the desired date.

b. To export log files, insert a USB memory device into the USB Port. Press either the 30 DAYS or ALL LOGS button, then remove the USB Memory device after the Intego PET Infusion System indicates the export is complete.

c. Bayer contact information is on the Service Configuration screen.


4.4 Security Settings

1. On the left side of the Configuration screen, press the SECURITY button to access the security options. A Manage Lock Code dialog box will appear.


2. Under Device Lock, press the ON button to enable the Device Lock feature. Press the OFF button to disable the Device Lock feature.

NOTE: Locking the device prevents unauthorized use or movement of the Intego PET Infusion System and requires a user to enter the default or a user-specific lock code after powering On.

3. To create a unique lock code, enter the current lock code by pressing the Current Code field. The factory set lock code is 7237. It can be used to override any user configured code. If the current user code is unknown, enter 7237 in the Current Code field.
4. Press the Enter New Code field. Enter a unique four digit lock code.
5. Press the Re-enter New Code field. Enter the same four digit lock code.
6. Press the ENTER button. The OK button will appear. Press OK to accept the new code.


4.5 Communication Settings

NOTE: Software features must be licensed in order to enable configuration options. Access to the configuration screens is available, but configuration options for unlicensed features are disabled.

NOTE: Users should not make changes to communication settings without the assistance of hospital IT staff.

On the left side of the Configuration screen, press the COMMUNICATIONS button to access the configuration options. Overview, Worklist, PACS and Reports will be displayed across the top.

1. Press the Overview tab to access the communications overview configuration options.
 a. Configure the AETitle by pressing the Intego AETitle field and entering the desired value on the keyboard.
 b. Configure the Institution by pressing the Institution field and entering the desired value on the keyboard.
 c. Configure the Department by pressing the Department field and entering the desired value on the keyboard.


2. Press the Worklist tab to access the worklist configuration options.
 a. Enable communication to the Worklist by pressing the ON button. Disable communication by pressing the OFF button.
 b. Configure the AETitle by pressing the AETitle field and entering the desired value on the keyboard.
 c. Configure the IP Address or Hostname by pressing either the IP Address or Hostname field and entering the desired value on the keyboard.
 d. Configure the Port by pressing the Port field and entering the desired value on the keyboard.
 e. Configure the Matching Station Name(s) by pressing the Matching Station Name(s) field and entering the desired value on the keyboard.
 f. Select a Modality Filter by pressing one of the Modality buttons or by pressing ALL for no filtering.

 NOTE: For Matching Station Name and AETitle, comma separation is required if adding multiple values.

g. Verify communication to the worklist by pressing the VERIFY button.

h. Configure the Matching AETitle(s) by pressing the Matching AETitle(s) field and entering the desired value on the keyboard.

3. Enable communication to the PACS by pressing the ON button. Disable communication by pressing the OFF button.

 a. Configure the PACS by pressing either the ON button or the OFF button.

 b. Verify communication to the PACS by pressing the VERIFY button.

 c. Configure the AETitle by pressing the AETitle field and entering the desired value on the keyboard.

 d. Configure the IP Address or Hostname by pressing the IP Address or Hostname field and entering the desired value on the keyboard.

 e. Configure the Port by pressing the Port field and entering the desired value on the keyboard.

4. Press the Reports tab to access the Reports configuration options.

 a. Configure the secondary capture reports that will be send to PACS by pressing the INFUSION RESULTS REPORT button, RADIATION DOSE REPORT button, or BOTH button.

5 Daily Setup

⚠WARNINGS
Patient injury could result from using improper components. ◆ Use only components and options provided by Bayer designed for Intego PET Infusion System.
Explosion hazard. DO NOT use the Intego PET Infusion System when flammable gases are present. ◆ Patient injury could result from using the system in the presence of flammable gases (such as anesthetics).
Trip hazard. ◆ Injury may occur. Do not operate Drive System if the power cable is plugged into an AC outlet.
No modification of this equipment is allowed.
Do not use the Intego PET Infusion System in an oxygen rich environment. ◆ Environment in which the concentration of oxygen is greater than 25% for ambient pressures up to 110kpa.


5.1 Placing the System into Clinical Mode

1. Power On the Intego PET Infusion System.


 NOTE: If the Device Lock feature has been enabled, the Intego PET Infusion System will prompt the user to enter the Lock Code. Press ENTER when complete.

2. The Introduction screen appears. Press the CLINICAL button


 NOTE: The screen displays the current mode in the lower right corner of the Introduction screen.

3. If the Intego PET Infusion System was in Training Mode, the system will display the message “Changing from Training Mode to Clinical Mode requires a restart. Continue?” Press OK to proceed.


4. The Intego PET Infusion System will display a progress bar while switching to Clinical Mode. When complete, the system will display a message prompting the user to restart. Press the RESTART button to Restart the system.

5. At the Introduction screen, press the CLINICAL button.

NOTE: Only screens for ^{18}F-FDG are shown here. When using ^{18}F-NaF, the text “FDG” changes to “NaF” and RP-specific screen elements change from orange to green (e.g., the RP button on the Navigation Bar).

6. When the System Preparation screen appears, press the OK button to proceed. Allow up to 15 minutes for the Dose Calibrator to warm up before proceeding; the Status Bar indicates the warm-up time remaining.


5.2 Using the Schedule

To access the Schedule screen, press the SCHEDULE button on the Navigation Bar. The schedule can either be entered manually or imported via the USB Port. The user may edit the schedule as it changes throughout the day.

Schedule appointments can also be imported as a modality worklist through the network interface if the worklist software feature has been licensed.

NOTE: The Intego PET Infusion System may highlight scheduled RP infusions. Refer to the "Monitoring the Vial" section of this document for more information.

5.2.1 Manually Entering Schedule Information


NOTE: The Patient ID, Patient Name, and Accession Number Fields are optional, but they are required to send Infusion Results or Radiation Dose reports to the PACS. If any of these fields are not specified, the Transmission Not Possible icon will be shown in the Schedule Monitor

1. Press the ADD APPOINTMENT button or a blank row in the schedule.
2. Press the Patient Name field to enter the patient name.

 NOTE: The Patient Name should be just the last name or entered as LAST NAME, FIRST NAME.

3. Press the Patient ID field to enter the patient ID.
4. Press the Time field to enter the planned infusion time.
5. Press the Requested Activity field to enter the planned dose activity.
6. Press the Accession Number field to enter the accession number.
7. Repeat steps 1 through 6 until schedule entry is complete.

 NOTE: More than five appointments can be added. Use the scroll bar to view appointments not on the Schedule screen.

5.2.2 Importing the Schedule Information via the USB Port

NOTE: Importing a schedule replaces the existing schedule.


1. Create an electronic version of the schedule named “schedule.csv” and copy it to a USB memory device. There are three ways to create an electronic schedule:
 a. Bayer’s MVP software tool (Recommended).
 b. Select third party nuclear medicine management software packages.
 c. Using a spreadsheet program or text editor.

 NOTE: Please contact Bayer for more information about creating an electronic schedule.

2. Insert the USB memory device into the USB Port on the Display.

3. Press the IMPORT DOSE PLAN button. The time, activity, and patient information automatically populate for all scheduled RP infusions

NOTE: Patient Names and Patient IDs that are imported from a dose plan may not be edited.


5.2.3 Importing the Schedule from the Modality Worklist

NOTE: This option is only available when users have licensed the worklist software feature.

NOTE: If the wireless connection is lost, the user will not be able to wirelessly import the worklist.

1. Press the IMPORT SCHEDULE button. The time, activity, and patient information automatically populate for all scheduled RP infusions.

NOTE: For worklist imports, all activity values are set to the system configured default activity.

5.2.4 Editing the Schedule

1. To edit a scheduled RP infusion, press the button with the patient name, ID, requested activity, and time to bring up the edit window.
 a. Press the Patient field to edit the Patient Name.

 NOTE: The Patient Name should be just the last name or entered as LAST NAME, FIRST NAME.

 b. Press the Patient ID field to edit the Patient ID.
 c. Press the Time field to edit the planned infusion time.
 d. Press the Requested Activity field to edit the planned dose activity.
 e. Press the Accession Number field to edit the accession number.

 NOTE: The Intego PET Infusion System automatically sorts the schedule by infusion time. Editing infusion time may reorder the schedule.


 NOTE: If users have licensed the software features and imported the schedule from a worklist, the Patient, Patient ID, and Accession Number fields are not editable.

2. To hide a scheduled infusion, press the HIDE/UNHIDE APPOINTMENT button. The button will then display the HIDDEN APPOINTMENT icon to indicate that the appointment has been hidden on the Dosing screens.
3. To unhide a hidden appointment, press the HIDE/UNHIDE APPOINTMENT button again. The button will then display the UNHIDDEN APPOINTMENT icon to indicate that the appointment is not hidden on the Dosing screens.

 NOTE: When a scheduled infusion is hidden or unhidden, the number following the HIDDEN APPOINTMENT icon at the top of the list updates to show the current number of hidden appointments on the schedule.

 NOTE: Hiding a scheduled infusion will not remove the infusion from the Schedule screens. It will only remove the infusion from the Dosing screens. A hidden infusion is signified on the Schedule screen by the HIDDEN APPOINTMENT icon on the HIDE/UNHIDE APPOINTMENT button.


4. To remove the entire schedule, press the CLEAR SCHEDULE button.

5.2.5 Exporting the Infusion History

The Intego PET Infusion System creates an infusion history as each infusion is completed. The infusion history consists of a .csv file for each of the last 31 days the Intego PET Infusion System was used.

1. To export the infusion history, insert a USB memory device in the USB Port and press the EXPORT HISTORY button.
2. Remove the USB memory device when the system indicates the export has been completed.

5.3 Performing Daily QC

⚠WARNING


Using the Intego PET Infusion System without performing the necessary QC processes could result in the patient receiving an incorrect dosage of the RP indicated for use.

- Be certain to perform all necessary QC checks at the recommended intervals.

The Daily QC confirms proper operation of the Dose Calibrator. Daily QC can only be performed without a Multi-Dose Vial or SAS installed. It is possible to perform RP infusions without valid Daily QC results, but this is not recommended. An Attention icon is displayed on the DAILY QC button when the system is being used without valid Daily QC results.


1. Access the System Preparation screen by pressing the SALINE button or RP button on the Navigation Bar.

2. Press the DAILY QC button.

NOTE: Nearby radiation sources may cause the zero adjustment to fail. Remove any nearby radiation sources before starting daily QC.


3. Press the START button. The Intego PET Infusion System will automatically proceed with all listed Daily QC tests and check the box next to each test as it completes.


NOTE: The Bias Adjustment may compensate for a failed Zero Adjustment. The Daily QC is valid if the Zero Adjustment fails but the Bias Adjustment is successful.

4. If the Zero Adjustment fails, check for and remove any nearby sources of radiation and repeat the Daily QC test from Step 3.
5. Daily QC will pause at the Constancy/Accuracy test and produce the prompt shown below.
 a. If the ^{137}Cs calibration source is not available, press the SKIP button to bypass the Constancy/Accuracy test and proceed to step 10.
 b. If the ^{137}Cs calibration source is available but the displayed assay information is incorrect, press the EDIT button, update the assay information, and proceed to step 6.

NOTE: When the ^{137}Cs calibration source lot number is changed, the Constancy check will establish a new reference activity for future Constancy checks instead of testing Constancy against the previous reference.


6. Grasp the Shielded Chamber Lid Latch to release it, and open the Shielded Chamber Lid to gain access to the Shielded Chamber.
7. Place the ^{137}Cs calibration source into the Dose Calibrator using the Calibration Source Holder.

8. Close the Shielded Chamber Lid until the Shielded Chamber Lid Latch engages, then press the OK button.
9. When prompted, open the Shielded Chamber Lid, remove the test source from the Dose Calibrator, slide the Shielded Chamber Lid closed until the Shielded Chamber Lid Latch engages, then press the OK button.
10. At the conclusion of the Daily QC, the Intego PET Infusion System displays a summary screen. Print this screen to keep with the system's records.

11. Press the SUMMARY button to return to the main System Preparation screen.

5.4 Entering RP Assay and Saline Information


⚠WARNING

Radiation exposure hazard.


- Do not insert more than 27.75 GBq (750 mCi) into the Intego PET Infusion System.

NOTE: Entering RP Assay and Saline Information may be performed before or after Installing the SAS or Installing the Multi-Dose Vial.

1. Access the System Preparation screen by pressing the SALINE or RP buttons on the Navigation Bar.
2. If a Multi-Dose Vial is loaded in the Intego PET Infusion System, remove the Vial Shield and SAS, then press the REMOVE button.
3. Press the NEW button to enter RP Assay information.


4. The Vial Assay Information screen appears. Press the respective fields to enter the Lot Number, Date, Reference Time, RP Type, Activity at Reference Time, and Volume. When all data is entered, press the OK button.


NOTE: If the level of radioactivity in the Multi-Dose Vial exceeds 27.75 GBq (750 mCi), the Intego PET Infusion System displays the message "Vial activity exceeds safe shielding level. Immediately return vial to protective container." If this message appears, remove the Vial Shield from the system and return to hot lab or other storage capable of providing adequate shielding for the contained activity.

5. The System Preparation screen reappears with the usable activity in the Multi-Dose Vial displayed within the vial symbol. The total Multi-Dose Vial volume and activity are displayed to the right of the vial symbol.

6. Press the REPLACE button in the Saline section.

7. Enter the volume of the Saline Container using the keypad.

 NOTE: The minimum recommended Saline Container volume is 750 ml.

5.5 Installing the SAS

⚠WARNINGS

Do not use the SAS with more than one multi-dose vial.

- Using the SAS with more than one multi-dose vial could result in injury to the patient or operator.

Severe operator injury can occur if care is not taken when opening or closing the Shielded Chamber Lid or the PAS Compartment Cover.

- Be certain that all fingers are clear before opening or closing the Shielded Chamber Lid or the PAS Compartment Cover.

Do not use if sterile package is opened or damaged.

- Patient or operator injury may result if package is opened or damaged, or if damaged components are used. Visually inspect contents and package before each use.

Patient Infection may result from the use of non-sterile components.

- Maintain sterility of all disposable components using aseptic techniques.

Only the fluid path of the set is sterile and non-pyrogenic.

- Do not use in a sterile or aseptic area without proper precautions.

Puncture hazard. The SAS contains needles.

- To prevent injury, take care when handling the SAS to control the needles. Dispose of the needles per individual facility guidelines.

Improper injection dosage or radiation spillage may result if the Shielded Chamber Lid is closed onto the tubing set.

- Ensure that the tubing set and all administration set components are properly installed before closing the Shielded Chamber Lid.

⚠CAUTION

Incorrect installation of the SAS may result in improper functioning of the Intego PET Infusion System.


NOTE: Always install a new SAS before installing a new Multi-Dose Vial.

NOTE: The Intego PET Infusion System is intended for use with ^{18}F-FDG and ^{18}F-NaF. Intego PET Infusion System disposables intended for use with ^{18}F-FDG may be used to deliver ^{18}F-FDG or ^{18}F-NaF. Please contact Bayer for more information.


NOTE: Installing the SAS or installing the Multi-Dose Vial may be performed before or after entering RP Assay and Saline Information.

NOTE: In the following steps, it will be helpful to refer to the diagrams of the Intego Infusion System Components, Intego Shielded Chamber Components, and the SAS Components.


1. Place the Saline Container onto the Saline Hook at the rear of the cart.
2. Open the PAS Compartment Cover.
3. Grasp the Shielded Chamber Lid Latch to release it, and slide the Shielded Chamber Lid towards the front of the cart to gain access to the Shielded Chamber.
4. Open the SAS Disposable Package.
5. Insert the Saline Spike into the Saline Container. Verify the Saline Spike is fully inserted into the container as illustrated below.


6. Open the Saline Pump door, install the Saline Tube, then close the Saline Pump door.
 a. Center the Saline Tube on the Saline Pump rollers. The Saline Pump will not pump fluid properly if the tubing falls below the rollers.

 b. Ensure the Saline Disconnect Luer rests against the edge of the Saline Pump.

 c. Center the Saline Tube in the Saline Pump door to prevent the tubing from being pinched when the door is closed.


7. Route the Saline Tube along the SAS Track, starting at the SAS Confluence Holder and moving towards the Saline Pump. Be sure to press the tubing firmly into the SAS Track.

8. Press the Saline Tube into the Saline Tube Holder between the Saline Pump and SAS Track.

9. Load the Needle Cartridge into the Vial Shield Compartment:
 a. Grasp the Needle Insertion Device's handle and lift the Needle Insertion Device straight up until it stops.
 b. Rotate the Needle Insertion Device counter-clockwise until it stops. Release the handle. The Needle Insertion Device will remain in position.
 c. Lift the hinged Needle Insertion Device cover to access the Needle Cartridge Holder.

d. Insert the Needle Cartridge into the holder so the RP Tube exits the cartridge away from the holder.

e. Stow the needles safely in the Vial Compartment until ready to install the Multi-Dose Vial by lowering the Needle Insertion Device cover, rotating the Needle Insertion Device clockwise until it stops, then pushing down on the Needle Insertion Device until it is back in its original position.

10. Insert the Confluence into the SAS Confluence Holder.

11. **Verify that the dose calibrator is empty before installing the coil.**

12. Insert the SAS Coil into the Dose Calibrator. Be sure the Pre-Coil Tube enters and the Post-Coil Tube exits on top of the SAS Coil and the SAS Coil is placed at the bottom of the Dose Calibrator.

13. Route the Pre-Coil Tube into the SAS track between the SAS confluence holder and the Dose Calibrator. Be sure to press the tubing firmly into the SAS track.


14. Unroll the Waste Container and insert it into the Waste Storage.

15. Place the T-Connector into the T-Connector Holder.


16. Place Waste Tube into the SAS track between the Waste Pinch Valve and Waste Storage. Be sure to press the tubing firmly into the SAS track.

17. Install the Waste Tube into the Waste Pinch Valve. Install the Patient Tube into the PAS Pinch Valve.


NOTE: To install tubing in a Pinch Valve, open the valve by pressing the button on the top. While holding the valve open and grasping the tubing on either side of the valve, pull the tubing in from the side. The tubing must be inserted completely as shown in the left panel below. Release the button on the valve to pinch the tubing in place.

18. Install the Waste Pinch Valve cover (if available) on the Waste Pinch Valve. Install the PAS Pinch Valve cover (if available) on the PAS Pinch Valve.


NOTE: The left panel in the picture below shows the correct layout of the SAS when the tubing has been properly installed into the SAS T-Connector Holder, Waste Pinch Valve, and PAS Pinch Valve.


19. Firmly press the Patient Tube exiting the PAS Pinch Valve into the Air Detector Holder and Air Detector. Insert the Swabbable Valve into the Swabbable Valve Holder.

NOTE: Make certain the tubing is fully inserted into the Air Detector as illustrated below.

20. Insert RP Tube in the RP Pump:
 a. Open the RP Pump door.
 b. Insert the RP Tube into the RP Pump as shown in the left panel below. For proper installation, the tubing collar must rest immediately outside the RP Pump on the side closest to the Vial Shield Compartment.

 c. Close the RP Pump door.
21. Once the SAS has been installed, ensure all components are laying flat prior to closing the lid.

5.6 Installing the Multi-Dose Vial

⚠WARNINGS

DO NOT use any RP other than ^{18}F-FDG or ^{18}F-NaF with the Intego PET Infusion System.

- Use of radiopharmaceuticals other than ^{18}F-FDG or ^{18}F-NaF may cause inaccurate dose measurement and/or excessive radioactive exposure..

Radiation exposure hazard. ^{18}F-FDG or ^{18}F-NaF emit gamma radiation. Operator should exercise extreme caution when handling radiopharmaceuticals.

- Doses may exceed maximum permissible dose if appropriate precautions are not properly exercised. Follow all individual facility guidelines for the handling and disposal of radioactive materials. DO NOT open the Shielded Chamber without first taking radiation precautions per facility procedures.

Severe operator injury can occur if care is not taken when opening or closing the Shielded Chamber Lid or the PAS Compartment Cover.

- Be certain that all fingers are clear before opening or closing the Shielded Chamber Lid or the PAS Compartment Cover.

Radiation exposure hazard. Using the Shielded Chamber Lid as a work space may expose the operator to excess radiation.

- Do not use the Shielded Chamber Lid as a work space while there are radioactive materials present within the Intego PET Infusion System. Using the Shielded Chamber Lid as a work space or placing items on the Shielded Chamber Lid while it is open may damage the Shielded Chamber Lid, hampering its ability to provide an adequate radiation shield.

Radiation exposure hazard.

- Do not insert more than 27.75 GBq (750 mCi) into the Intego PET Infusion System.

Puncture hazard. The SAS contains needles.

- To prevent injury, take care when handling the SAS to control the needles. Dispose of the needles per individual facility guidelines.

Radiation exposure hazard.

- Radiation activity exposure is increased once the Access Cap is removed. The user should avoid working directly over the open Vial Shield to minimize exposure.
- Once RP priming commences, RP is in the SAS and no longer contained within the Vial Shield. To minimize radiation exposure, opening the Shielded Chamber Lid should be avoided once RP priming commences.

Biological contamination may occur if aseptic techniques are not followed.

- Follow facility guidelines to disinfect the vial septum.


NOTE: Entering RP Assay and Saline Information may be completed before or after installing the Multi-Dose Vial. Installing the SAS should be completed prior to installing the Multi-Dose Vial.

NOTE: In the following steps, it will be helpful to refer to the diagrams of the Intego Infusion System Components, Intego Shielded Chamber Components, and the SAS Components.

1. Grasp the Needle Insertion Device's handle, lift straight up until it stops, then rotate counter-clockwise until it stops. Release the handle; the Needle Insertion Device will remain in position providing access to the Vial Shield Compartment.
2. Insert the Vial Shield containing a Multi-Dose Vial of RP into the Vial Shield Compartment. DO NOT remove the Vial Shield's Access Cap. Leave the carrying handle on the Vial Shield until ready to insert the SAS Needles into the vial.
3. When ready to insert the SAS Needles into the Multi-Dose Vial, remove the Access Cap: grasp the Vial Shield Carrying Handle and simultaneously push down and rotate counter-clockwise 1/8 of a turn. After completing the turn, lift the handle straight up. The Access Cap will lift from the Vial Shield, providing access to the Multi-Dose Vial septum.

4. Insert SAS Needles into the Multi-Dose Vial:
 a. Grasp the Needle Insertion Device's handle and rotate it clockwise until it stops; the SAS Needles should be oriented over the Multi-Dose Vial septum.
 b. Press down on the Needle Insertion Device's handle until the Needle Insertion Device is in its original closed position as shown in the Intego Shielded Chamber Components diagram. The SAS Needles will penetrate the Multi-Dose Vial as the Needle Insertion Device is lowered.
 c. Lower the Needle Insertion Device handle to rest on the cover.
5. Close the Shielded Chamber Lid until the Shielded Chamber Lid Latch engages and close the PAS Compartment Cover over the PAS Compartment.

 The Vial Shield is keyed such that it can be properly inserted into the Vial Shield Compartment only one way. The posts in the Vial Shield Compartment must rest in the keyed areas located on one flat surface of the Vial Shield.

5.7 Installing the PAS

⚠WARNINGS

Biological contamination can result from reusing the PAS or failure to follow sterile technique.

- DO NOT reuse or attempt to resterilize. Properly discard disposable items after use. If there is any possibility that contamination may have occurred during set-up or use, disassemble and set-up a new sterile product.

This product contains DEHP, a chemical known to the state of California to cause birth defects or other reproductive harm.

- Please visit www.medrad.com/DEHP for more information regarding DEHP in Bayer products.

Do not use if sterile package is opened or damaged.

- Patient or operator injury may result if package is opened or damaged, or if damaged components are used. Visually inspect contents and package before each use.

For devices labeled for single use, please note: This product is intended for single use only. Do not resterilize, reprocess or reuse. The disposable devices have been designed and validated for single use only.

- Re-use of the single use disposable devices pose risks of device failure and risks to the patient. Potential device failure includes significant component deterioration with extended use, component malfunction, and system failure. Potential risks to the patient include injury due to device malfunction or infection as the device has not been validated to be cleaned or re-sterilized.

Patient infection may result from the use of non-sterile components.

- Maintain sterility of all disposable components using aseptic techniques.

Component damage may occur if not installed properly.

- Ensure all connections are secure. Do not over-tighten.

⚠WARNINGS

Severe operator injury can occur if care is not taken when opening or closing the Shielded Chamber Lid or the PAS Compartment Cover.

- Be certain that all fingers are clear before opening or closing the Shielded Chamber Lid or the PAS Compartment Cover.

NOTE: A PAS must be installed to prime the SAS and to infuse a patient. A new PAS must be installed for each patient. The PAS installed to prime the SAS may be used for the first patient after SAS priming.

NOTE: In the following steps, it will be helpful to refer to the diagrams of the Intego Infusion System Components, Intego Shielded Chamber Components, and the SAS Components.

1. Open the PAS Compartment Cover to expose the Swabbable Valve.
2. Clean the SAS Swabbable Valve following aseptic techniques described in facility guidelines.
3. Open the PAS disposable package, being sure not to remove the Removable Prime Tube.
4. Attach the PAS Luer Connector to the SAS Swabbable Valve.
5. Route the PAS from the PAS Compartment.


 NOTE: Route the PAS out the left, right, or rear of the PAS Compartment, whichever is most convenient for connecting the PAS to the patient.

6. Close the PAS Compartment Cover back over the PAS Compartment. Make sure the PAS Compartment Cover does not pinch the PAS.

5.8 Priming the SAS

NOTE: Entering RP Assay and Saline Information, Installing the SAS, Installing the Multi-Dose Vial, and Installing a PAS must be completed prior to Priming the SAS. To prime the SAS:

1. Access the System Preparation screen by pressing the SALINE or RP buttons on the Navigation Bar.
2. Press the PRIME SAS button.

3. A dialog box appears with the following message: “Ensure that a PAS is connected to the SAS. Continue?” Install a PAS if necessary, then press the YES button.
4. SAS priming proceeds in three phases.

a. Priming Saline - The Intego PET Infusion System fills the Saline Tube, Pre-Coil Tube, SAS Coil, Post-Coil Tube, Waste Tube, and Patient Tube with saline. All RP is still contained within the Vial Shield at this point, allowing safe correction of SAS installation errors detected during this phase.

b. Priming RP - The Intego PET Infusion System fills the SAS Needles and RP Tube with RP.

c. Verifying Concentration - The Intego PET Infusion System measures RP activity concentration to ensure RP assay information has been entered accurately.

NOTE: During each phase of SAS priming, the Status Bar indicates the phase being performed and displays a progress bar showing progress to completion.

5. When SAS Priming is complete, the Intego PET Infusion System will beep and the Status Bar will state “System is Ready.”

NOTE: The attached PAS is not primed during this process; it is primed when preparing for the first patient infusion. Excess fluid used during SAS priming is collected within the Waste Container inside the Shielded Chamber.

NOTE: To resolve priming issues, refer to the "Troubleshooting Tips" section of this manual.

6 Patient Infusion

⚠WARNING

Explosion hazard. DO NOT use the Intego PET Infusion System when flammable gases are present.


- Patient injury could result from using the system in the presence of flammable gases (such as anesthetics).

6.1 Adding a Patient / Editing Patient Information

If a schedule has not been loaded into the system, the user may add a patient to the system from the Dosing screen. This adds that patient to the schedule and automatically selects that patient for dosing. If a patient has already been selected, the user may edit the selected patient's information until a dose has been extracted. To add a patient or edit patient information, do the following:

1. Press the DOSING button on the Navigation Bar to go to the Dosing screen.
2. To enter patient information, press the ADD button located on the Patient Information panel.


 NOTE: If a patient has not been selected, the Patient Information panel will show an ADD button. Press the ADD button to add a patient and have them automatically selected. If a patient has already been selected, the Patient Information panel will show an EDIT button. Press the EDIT button to edit the patient information for the selected patient.

3. A patient information entry screen will appear.

 NOTE: The patient information entry screen may appear with some or all patient information populated with valid data. Pre-populated patient information can be provided with an imported Multi-Dose Vial schedule or if retrying a patient for which patient information was already entered. Pre-populating patient information is a convenience feature to reduce the amount of hand-entry a user must perform.

 NOTE: The Intego PET Infusion System will prevent initiation of a patient infusion with less than three minutes remaining on the battery meter.

4. Enter patient and operator information and update any information that is incorrect from the patient information entry screen.

 a. To enter a Patient Name, Date of Birth, Patient ID, Accession Number, or Operator ID, touch the corresponding data field and an on-screen keyboard will appear.

 b. To enter a Gender, press the appropriate FEMALE or MALE button.

 c. To enter an Infusion Site, press the button associated with the desired Infusion Site. If "Other" is entered for the Infusion Site, the operator can add additional detail in the Detail field by touching the data field and using the on-screen keyboard.

 d. To enter a Blood Glucose, press the corresponding data field and an on-screen numeric keypad will appear.

 e. To enter a Flow Rate, press the corresponding data field and press the desired 1.0 ML/SEC or the 0.5 ML/SEC button.

 f. Press OK when finished entering patient information.

 NOTE: The Intego PET Infusion System requires confirmation that the Patient Entry field has been reviewed, even if information was not entered. To proceed without entering information, press ADD, then press OK.

6.2 Selecting a Patient from the Schedule

If a schedule has been added by manually adding appointments, importing a dose plan, or importing a worklist through the worklist software, the user should select the patient to dose from the schedule as follows:

1. Press the DOSING button on the Navigation Bar to go to the Dosing screen.
2. Press the SELECT button located on the Patient Information panel. A patient schedule will appear.
3. Search for the desired patient in the list. If the patient appears in the list, go to step 6.
4. If the patient does not appear in the list, press the REFRESH button to update the list from the worklist. If the patient appears in the list, go to step 6.

5. If after pressing the REFRESH button the desired patient still doesn't appear in the list, press the ADD PATIENT button. Proceed to step 3 in the "Adding a Patient/Editing Patient Information" on page 71.
6. Press the PATIENT button to select the desired patient. A patient information screen will appear.

 NOTE: Patient and operator information is used for dose infusion record-keeping. All patient and operator information is optional. However, to send images to the PACS, the Patient Name, Patient ID, and Accession Number must be defined. Fields that are not entered are left blank in the final dose infusion record.

 NOTE: Patient information fields will be populated with information that was manually entered via the schedule screen, imported from a dose plan, or imported from the worklist if the user has licensed the worklist software feature.

 NOTE: Patient Name, Patient ID, Date of Birth, Gender, and Accession Number fields will not be editable if users have imported the schedule from the worklist or dose plan.
7. Users can enter any of the remaining patient information from this screen by following the steps in "Selecting a Patient from the Schedule".
8. Press OK. The selected patient is highlighted in the Schedule Monitor.

6.3 Priming the PAS

⚠WARNINGS

Air embolism hazard. Injury or death can result from infusing air.

- Ensure air is expelled from the disposable sets before beginning infusion.

Biological contamination can result from reusing the PAS or failure to follow sterile technique.

- DO NOT reuse or attempt to resterilize. Properly discard disposable items after use. If there is any possibility that contamination may have occurred during set-up or use, disassemble and set-up a new sterile product.

Patient infection may result from the use of non-sterile components.

- Maintain sterility of all disposable components using aseptic techniques.

Radiation exposure hazard. Higher residual radiation may occur when using an extension tube or additional tubing between the catheter and the PAS.

- The flush volume needs to be increased to account for additional tubing.

Injury to patient may occur if system is moved without disconnecting the PAS from the patient.

- Disconnect the PAS from the patient prior to moving the Intego PET Infusion System.

The system is not to be serviced or maintained while in use with a patient.

⚠CAUTION

Catheters no more restrictive than 24 gauge x 0.75 inches (19 mm) long, and steel needles no more restrictive than 23 gauge x 0.5 inches (12 mm) long, should be used with the Intego PET Infusion System.

- The system performance may not be achievable if more restrictive catheters or needles are used.

Only use the Intego PET Infusion System PAS to connect to the patient catheter or needle.

- The system performance may not be achievable if tubing extension sets, additional tubing, or tubing other than an Intego PET Infusion System PAS is used.

1. Install a PAS, if not already installed.

2. Hold the Removable Prime Tube on the PAS vertically so that the filter (white tip) is pointed upward.
3. Press the PRIME button.

 NOTE: The Intego PET Infusion System will prevent priming with less than three minutes remaining on the battery meter.

4. When priming is complete, the Intego PET Infusion System will beep. Confirm that fluid is in the Removable Prime Tube and visually inspect the PAS for air.

 NOTE: If air is present within the tubing set after PAS priming is complete, evaluate the risk to the patient and reprime as needed by pressing the PRIME button. Repriming the PAS may overflow the Removable Prime Tube. Before repriming the PAS, remove the Removable Prime Tube and place the end of the PAS in an appropriate container for the excess fluid. If such a reservoir is not available, replace the PAS with a new one. A dialog box is displayed stating that PAS priming is already complete and asks whether the clinician wishes to prime again. Press the YES button.

5. Remove the Removable Prime Tube from the end of the PAS and discard.
6. Connect the PAS to the patient.

6.4 Selecting the Flow Rate

1. The currently configured flow rate is displayed. If this flow rate is not the desired flow rate, press the button and select the desired flow rate from the displayed list.

2. The selected flow rate will be displayed.

6.5 Performing a Saline Test Injection

⚠WARNING


Patient injury could result from not following standard extravasation minimizing techniques.

- To reduce the risks of extravasation conditions, follow all standard extravasation risk minimizing techniques.


NOTE: The Test Inject is optional.

NOTE: The Intego PET Infusion System will prevent a Saline Test Inject with less than three minutes remaining on the battery meter.

1. On the Dosing screen, press the TEST INJECT button.

2. The Intego PET Infusion System displays a message asking the user to confirm that a check for air has been completed. Inspect the PAS Tubing for air.
 a. If the PAS contains air, press NO. Disconnect the PAS from the patient and return to the "Priming the PAS" section of the manual.
 b. If the PAS does not contain air, press YES. The Intego PET Infusion System will infuse the patient with a small, preset volume of saline. See "System Configuration".

6.6 Entering the Requested Dose Activity

The requested patient dose is displayed in the Requested Activity field. The requested patient dose may be entered by using a personalized dosing formula or manually.


6.6.1 Personalized Dose Entry

To select Personalized Dose Entry, press the dosing method button until it reads "Dosing: Personal." The requested patient dose is then established as follows:

1. The FORMULA button displays the currently selected formula. If this formula is not the desired formula, press the button and select the desired formula from the displayed list.

 NOTE: See "System Configuration" to create formulas.


2. Press the Patient Weight field, enter the patient's weight on the keypad, and then press ENTER.
3. The Intego PET Infusion System automatically calculates the requested patient dose using the entered weight and selected formula.

6.6.2 Manual Dose Entry

To select Manual Dose Entry, press the dosing method button until it reads "Dosing: Manual." The requested patient dose is established via one of the following:

- Default Dose: When there are no scheduled RP infusions, the requested patient dose is automatically set to the default value when initiating dosing for a new patient or when the RESET button is pressed.
- Scheduled Dose: If there are scheduled RP infusions, the requested patient dose is automatically set to the current scheduled patient dose when initiating dosing for that patient or when the RESET button is pressed.
- Direct Entry: Direct entry may be used to override either a default patient dose or a scheduled patient dose. Press the Requested Activity field and use the keypad to enter the activity value. The Intego PET Infusion System will only accept a value within the displayed Valid Range (the range of patient doses the system is capable of delivering at the time).


6.6.3 Requested Activity and the Activity Bar

Once the requested patient dose is established, it is displayed graphically on the Activity Bar below the Requested Activity field. The Activity Bar shows the requested patient dose, valid dosing range, and personalized dosing limits defined by the user selected formula within the total patient dosing range of the Intego PET Infusion System.

NOTE: The filled color portion covers the valid range of doses the Intego PET Infusion System can deliver under current conditions.

NOTE: The white dashed lines and pointers mark the upper and lower dosing limits imposed by the currently selected formula (only visible when using Personalized Dose Entry).

NOTE: If the requested dose activity is within the valid range, it is marked with a black line and pointer.

NOTE: If the requested dose activity is not within the valid range, it is marked with a yellow line and pointer. The valid range is also displayed within the Requested Activity field under the requested dose activity value.

6.7 Preparing an RP Dose


⚠WARNING

Radiation exposure hazard. ^{18}F-FDG or ^{18}F-NaF emit gamma radiation. Operator should exercise extreme caution when handling radiopharmaceuticals.

- Doses may exceed maximum permissible dose if appropriate precautions are not properly exercised. Follow all individual facility guidelines for the handling and disposal of radioactive materials. DO NOT open the Shielded Chamber without first taking radiation precautions per facility procedures.

NOTE: When operating on battery power, the Intego PET Infusion System will prevent Preparing an RP Dose with less than three minutes remaining on the battery meter.

1. Press the PREPARE button to prepare the dose of RP to be infused.
2. If a test injection was not already performed, the Intego PET Infusion System displays a message asking the user to confirm that the PAS has been inspected for air (see the "Performing a Saline Test Injection" section). Inspect the PAS Tubing for air:
 a. If the PAS contains air, press NO, disconnect the PAS from the patient, and return to the "Priming the PAS" section of the manual.
 b. If the PAS does not contain air, press YES. The Intego PET Infusion System will prepare the requested dose.


NOTE: Once a patient dose is prepared, the Requested Activity field displays the patient dose's actual activity as measured by the Dose Calibrator. The Activity Bar changes to show the percent deviation of the measured dose with respect to the requested dose on a scale of +/- 10%. If the measured dose is beyond +/- 10% of the requested dose, the Intego PET Infusion System notifies the user by displaying the triangle at the end of the Activity Bar as color-filled. The Requested Activity field and Activity Bar continuously update as the prepared dose decays.

6.8 Infusing or Discarding the Dose

After the dose is prepared, the user may infuse the patient with the dose or discard the dose to the Waste Container.

If the wireless feature is enabled, the following transmission status icons will appear on the dosing screen to alert the user of the status of the wireless transmission:

Transmission in progress icon. Used to indicate secondary captures for the patient are in transmission to PACS.

Transmission successful icon. Used to indicate secondary captures for the patient have been successfully transmitted to PACS.

Transmission error icon. Used to indicate an error in the transmission of secondary captures to PACS for the patient.
The system will retry information transmission three times. There is a 30 second delay between tries. If the system cannot upload the data after three attempts, this icon will display.

Transmission not possible icon. Used to indicate transmission of secondary captures to PACS for the patient is currently not possible.
Transferring to PACS will only be attempted if Patient Name, Patient ID, and Accession Number are entered. If one of these is not entered, this icon will display.


6.8.1 Start and Monitor Infusion

⚠WARNING


In case of occlusion, do not disconnect disposables before discarding dose, as liquid leakage can occur.

NOTE: When operating on battery power, the Intego PET Infusion System will prevent Infusing or Discarding the Dose with less than three minutes remaining on the battery meter.

1. To infuse the dose, press the INFUSE button. To discard the dose, press the DISCARD button.

2. While infusing or discarding the dose, the Intego PET Infusion System displays a graph showing the percentage of the prepared dose activity remaining in the Dose Calibrator over the infusion or discard. The graph indicates whether activity is moving as expected or if there is a possible occlusion restricting the activity being delivered to the patient or Waste Container.

NOTE: The Intego PET Infusion System may detect an occluded line to the patient during an infusion. Appendix C, PAS Recovery, provides instructions for handling occlusions detected during an infusion.

6.8.2 Infusion Completion


After an infusion completes, the Intego PET Infusion System displays a summary screen that shows the amount of RP delivered (both volume and activity), total fluid volume delivered, and any patient information provided.

1. Press the PRINT button to print this summary.
2. Press the OK button to update the infusion history and close the summary window.


NOTE: The Intego PET Infusion System should not be put in Standby before pressing the OK button.

NOTE: The printed label does not have all the information included on the Display.

NOTE: After the infusion completes, disconnect the PAS from the patient and Intego PET Infusion System. Although PAS residuals are negligible for infusions without occlusions, the PAS should be disposed of using the facility procedure for disposal of radioactive material.

NOTE: After pressing the OK button, the main Dosing screen appears.

NOTE: The Schedule Monitor will replace the last RP infusion schedule information with the actual infusion results and marked with a Complete Infusion Delivered icon to indicate a completed infusion.

NOTE: The Schedule Monitor on the left side of the screen will have moved its highlight to the next RP infusion in the Multi-Dose Vial schedule.

NOTE: On the right side, only the patient information panel is visible. To start the next patient, install a PAS, then proceed with "Adding a Patient / Editing Patient Information".

6.8.3 Discard Completion

After a discard completes, the Intego PET Infusion System displays a summary screen that shows the RP activity delivered, RP activity discarded, and any patient information provided.

1. To retry infusing the same patient, press the RETRY button. The main Dosing screen appears with the patient information panel visible on the right side. To start dosing the patient, install a PAS, then proceed with entering the patient information (all information will be pre-populated with previously entered patient information).
2. If the patient must be rescheduled or removed from the Multi-Dose Vial schedule, press the RESCHEDULE button. The Schedule screen appears, from which the scheduled infusion time can be changed or the RP infusion can be deleted.
 a. To delete a scheduled infusion, press the X button.
 b. To edit a scheduled RP infusion, press the Time field or Activity field and enter the new value.
3. To start the next patient, install a PAS, go to the Dosing screen, and proceed with ""Adding a Patient / Editing Patient Information""

6.9 Sending Secondary Capture Images to PACS

NOTE: This option is only available if users have licensed the PACS software features.


1. If users have licensed the Infusions Results Report and have enabled the report in the configuration settings, a secondary capture of the infusion results will be transmitted to PACS after each patient infusion.
2. If users have licensed the Radiation Dose Report and have enabled the report in the configuration settings, a secondary capture of the radiation dose information will be transmitted to PACS after each patient infusion. Radiation dose report provides estimated radiation dose by organ, as well as the estimated total effective dose calculated using NUREG/CR-6345.

6.10 Monitoring the Vial


The Intego PET Infusion System continuously analyzes if the activity and volume of the currently loaded Multi-Dose Vial is sufficient to meet the dosing requirements of the Multi-Dose Vial schedule. When the Multi-Dose Vial schedule cannot be satisfied with the Multi-Dose Vial, the system alerts the user as follows:

NOTE: Hidden items are not considered in these calculations. Users should hide any patients they do not want to be considered in the vial monitoring.

1. On the Dosing screen: By displaying an Attention icon and message in the Status Bar, by displaying an Attention icon at the top of the Schedule Monitor, and by highlighting the unachievable RP infusions in yellow.

2. On the Schedule screen: By displaying an Attention icon and message in the Status Bar, by highlighting the unachievable RP infusions in yellow, and by displaying the valid dose range under the planned dose for the first unachievable RP infusions in the schedule.

While the Multi-Dose Vial is intended to meet the original Multi-Dose Vial schedule, there are several reasons a schedule may become unachievable:

- Scheduled infusions were delayed.
- More activity is used per RP infusion than originally planned.
- Unplanned RP infusions are added during the day.

Users may temporarily enter schedule changes under consideration to determine if sufficient RP activity is in the Multi-Dose Vial to support the new plan. Users should only commit to schedule changes if the Vial Monitor does not raise an alert. If changes that create an unachievable schedule cannot be avoided, the Vial Monitoring early warning alerts seek to provide information and time for the clinician to plan for the potential Multi-Dose Vial activity shortfall.

6.11 Replacing the Saline Bag

⚠WARNINGS

Air embolism hazard. Injury or death may result from infusing air.

- Ensure that air is expelled from the disposable sets before beginning infusion.

Biological contamination can result from reusing the PAS or failure to follow sterile technique.

- Do not reuse or attempt to re-sterilize. Properly discard disposable items after use. If there is any possibility that contamination may have occurred during set up or use, disassemble and setup a new sterile product.


NOTE: Follow your institution's policy for saline administration and aseptic technique.

NOTE: Good clinical practice must be exercised at all time when using the device.

1. Select the Saline tab to go to the System Preparation screen.

2. Press the Replace button.

3. Remove the used saline bag from the hook.

 NOTE: To minimize the introduction of air, perform steps **4** and **5** while holding the used saline bag and spike at a level below the saline pump.

4. Remove the SAS spike from the used saline bag.
5. Insert the SAS spike into the new saline bag.

 NOTE: If air is introduced into the line, re-prime the SAS.

6. Place the saline bag onto the hook at the rear of the cart.
7. Enter the volume of the saline bag using the keypad on the display.
8. Press the Dosing button to return to the Dosing page.

7 Vial Shield and SAS Removal

⚠WARNING

Radiation exposure hazard. ^{18}F-FDG or ^{18}F-NaF emit gamma radiation. Operator should exercise extreme caution when handling radiopharmaceuticals.

- Doses may exceed maximum permissible dose if appropriate precautions are not properly exercised. Follow all individual facility guidelines for the handling and disposal of radioactive materials. DO NOT open the Shielded Chamber without first taking radiation precautions per facility procedures.

Puncture hazard. The SAS contains needles.

- To prevent injury, take care when handling the SAS to control the needles. Dispose of the needles per individual facility guidelines.

Puncture hazard.

- Do not attempt to recap needles.

Radiation exposure hazard. The Intego PET Infusion System may contain radioactive materials.

- Do not reach under cart while the system contains activity. Do not remove access panels while the system contains activity.

Severe operator injury can occur if care is not taken when opening or closing the lids over the Shielded Chamber or the Patient Administration Set compartment.

- Be certain that all fingers are clear before opening or closing the Shielded Chamber Lid or the PAS Compartment Cover.

Both the SAS and the Multi-Dose Vial may contain radiation even after the end-of-the-use period. If possible, allow the used SAS and Multi-Dose Vial to decay to background levels (at least 10 half-lives) before removing. If removal of the used SAS is required before decay to background has been achieved, take radiation precautions per facility radiation safety procedures.

NOTE: The Multi-Dose Vial and used SAS may be left within the Shielded Chamber until the radiation level has decayed sufficiently for safe removal.

7.1 Vial Shield and SAS Removal

1. Remove the Saline Tube from the Saline Pump and disconnect SAS from the Saline Container. The Saline Container and/or the Saline Tube may leak when disconnected. Take precautions to prevent spills.
2. Remove the Swabbable Valve and the Patient Tube from the Swabbable Valve Holder and Air Detector.
3. Open the PAS Compartment Cover. Release the Shielded Chamber Lid Latch and slide the Shielded Chamber Lid towards the front of the cart far enough to expose the two Pinch Valves.
4. Remove the pinch valve covers from the pinch valves.
5. Remove the Waste Tube from the Waste Pinch Valve and the Patient Tube from the PAS Pinch Valve. Press the button on top of the Pinch Valve and pull the SAS out.
6. Completely open the Shielded Chamber Lid and remove the RP Tube from the RP Pump.
7. To remove the SAS Needles from the Multi-Dose Vial, grasp the Needle Insertion Device's handle, lift straight up until it stops, then rotate counter-clockwise until it stops. Release the handle and the Needle Insertion Device will remain in position.
8. Using the Vial Shield Carrying Handle, insert and attach the Vial Shield Access Cap by pushing down and rotating clockwise. Ensure that the Access Cap is fully engaged with the Vial Shield. Do not remove the Vial Shield at this time.
9. Lift the hinged Needle Insertion Device cover to access the Needle Cartridge holder. Remove the Needle Cartridge from the Needle Insertion Device.

10. Remove the SAS Coil from the Dose Calibrator.
11. Using tongs, lift the Waste Container from the Waste Storage and remove the SAS from the Shielded Chamber.
12. Discard the SAS in a designated hot waste disposal container per facility guidelines.
13. Using the Carrying Handle, remove the capped Vial Shield from the Shielded Chamber and return it to the hot lab or radiopharmacy.

8 Training Mode

⚠WARNINGS

Radiation exposure hazard. DO NOT install RP into the Intego PET Infusion System while in Training Mode.

- If the system detects radiation in the Dose Calibrator while operating in Training Mode, it will indicate an error and operation will stop.

Explosion hazard. DO NOT use the Intego PET Infusion System when flammable gases are present.

- Patient injury could result from using the system in the presence of flammable gases (such as anesthetics).

While in Training Mode, the Intego PET Infusion System should only be used for training purposes.

- Do not use the system to infuse fluid into a patient.

NOTE: The Intego PET Infusion System will simulate normal functionality in Training Mode. The user can perform all of the tasks detailed in the "Daily Setup" and "Patient Infusion" chapters of this manual. Saline or other non-radioactive fluid must be substituted for an RP to simulate an infusion.


NOTE: The Intego PET Infusion System cannot perform the Daily QC Constancy and Accuracy Test, Dose Calibrator Check, Linearity check, or Dose Calibrator Calibration while in Training Mode. Bayer recommends that the user performs these tasks prior to using the system in Training Mode.

NOTE: In Training Mode, the pinch valve calibration is simulated and not stored when the system is powered Off.

NOTE: The Intego PET Infusion System cannot simulate occlusions in Training Mode.


8.1 Training Mode

1. Power on the Intego PET Infusion System.
2. The Introduction screen appears. Press the TRAINING button.


NOTE: The screen displays the current mode in the lower right corner of the Introduction screen.

3. If the Intego PET Infusion System was in Clinical Mode, the system will display the message “Changing from Clinical Mode to Training Mode requires a restart. Continue?” Press OK to proceed.


4. The Intego PET Infusion System will display a progress bar while switching to Training Mode. When complete, the system will display a message prompting the user to restart. Press the RESTART button to Restart the system.

5. At the Introduction screen, press the TRAINING button.

 NOTE: When the Intego PET Infusion System is in Training Mode, the Title Bar is colored red and the word "Training" is displayed.

6. At the conclusion of the training exercise, return the Intego PET Infusion System to Clinical Mode (See "Daily Setup"). Remove and discard the used disposable sets per site guidelines.

A Cleaning and Maintenance

⚠WARNINGS

Shock hazard. Patient or operator injury could result from worn cabling or unit disassembly.

- To avoid exposure to potentially hazardous voltages, DO NOT disassemble the Intego PET Infusion System in any way. Worn cabling also creates hazards. If any worn or damaged cables are detected, DO NOT use the system. Contact Bayer for service or replacement.

Severe operator injury can occur if care is not taken when opening or closing the lids over the Shielded Chamber or the Patient Administration Set compartment.

- Be certain that all fingers are clear before opening or closing the Shielded Chamber Lid or the PAS Compartment Cover.

Radiation exposure hazard. The Intego PET Infusion System may contain radioactive materials.

- Do not reach under cart while the system contains activity. Do not remove access panels while system contains activity.

Radiation exposure hazard. Radiation activity exposure is increased once the Access Cap is removed.

- The user should avoid working directly over the open Vial Shield Access Cap to minimize exposure.

Radiation exposure hazard. Using the Shielded Chamber Lid as a work space may expose the operator to excess radiation.

- Do not use the Shielded Chamber Lid as a work space while there are radioactive materials present within the Intego PET Infusion System. Using the Shielded Chamber Lid as a work space or placing items on the Shielded Chamber Lid while it is open may damage the Shielded Chamber Lid, hampering its ability to provide an adequate radiation shield.

Radiation exposure hazard. ^{18}F-FDG or ^{18}F-NaF emit gamma radiation. Operator should exercise extreme caution when handling radiopharmaceuticals.

- Doses may exceed maximum permissible dose if appropriate precautions are not properly exercised. Follow all individual facility guidelines for the handling and disposal of radioactive materials. DO NOT open the Shielded Chamber without first taking radiation precautions per facility procedures.

Puncture hazard. The SAS contains needles.

- To prevent injury, take care when handling the SAS to control the needles. Dispose of the needles per individual facility guidelines.

Operator injury. Batteries are not user replaceable and must be replaced by Bayer personnel.

- Contact Bayer.

⚠CAUTION

Do not remove any covers or disassemble the Intego PET Infusion System.

- Periodically inspect for loose or frayed cables, loose covers, cracks, dents, or loose hardware. Contact Bayer for repairs.

A.1 Cleaning Guidelines

⚠CAUTIONS

Do not use strong cleaning agents and solvents.

- Warm water and a mild soap solution are all that are required to clean the Intego PET Infusion System. Allow all components to dry thoroughly before using.

DO NOT soak or immerse any part of the Intego PET Infusion System or expose to excessive amounts of water or cleaning solutions while cleaning.

- Avoid allowing any fluids to leak inside system components. Improper or careless cleaning methods may result in equipment damage. Wipe components with a gauze pad dampened with a mild soap and water solution.

NOTE: Follow facility guidelines for the clean-up of spilled radioactive pharmaceuticals.

NOTE: For all body fluid spills, follow institutional decontamination procedures. If fluid has leaked inside any component of the Intego PET Infusion System, the affected subassembly should be disassembled and cleaned by Bayer personnel or returned to Bayer.

NOTE: In the event of fluid ingress or spillage on the injection system ensure all equipment and accessory connections are removed, dried, and inspected. Follow hospital policies and procedures or contact Bayer for performing appropriate electrical safety and operational checks prior to use.

Before beginning, disconnect the Intego PET Infusion System from AC power. If necessary, remove the SAS and the Multi-Dose Vial from the system by following the instructions outlined in the "Vial Shield and SAS Removal" section of the manual.

Equipment	Cleaning Guidelines	Cleaning Tools
Interior of the Shielded Chamber	Wipe clean and dry thoroughly.	A soft cloth or paper towel dampened with mild cleaning solution.
Exterior of the Cart	Wipe clean and dry thoroughly.	A soft cloth or paper towel dampened with mild cleaning solution.
Touch Screen	Wipe clean.	A soft, non-abrasive cloth or paper towel dampened with cleaning solution.
Vial Shield	Wipe clean and dry thoroughly. Ensure the vial shield exterior is free from debris.	A soft cloth or paper towel dampened with cleaning solution or warm water.
Needle Insertion Guide	Ensure that the needle insertion guide is free from debris.	A soft cloth or paper towel dampened with cleaning solution or warm water.

A.2 Recommended Maintenance Schedule

⚠WARNING

Using the Intego PET Infusion System without performing the necessary QC processes could result in the patient receiving an incorrect dosage of the RP indicated for use.

- Be certain to perform all necessary QC checks at the recommended intervals.

⚠CAUTION

Intego PET Infusion System malfunction may be caused by failure to perform regular maintenance.

- Regular preventive maintenance is recommended to ensure that the system functions properly. Refer to this manual or contact Bayer for additional information.

NOTE: Remove all radioactive fluids from the Intego PET Infusion System in preparation for a service call.

NOTE: Any problems detected during these or any other procedure should be corrected before using the Intego PET Infusion System in patient procedures.

NOTE: Local regulations or hospital protocol may require electrical leakage checks at more frequent intervals. If this applies, local regulations for leakage must be followed.

NOTE: Failures that occur due to lack of proper maintenance or abuse are not covered under warranty.

NOTE: Bayer makes on-site consulting or consulting references available upon request.

The Intego PET Infusion System must be properly maintained to ensure that it is meeting product requirements. The individual maintenance system and schedule depends upon how each system is used, the type of procedures performed, and the frequency of use. The following maintenance schedule is recommended for the system.

Frequency	Guidelines
Daily	1. Perform the Daily QC procedure before the first use of the day. Refer to the "Daily Setup" chapter of this manual. 2. With the Intego PET Infusion System powered Off, check that the Forward/Reverse drive controller and drive engage switch are not broken or damaged and that they return to the rest position when pushed and released for both the forward and reverse directions. If there is a problem, contact Bayer. 3. Thoroughly inspect the Intego PET Infusion System for any external damage, which could indicate damage to the radiation shielding. If severe damage is found, contact Bayer.
Weekly	With the Intego PET Infusion System powered Off, put the Drive Controller to the full-speed forward position and power On the system. The system should not move. If the system moves, contact Bayer.
Quarterly	Perform the Dose Calibrator Linearity check. Refer to the "Dose Calibrator Linearity Check" section of this manual.
Annually	As part of an annual maintenance program performed by Bayer or an authorized dealer, Bayer recommends that the following checks and procedures be performed. 1. Perform the Electrical Leakage check. 2. Perform the Ground Continuity check. 3. Perform the Dose Calibrator Calibration check (Refer to the "Dose Calibrator Calibration" section of this manual). 4. Calibrate the Display. 5. Perform a complete system performance checkout. Contact Bayer for complete details.
As Needed	Perform Pinch Valve calibration. Refer to the "Calibrating the Pinch Valves" section of this manual.

A.3 Dose Calibrator Linearity Check

NOTE: The Dose Calibrator Linearity check is a required procedure to be performed during the quarterly maintenance schedule.


NOTE: Use ^{18}F to perform the Dose Calibrator Linearity check. The activity should be higher than 925 MBq (25 mCi) to ensure the full range of programmable doses are checked.

NOTE: The Linearity check proceeds until the ^{18}F activity decays to 3.7 MBq (100 µCi). If required, this limit can be configured by Bayer to be as low as 370 kBq (10 µCi).


NOTE: If the lower limit on the Linearity check is set to less than 3.7 MBq (100 µCi), the decay time will increase.

To perform the Dose Calibrator Linearity Check, follow these steps:


1. Press the CONFIGURATION button on the Navigation Bar to access the Configuration screen.

2. The Configuration screen appears. Press the MAINTENANCE button and then press the Linearity tab.

3. Press the NEW MEASUREMENT button to begin a new Linearity check.


4. Enter the activity.

NOTE: The Intego PET Infusion System automatically measures the activity level. If the clinician enters a value for activity, that value will be retained for future reference and reporting.


a. Press the Activity field.

b. Using the keypad, enter the activity and reference time from the assay information. Press the ENTER button when complete.


5. Using the Calibration Source Holder, place vial of ^{18}F into the Dose Calibrator and close the Shielded Chamber Lid.

6. Press the BEGIN MEASUREMENT button.


7. While the check is in progress, the Intego PET Infusion System provides a progress screen. To stop the Linearity check before it is complete, press the ABORT button.


NOTE: During the Linearity check, the system displays the measured activity every 10 minutes.

8. When the check is complete, the Intego PET Infusion System displays a summary screen (the screen shown is an example). Press the PRINT button for a hard copy of the results. Press the DONE button to exit.

9. To export the Linearity data, insert a USB memory device into the USB Port on the rear of the Display and then press the EXPORT button.

A.4 Dose Calibrator Calibration

⚠WARNING

Using the Intego PET Infusion System without performing the necessary QC processes could result in the patient receiving an incorrect dosage of the RP indicated for use.


- Be certain to perform all necessary QC checks at the recommended intervals.

NOTE: Dose Calibrator calibration requires ^{57}Co and ^{137}Cs sources.


NOTE: DO NOT use a rod calibration source to calibrate the system. Only a vial source should be used to calibrate the Dose Calibrator.

NOTE: Bayer recommends that the calibration source is greater than 3.7 MBq (100 µCi) with a 3% source uncertainty. If a 3% or better source is not available, contact Bayer.


1. Press the CONFIGURATION button on the Navigation Bar to access the Configuration screen.
2. The Configuration screen appears. Press the MAINTENANCE button and then press the Calibration tab.


3. Press the BEGIN CALIBRATION button.

4. The Intego PET Infusion System displays a dialog box prompting the clinician to insert the ^{57}Co source into the Dose Calibrator. Using the Calibration Source Holder, insert the ^{57}Co source into the Dose Calibrator and close the Shielded Chamber Lid.


5. Confirm the assay information listed for the first source is correct. If it is not, press the EDIT button to make changes. Once the information is correct, press the OK button.

6. The Intego PET Infusion System measures the source and displays a status window for the ^{57}Co. Adjust the gain to achieve an error less than 1% by clicking in the GAIN cell and entering the desired value using the keypad or by pressing the + and - buttons. The system will re-measure the source. Continue to adjust the gain until an acceptable (<1%) error is achieved.


 NOTE: Bayer recommends using the factory calibration gain located on the Dose Calibrator certificate as the initial gain value.

7. Remove the ^{57}Co source.


8. Press the Cs-137 tab.

9. The Intego PET Infusion System then displays a dialog box prompting the user to insert the ^{137}Cs source into the Dose Calibrator. Using the Calibration Source Holder, insert the ^{137}Cs source into the Dose Calibrator. Confirm the assay information listed for the second source is correct. Press the EDIT button to make changes. Once the information is correct, press the OK button.

10. The Intego PET Infusion System measures the source and displays a status window for the ^{137}Cs source. Adjust the gain to achieve an error less than 1% by clicking in the GAIN field and entering the desired value using the keypad or by pressing the + and - buttons. The system will re-measure the source. Continue to adjust the gain until an acceptable (<1%) error is achieved.

 NOTE: Bayer recommends using the factory calibration gain located on the Dose Calibrator certificate as the initial gain value.

11. If gain adjustments were made on either source, the user must re-measure both sources by performing the previous steps until no gain adjustments are needed to achieve errors of less than 1% for both sources.
12. Once no gain adjustments are needed on either source, press DONE to complete and exit the calibration.

A.5 Geometry Check

Acceptable geometry testing was demonstrated for the Intego Dose Calibrator in accordance with IEC61145: Calibration and Usage of Dose Chamber Systems for Assay of Radionuclides, and NRC Information Notice No. 93-10: Dose Calibrator Quality Control, Section 4.


The test report indicates that the licensee need not repeat geometry testing due to the following:

- The intended use of the Intego Dose Calibrator is for a single isotope: ^{18}F.
- The location of the isotope in the Dose Calibrator is fixed based on the design of the Intego SAS.


A.6 Calibrating the Pinch Valves

1. Press the CONFIGURATION button on the Navigation Bar to access the Configuration Screen.
2. Press the MAINTENANCE button and then press the Calibration tab.


3. Press the VALVE CALIBRATION button.

4. The system displays a message prompting the user to remove all tubing from the Pinch Valves. Ensure that the tubing has been removed and Press OK.

5. When the calibration procedure is complete the system will display a message that calibration is complete. Press OK.

A.7 Replacing the Printer Paper


1. Depress the two locking catches on the sides of the printer (one on either side) and pull the printer cover away from the cart. The cover pivots forward and downward. Spread the two arms holding the roll of paper and remove the spent roll.
2. Place the new roll onto the arms so that the loose end of the paper rolls from on top.

3. Insert the loose end of the roll into the paper guide.

4. Close the printer cover by pressing the release catch while lifting the printer cover back into place. The buttons on either side of the printer cover should both snap back into place.

A.8 Battery Maintenance and System Storage

For optimal battery life, keep the Intego PET Infusion System batteries fully charged by plugging the system into AC power whenever possible. This would include the following:

- Plug the system into AC power overnight.
- Plug the system into AC power immediately after it is driven every time.
- Plug the system into AC power immediately if the Drive System Status (Blue) Status Indicator begins flashing.
- Only drive the cart if the Drive System Status Indicator (Blue) is illuminated continuously.

NOTE: If the Intego PET Infusion System will not be used and cannot be powered for more than 1 week, the user should contact Bayer to prepare the system for storage.


A.9 Disposal of Equipment


When the Intego PET Infusion System is at the end of its usable life cycle, certain precautions should be taken prior to the removal of the system:

- Ensure that all calibration sources have been removed and properly stored.
- Interior and exterior wipe tests and area surveys should be performed, using the appropriate radiation detectors to determine any possible contamination of radioactive materials.
- The Dose Calibrator should be removed from the system and disposed of per manufacturer's instructions.
- Lead is used as the primary shielding material for the Intego PET Infusion System. Lead is a hazardous waste regulated by the U.S. Environmental Protection Agency (EPA) and cannot be discarded in the regular trash or with the radioactive waste. Disposal should be made in accordance with federal, state, and local regulations. Bayer recommends end of life units be sent to reputable firms for electronic recycling in lieu of disposal.
- For equipment disposal information in accordance with European WEEE directives, please visit www.weee.bayer.com
- Remember to wear gloves when handling potentially contaminated items and wash hands thoroughly after handling lead items.
- Contact facility manager for instructions on proper removal or recycling of lead items and batteries.

B Specifications

B.1 Mechanical

	Nominal Value
Weight, INT SYS 200	355 kg (783 lb)
Weight, Vial Shield	8 kg (18lb)

B.2 Power Cable Specifications

The specifications required by the MEDRAD® Intego PET Infusion System relative to the power cable (plug, receptacle, and cord) are:

- Operating Temperature: 60° C minimum
- Receptacle Type: IEC-60320 C13
- Normal Cord Voltage: 300 VAC minimum
- Wire Gauge: 1.00 mm^2 minimum
- Cord Type: IEC 60245-1, Annex A, Designation 53, or IEC 60227-1, Annex A, Designation 53 Certified
- Cord Length: 3 m maximum

The power cable must meet applicable plug, cord, and receptacle specifications including type, voltage, current, and safety approval markings for the country in which the power cable is being used.

B.3 MEDRAD® Intego PET Infusion System to IT Network Connections

Connecting the system to an IT-Network that includes other equipment could result in unidentified risks to patients, operators, or third parties.

The organization responsible for managing the network should identify, analyze, evaluate, and control risks associated with connecting the equipment to the IT-Network.


Subsequent changes to the IT-Network could introduce new risks and require additional analysis. For example:

- Changes in the IT-network configuration
- Connecting additional items to the IT-Network
- Disconnecting items from the IT-Network
- Updating equipment connected to the IT-Network
- Upgrading equipment connected to the IT-Network

B.4 Radiation Shielding Profile

B.4.1 Vial Shield Radiation Profile

Exposure Rate at the surface with 27.75 GBg (750 mCi) in 30 ml saline in the vial.


B.4.2 Radiation Exposure Reduction Specification and Testing Results

Meets the NRC Code of Federal Regulations for Public Exposure, which limits radiation exposure to 0.02mSv (2mrem) in any one hour period. With the maximum activity loaded, the maximum radiation exposure at the front, rear, left, right, and top of the cart does not exceed 0.02mSv/hr (2mrem/hr) when the system is idle. The maximum radiation exposure at the front, rear, left, right, and top of the cart does not exceed 0.05mSv/hr (5mrem/hr) while holding a prepared dose prior to delivery to the patient.

Measured Exposure:

- Peak rate equal to or less than 0.014 mSv/hr (1.4 mRem/hr) at 30.5 cm (12 in) from any surface of the cart with 27.75 GBq (750 mCi) in the vial.
- Peak rate less than 0.001 mSv/hr (0.1 mRem/hr) at 30.5 cm (12 in) from the surface in the Operator Position with 27.75 GBq (750 mCi) in the vial.

Operator Exposure:

B.5 Environmental

	Nominal Value	Notes
Non-Operating Conditions		
Temperature, maximum ambient	60°C (140°F)	--
Temperature, minimum ambient	-20°C (-4°F)	--
Humidity, maximum	100% R.H.	Non-condensing
Humidity, minimum	5% R.H.	Non-condensing
Pressure, maximum barometric	57 kPa (8 psi)	--
Pressure, minimum barometric	106 kPa (15 psi)	--
Vial Shield Temperature, maximum ambient	70°C (158°F)	--
Vial Shield Temperature, minimum ambient	-40°C (-40°F)	--
Operating Instructions		
Temperature, maximum ambient	40°C (104°F)	--
Temperature, minimum ambient	10°C (50°F)	--
Humidity, maximum	90% R.H.	--
Humidity, minimum	20% R.H.	--
Pressure, maximum barometric	70 kPa (10 psi)	--
Pressure, minimum barometric	106 kPa (15 psi)	--

B.6 Electrical

	Nominal Value	Notes
Mains Power Connection		
Mains supply voltage, maximum	240 VAC	--
Mains supply voltage, minimum	100 VAC	--
Mains supply frequency, maximum	60 Hz	--
Mains supply frequency, minimum	50 Hz	--
Power consumption, maximum	300 VA	--
AC Power Line Fuse	3.0A, 250V 5X20MM, IEC TYPE T, HIGH	
Safety		
Electrical leakage current, earth	<500 µA (Normal Condition)	See note below
Electrical leakage current, chassis (touch)	<100 µA (Normal Condition)	See note below
Electrical leakage current, patient connection	< 100 µA, a.c. (Normal Condition)	See note below
Ground continuity	< 0.2 Ohms	Measured from power cable ground terminal to the system power and grounding points. See note below.
NOTE: Per IEC 60601-1, the Intego PET Infusion System is designated as a Class 1 medical device with a type BF applied part.		

B.6.1 Protection Against the Ingress of Fluids

Per IEC 60529, the Intego PET Infusion System has been classified as drip proof. This is indicated by the IPX1 designation.


NOTE: In the event of fluid ingress or spillage on the injection system ensure all equipment and accessory connections are removed, dried, and inspected. Follow hospital policies and procedures or contact Bayer for performing appropriate electrical safety and operational checks prior to use.

B.6.2 Mode of Operation

Per IEC 60601-1, the mode of operation for the Intego PET Infusion System is continuous operation. The system is capable of operation under normal load for an unlimited period, without excessive temperature being developed.

B.7 Fluid Delivery

	Nominal Value	Notes
RP Vial		
Maximum vial activity for shielding effectiveness	27.8 Gbq (750 mCi)	--
Maximum vial activity for dose preparation	25.9 GBq (700 mCi)	System will not move fluid if vial activity exceeds this level
Maximum RP vial volume	30 ml	30H vial. See note below.
Unextractable RP volume	0.64 ml	30H vial
Patient doses/vial	11	Based on 25.9 GBq (700 mCi) in vial, 555 MBq (15 mCi) patient doses and 30 minute interval between doses.
NOTE: A wide variety of vials from 10 ml to 30 ml are compatible with the system, with use of one of six different vial shields. See "Vials and Vial Shields"..		
Concentration		
Maximum RP concentration	3.7 GBq/ml (100mCi/ml)	--
Minimum RP concentrations	12.3 MBq/ml (0.33mCi/ml)	See note below

	Nominal Value	Notes
Dose Activity		
Maximum dose activity	925 MBq (25 mCi)	Maximum Dose Activity is reduced for RP concentrations ≤ 308 MBq/ml (8.3 mCi/ml). See Concentration-Based Dose limits graph.
Minimum dose activity	37 MBq (1 mCi)	Minimum Dose Activity is increased for RP concentrations ≥ 296 MBq/ml (8 mCi/ml). See Concentration-Based Dose limits graph.
Dose activity accuracy	± 2%	Exclusive of Dose Calibrator calibration factor (if using Dose Calibrator accuracy of ± 3%, then able to deliver ^{18}F-RP within ± 4.7% of the measured dose). See note below.
Half-life (^{18}F)	109.74 min	Value used for all decay calculations.
Dose Volume		
Maximum RP dose volume	3.0 ml	--
Minimum RP dose volume	0.125 ml	--
Minimum dose fluid volume	35 ml	Includes RP and saline flush. Additional Flush volume is added to this value. See note below.
Test Injections - Optional and user configurable		
Minimum saline volume	5 ml	--
Maximum saline volume	30 ml	--
Additional Flush - Optional and user configurable		
Minimum saline volume	0 ml	--
Maximum saline volume	30 ml	--

	Nominal Value	Notes
Fluid to Patient		
Flow rate	0.5 ml/s	See note below.
	1.0 ml/s	See note below.
NOTE: Dose activity and fluid delivery accuracies are only guaranteed if Intego PET Infusion System disposables are used and if the patient connection uses no attachment more restrictive than a 24 gauge x 0.75 inch (19 mm) long catheter or a 23 gauge x 0.5 inch (12 mm) long steel needle. Use of more restrictive catheters or needles will affect flow rate and delivered volume and may affect delivered dose activity.		
PAS		
Extended length	155 cm (60 in.)	Approximate length
SAS		
Waste container fluid capacity	375 ml	--
Waste container activity capacity	3.7 GBq (100 mCi)	--

B.8 Mobility

	Nominal Value	Notes
Speed, maximum forward, high range	0.67 m/s (2.2 ft/s)	--
Speed, maximum forward, low range	0.4 m/s (1.3 ft/s)	--
Speed, maximum reverse	0.4 m/s (1.3 ft/s)	--
Recharging time	8 h	See note below.
Braking	Automatic electric	Brakes are applied by the mobility controller when the driver controls are released.
Manual Override	--	Pulling the handle releases the drives, allowing the system to be moved manually.
Caster Brakes	Manual	User-applied to maintain position.
NOTE: Mobility battery recharge time is approximately 3x the discharge time. Users are advised to connect the Intego PET Infusion System to mains power whenever it is not moving.		

B.9 MEDRAD® Intego Wireless Connection Feature

The wireless connection feature on the system (INT SYS 200 with POC 127) provides a connection to a hospital network for the purposes of:

- importing a schedule from a modality worklist
- providing infusion results in the form of a graphic into a PACS system

⚠CAUTIONS

Connection of the MEDRAD® Intego system to an IT network that includes equipment other than what has been validated could result in a previously unidentified risk to patients, operators, and/or third parties.

- The hospital IT department responsible for the network should identify, analyze, evaluate, and control these risks. (IEC 80001-1:2010 provides guidance to address these.)

Any subsequent changes to the IT network may introduce new risks and require additional analysis. Such changes may include:

- Network configuration changes
- Connection of additional products or items to the network
- Disconnecting items from the network
- Updating and/or upgrading equipment connected to the network

Failure to adhere to the identified network standards may lead to data corruption while transferring dose infusion records or dose absorption reports from MEDRAD® Intego to PACS, resulting in incorrect patient records.

B.9.1 Required Characteristics

In order for the system's wireless connection feature to operate, the hospital IT network must have the following characteristics:


1. A wireless network conforming to 802.11 b/g specifications

 NOTE: Default support is for 2.4 GHz band with an option to enable 5 GHz band.

2. A network incorporating wireless security based on WPA and WPA2 encryptions
 a. It is recommended that the MEDRAD® Intego wireless communication is operated over a secure, NOT open, network.
 b. Recommended wireless security measures:
 - WPA2 PSK - AES
 - WPA PSK - TKIP

B.9.2 Required Configuration

The configuration of the system into the hospital network is managed by the hospital IT staff. A wireless router must be configured by the hospital IT staff to allow connection of the MEDRAD® Intego system. The system must have configured license files and appropriate network settings to allow interoperability on the network connection. Proper network security must be employed. (See “Required Characteristics” on page B - 106. for details.)

B.9.3 Intended Information Flow

- **Modality Worklist:** This enables a piece of imaging equipment (a modality) to obtain details of patients and scheduled examinations electronically.
- **PACS:** Picture archiving and communication system (PACS) is a medical imaging technology that provides economical storage and convenient access to images from multiple modalities. Electronic images and reports are transmitted digitally via PACS, eliminating the need to manually file, retrieve, or transport film jackets.
- **DICOM:** Digital Imaging and Communications in Medicine (DICOM) is a standard for handling, storing, printing, and transmitting information in medical imaging. It includes a file format definition and a network communications protocol. The communication protocol is an application protocol that uses TCP/IP to communicate between systems. DICOM files can be exchanged between two entities that are capable of receiving image and patient data in DICOM format.

B.9.4 Network Communication Technical Specifications

The wireless radio specifications used in MEDRAD® Intego are:

Standard	IEEE 802.11 b/g/n Standard ◆ Supports 2.4 GHz by default ◆ If enabled, can support 5 GHz band
Frequency Range	IEEE 802.11b/g, 802.11g/n HT20: 2412MHZ ~ 2472 MHz
	IEEE 802.11g/n HT40: 2422MHz ~ 2462 MHz
Transmit Power (EIRP)	IEEE 802.11b: 12.38 dBm
	IEEE 802.11g: 12.78 dBm
	IEEE 802.11g/n HT20: 16.65 dBm
	IEEE 802.11g/n HT40: 16.42 dBm
Channel Number	IEEE 802.11b/g, 802.11g/n HT20: 13 channels
	IEEE 802.11g/n HT40: 9 channels
Security	WPA PSK- TKIP
	WPA2 PSK - AES (Recommended)

B.9.5 FCC RF Radiation Exposure Statement

1. This Transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.
2. This equipment complies with FCC RF radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with a minimum distance of 20 centimeters between the radiator and your body.

C Compliance to IEC 60601-1-2 / 2nd, 3rd, and 4th Editions

US: The Intego PET Infusion System, although exempt per 47CFR15.103(c)(e), has been verified to conform to the requirements of 47CFR15 as a Class B digital device.

EU: The Intego PET Infusion System has been verified to conform to the requirements of the Medical Device Directive (93/42/EEC).

The MEDRAD® Intego PET Infusion System complies with the requirements of:

IEC 60601-1-2: Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral standard: Electromagnetic compatibility – Requirements and tests

CISPR 11: Industrial, scientific and medical (ISM) radio-frequency equipment- Electromagnetic disturbance characteristics – Limits and methods of measurement

IEC 61000-3-2: Electromagnetic compatibility (EMC) – Part 3-2: Limits – Limits for harmonic current emissions (equipment input current ≤ 16 A per phase) (This does not apply to Class A equipment.)

IEC 61000-3-3: Electromagnetic compatibility (EMC)- Part 3-3: Limits - Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems, for equipment with rated current ≤ 16 A per phase and not subject to conditional connections) (This does not apply to Class A equipment.)

IEC 61000-4-2: Electromagnetic compatibility (EMC) – Part 4-2: Testing and measurement techniques – Electrostatic discharge immunity test

IEC 61000-4-3: Electromagnetic compatibility (EMC) – Part 4-3: Testing and measurement techniques – Radiated, radio-frequency, electromagnetic field immunity test

IEC 61000-4-4: Electromagnetic compatibility (EMC) – Part 4-4: Testing and measurement techniques – Electrical fast transient/burst immunity test

IEC 61000-4-5: Electromagnetic compatibility (EMC) – Part 4-5: Testing and measurement techniques – Surge immunity test

IEC 61000-4-6: Electromagnetic compatibility (EMC) – Part 4-6: Testing and measurement techniques – Immunity to conducted disturbances, induced by radio frequency fields

IEC 61000-4-8: Electromagnetic compatibility (EMC) – Part 4-8: Testing and measurement techniques – Power frequency magnetic field immunity tests

IEC 61000-4-11: Electromagnetic compatibility (EMC) – Part 4-11: Testing and measurement techniques – Voltage dips, short interruptions and voltage variations immunity tests

This system is in compliance to IEC-60601-1-2 / 2nd, 3rd, and 4th edition standards. Special precautions regarding Electromagnetic Compatibility (EMC), are required for installation and use of this system. Detailed EMC information contained in this addendum is intended to reflect conformance to IEC-60601-1-2 / 2nd, 3rd, and 4th edition standards.

⚠WARNINGS

For proper operation, use only accessories and options provided by Bayer that are designed specifically for the system. Other non-Bayer approved accessories or options may cause equipment damage or may result in increased emissions or decreased immunity of the system. System accessories listed in the operation manual comply with the requirements of electromagnetic emissions and immunity standards IEC-60601-1-2 / 2nd, 3rd, and 4th edition.

Do not use system adjacent to or stacked with other equipment. Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If adjacent or stacked use is necessary, the system and the other equipment should be observed to verify normal operation in the configuration in which it will be used.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the injector system unless a greater separation distance is required as indicated by the equation. Otherwise, degradation of the performance of this equipment could result.

⚠CAUTION

System may disarm or fail to operate when exposed to high magnetic fields. Portable and mobile RF communications equipment can affect the system.

Recommended separation distances between portable and mobile RF communications equipment and the system

The System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the system as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 KHz to 80 MHz $d = [3.5/V_1]\ \sqrt{p}$	**80 MHz to 800 MHz** $d = [3.5/V_1]\ \sqrt{p}$	**800 MHz to 2.7 GHz** $d = [7/V_1]\ \sqrt{p}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

THE SYSTEM REQUIRES SPECIAL PRECAUTIONS REGARDING EMC. Install and put into service according to the EMC information provided below:

Guidance and manufacturer's declaration - electromagnetic emissions

The system is intended for use in the electromagnetic environment specified below. The customer or user of the system should assure that it is used in such an environment.

Emissions test	**Compliance**	**Electromagnetic environment - guidance**
RF emissions CISPR 11	Group 1	The System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The emission characteristics of this system make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). If the system is used in a residential environment (for which CISPR 11 Class B is normally required), this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or reorienting the equipment.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable	

Guidance and manufacturer's declaration - electromagnetic immunity		
The system is intended for use in the electromagnetic environment specified below. The customer or user of the system should assure that it is used in such an environment.		
Immunity test	**IEC 60601 Test Compliance Level**	**Electromagnetic environment - guidance**
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2, ±4, ±8, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with a synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for a.c. mains ±1 kV for I/O ports	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± -0.5 kV, ± -1 kV, ± -2 kV line to ground ± -0.5 kV, ± -1 kV line to line	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips IEC 61000-4-11	100% Vac for 0.5 cycles at 0°, 45°, 90°,135°,180°, 225°, 270°, 315	Mains power quality should be that of a typical commercial or hospital environment. If the user of the system requires continuous operation during power mains interruptions, it is recommended the system be powered from an uninterruptible power supply or battery.
	100% Vac for 1.0 cycles at 0°	
	30% Vac for 30 cycles at 0°	
	100% Vac for 250 (50Hz) cycles or 300 (60Hz) cycles at 0°	
Voltage interruptions IEC 61000-4-11	0% a.c. 250(50 Hz) or 300(60 Hz) at 0°	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Guidance and manufacturer's declaration - electromagnetic immunity

The system is intended for use in the electromagnetic environment specified below. The customer or user of the system should assure that it is used in such an environment.

<table>
<tr><th>Immunity test</th><th>IEC 60601 Test Compliance Level</th><th>Electromagnetic environment - guidance</th></tr>
<tr><td>Conducted RF
IEC 61000-4-6</td><td>3Vrms from 150kHz to 80 MHz at 80% AM 1kHz
6Vrm, 80% AM 1kHz at ISM frequencies listed below:
<table><tr><th>Frequency (MHz-ISM List)</th><th>Test Level (Vrms)</th></tr><tr><td>1.8 - 2.0</td><td>6</td></tr><tr><td>3.5 - 4.0</td><td>6</td></tr><tr><td>5.3 - 5.4</td><td>6</td></tr><tr><td>6.765 - 6.795</td><td>6</td></tr><tr><td>7.0 - 7.3</td><td>6</td></tr><tr><td>10.1- 10.15</td><td>6</td></tr><tr><td>13.553 - 13.567</td><td>6</td></tr><tr><td>14.0 - 14.2</td><td>6</td></tr><tr><td>18.07 - 18.17</td><td>6</td></tr><tr><td>21.0 - 21.4</td><td>6</td></tr><tr><td>24.89 - 24.99</td><td>6</td></tr><tr><td>26.957 - 27.283</td><td>6</td></tr><tr><td>28.0 - 29.7</td><td>6</td></tr><tr><td>40.66 - 40.70</td><td>6</td></tr><tr><td>50.0 - 54.0</td><td>6</td></tr></table></td><td>WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the injector system unless a greater separation distance is required as indicated by the equation. Otherwise, degradation of the performance of this equipment could result.

Recommended separation distance
$d = 1.17\sqrt{p}$</td></tr>
<tr><td>Radiated RF
IEC 61000-4-3</td><td>3Vrms from 80 MHz to 2.7 GHz at 80% AM 1kHz
and specific ISM bands listed below:
<table><tr><th>Frequency (MHz)</th><th>Modulation Type</th><th>Modulation Frequency</th><th>Field Strength (Volts/meter)</th></tr><tr><td>385</td><td>Pulse</td><td>18 Hz</td><td>27</td></tr><tr><td>450</td><td>Pulse</td><td>18 Hz</td><td>28</td></tr><tr><td>710</td><td>Pulse</td><td>217 Hz</td><td>9</td></tr><tr><td>745</td><td>Pulse</td><td>217 Hz</td><td>9</td></tr><tr><td>780</td><td>Pulse</td><td>217 Hz</td><td>9</td></tr><tr><td>810</td><td>Pulse</td><td>18 Hz</td><td>28</td></tr><tr><td>870</td><td>Pulse</td><td>18 Hz</td><td>28</td></tr><tr><td>930</td><td>Pulse</td><td>18 Hz</td><td>28</td></tr><tr><td>1720</td><td>Pulse</td><td>217 Hz</td><td>28</td></tr><tr><td>1845</td><td>Pulse</td><td>217 Hz</td><td>28</td></tr><tr><td>1970</td><td>Pulse</td><td>217 Hz</td><td>28</td></tr><tr><td>2450</td><td>Pulse</td><td>217 Hz</td><td>28</td></tr><tr><td>5240</td><td>Pulse</td><td>217 Hz</td><td>9</td></tr><tr><td>5500</td><td>Pulse</td><td>217 Hz</td><td>9</td></tr><tr><td>5785</td><td>Pulse</td><td>217 Hz</td><td>9</td></tr></table></td><td>$d = 1.17\sqrt{p}$ 80 MHz to 800 MHz

$d = 2.33\sqrt{p}$ 800 MHz to 2.7 GHz

Where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).

Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,[a] should be less than the compliance level in each frequency range.[b]

Interference may occur in the vicinity of equipment marked with the following symbol:

 Non-ionizing Radiation Symbol (IEC TR 60878, 5140)</td></tr>
</table>

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the system is used exceeds the applicable RF compliance level above, the system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the system.
b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

D Troubleshooting Tips

D.1 Recovering from Priming Issues

⚠WARNING

Radiation exposure hazard. ^{18}F-FDG or ^{18}F-NaF emit gamma radiation. Operator should exercise extreme caution when handling radiopharmaceuticals.

- Doses may exceed maximum permissible dose if appropriate precautions are not properly exercised. Follow all individual facility guidelines for the handling and disposal of radioactive materials. DO NOT open the Shielded Chamber without first taking radiation precautions per facility procedures.

⚠CAUTION

Incorrect installation of the SAS may result in improper functioning of the Intego PET Infusion System.

- Ensure the Saline Container is sufficiently spiked.
- Confirm the tubing is placed correctly in the Saline Pump.
- Confirm the tubing is oriented on the Saline Pump rollers properly.
- Ensure the tubing is not pinched in the Saline Pump door.
- Ensure the Saline Disconnect Luer is against the Saline Pump.
- Confirm that the Needle Cartridge Clip is inserted with proper orientation into the Needle Insertion Device.
- Confirm the collar on the tubing is positioned properly in the RP Pump.
- Confirm the tubing is inserted into the Pinch Valves correctly.
- Confirm the tubing is fully inserted into the Air Detector.
- Confirm the RP Pump is closed and the tubing is in the pump.
- Confirm that all tubing is laying flat and not being pinched by the lid.

D.2 Recovering from Activation of Air Detector

⚠WARNING

Radiation exposure hazard. ^{18}F-FDG or ^{18}F-NaF emit gamma radiation. Operator should exercise extreme caution when handling radiopharmaceuticals.

- Doses may exceed maximum permissible dose if appropriate precautions are not properly exercised. Follow all individual facility guidelines for the handling and disposal of radioactive materials. DO NOT open the Shielded Chamber without first taking radiation precautions per facility procedures.

D.2.1 During Priming

1. Confirm that the SAS is properly inserted into the Air Detector. If necessary, slide the Shielded Chamber Lid open to expose the Air Detector and properly insert the SAS.
2. Re-prime the system.
3. The Display displays a message stating that the SAS has already been primed. Press the YES button to continue.

D.2.2 During a Saline Test Inject

1. Confirm that the SAS is properly inserted into the Air Detector. If necessary, slide the Shielded Chamber Lid open to expose the Air Detector and properly insert the SAS.
2. Check the PAS for air, including the connection points with the patient catheter and the SAS.
3. If an air pocket is detected, disconnect the PAS from the patient catheter and prime the tubing set again.
4. Reconnect the PAS to the patient catheter and initiate the saline test infusion.

D.2.3 During an Infusion

1. Confirm that the SAS is properly inserted into the Air Detector. If necessary, slide the lid to the Shielded Chamber open only far enough to expose the Air Detector and properly insert the SAS.
2. Check the PAS for air, including the connection points with the patient catheter and the SAS.
3. If an air pocket is detected disconnect the PAS from the patient catheter and discard.
4. Press the DISCARD button, which discards the current dose.
5. Attach a new PAS.
6. Prime the new PAS.
7. Resume normal operation.

D.3 Recovering a Dose Due to a System Failure

⚠WARNING

Radiation exposure hazard. ^{18}F-FDG or ^{18}F-NaF emit gamma radiation. Operator should exercise extreme caution when handling radiopharmaceuticals.

- Doses may exceed maximum permissible dose if appropriate precautions are not properly exercised. Follow all individual facility guidelines for the handling and disposal of radioactive materials. DO NOT open the Shielded Chamber without first taking radiation precautions per facility procedures.

D.3.1 With a Dose of RP in the Dose Calibrator

1. Restart the system. If the error condition does not persist, proceed with the infusion as normal.
2. If the error condition persists after the restart, the dose in the Dose Calibrator can be infused manually.
 a. Open the Saline Pump and remove the tubing from the pump.
 b. Disconnect the Saline Container from the SAS at the Saline Disconnect Luer adjacent to the Saline Pump.
 c. Insert the Saline Container tubing set into the Saline Pump to prevent saline from leaking onto the floor.
 d. Open the Shielded Chamber Lid only far enough to expose the PAS Pinch Valve and remove the SAS from the Pinch Valve.
 e. Close the Shielded Chamber Lid.
 f. Connect a syringe filled with 40 ml of saline to the free end of the SAS.
 g. Manually infuse the 40 ml of saline into the SAS using the syringe to push the RP dose through the system and into the patient.
 h. After the infusion, disconnect the patient from the tubing set.
 i. If necessary, remove the SAS and the Multi-Dose Vial from the Shielded Chamber following the instructions found in the "Vial Shield and SAS Removal" section.
 j. Secure the Intego PET Infusion System and contact Bayer.

D.3.2 With RP in the Vial Shield

⚠CAUTION

Puncture hazard. The SAS contains needles.

- To prevent injury, take care when handling the SAS to control the needles. Dispose of the needs per individual facility guidelines.

1. Restart the system. If the error condition does not persist, proceed with the infusion as normal.
2. If the error condition persists after the restart, the RP in the Vial Shield can be used for hand injection.

Hand Injection Method - Option 1

a. Make certain that the Vial Access Cap is ready to be reinstalled onto the Vial Shield.

b. Open the Shielded Chamber Lid.

c. Using the handle of the Needle Insertion Device, immediately pull the needles directly out of the Vial Shield. DO NOT raise the shield lid over the Needle Insertion Device.

d. Rotate the Needle Insertion Device counter-clockwise and immediately insert the Vial Access Cap onto the Vial Shield using the handle.

e. Remove the Vial Shield from the Vial Shield Compartment using the Vial Shield handle.

f. Rotate the Needle Insertion Device clockwise and lower it into place. Close the Shielded Chamber Lid.

g. Return the Vial Shield containing the RP directly to the facility hot lab where individual patient doses can be prepared.

Hand Injection Method - Option 2


a. Remove the PAS from the SAS and discard in the appropriate receptacle.


b. Clean the Swabbable Valve on the exposed end of the SAS using proper sterile technique.

c. Attach a sterile 10 ml syringe in a syringe shield to the Swabbable Valve.

d. Open the Shielded Chamber Lid.

e. Clamp off the Saline Line directly before the confluence using a slide clamp. (Some examples of appropriate slide clamps are Qosina 12053 and Qosina 12044).

f. Remove the tubing from the RP Pump and close the RP Pump lid.

g. Remove the tubing from the PAS Pinch Valve.

h. Close the Shielded Chamber Lid.

i. Using the syringe, extract 7 ml from the SAS. Assay the syringe. It should have only very low activity. Discard the syringe according to site policy. If the measured activity exceeds 740 MBq (20 mCi), skip to step m.

j. Load an empty, sterile 10 ml syringe into a syringe shield. Wipe the Swabbable Valve, connect the shielded 10 ml syringe to the Swabbable Valve, and extract an additional 7 ml. Assay the syringe. If the measured activity exceeded 740 MBq (20 mCi), skip to step m. Otherwise, discard the syringe according to site policy.

k. Load an empty, sterile 10 ml syringe into a syringe shield. Wipe the Swabbable Valve, connect the shielded 10 ml syringe to the Swabbable Valve, and extract 0.5 ml. (If it is not possible to measure 0.5 ml accurately near the start of the syringe, preload 5 ml of saline or WFI into the syringe before attaching to the Swabbable Valve.) Assay the syringe. If the measured activity exceeded 740 MBq (20 mCi), skip to step m. Otherwise, discard the syringe according to site policy.

l. Repeat step k. up to 3 more times (for a total of 4 x 0.5 ml). If significant activity has not been measured in an assay, there is a problem with the SAS or with the needles in the Multi-Dose Vial. In this instance, revert to “Hand Injection Method - Option 1."

m. Once significant activity has been detected, individual doses of RP can be extracted by hand. Attach a shielded, sterile syringe to the Swabbable Valve.

n. Extract the desired volume of RP and remove the syringe. Assay the syringe in a separate Dose Calibrator before injection.

o. Clean the Swabbable Valve and repeat steps m and n as necessary.

D.4 PAS Occlusion Recovery

⚠WARNINGS

Radiation exposure hazard. ^{18}F-FDG or ^{18}F-NaF emit gamma radiation. Operator should exercise extreme caution when handling radiopharmaceuticals.

- Doses may exceed maximum permissible dose if appropriate precautions are not properly exercised. Follow all individual facility guidelines for the handling and disposal of radioactive materials. DO NOT open the Shielded Chamber without first taking radiation precautions per facility procedures.

In case of occlusion, do not disconnect disposables before discarding dose, as liquid leakage can occur.


This section describes how to identify a potential occlusion and details the method to recover from an occluded PAS during an infusion.

1. At the end of the infusion, if more than 2% of the dose activity remains in the Dose Calibrator, the system will allow the user to resume the infusion or discard the dose. Bayer recommends that the user check the lines for occlusions (clamped tubing, kinked crushed tubing, or high-gauge needles or catheters, closed valves between the PAS and the patient) prior to resuming.

NOTE: During PAS Occlusion Recovery, the system will allow the clinician to attempt to infuse the remaining RP a total of 4 times. If the activity remains in the Dose Calibrator following 4 attempts to infuse, the RP will be discarded.

2. If the user selects the DISCARD option or if activity remains in the Dose Calibrator after all attempts to resume, the remaining activity will be discarded into the Waste Container. The status bar will display any messages relating to the infusion.

3. The history will update and an icon will appear, indicating a partial infusion.
4. If the occlusion remains after the user attempts to resume the infusion or discard the dose, the user should check for occlusions in the Saline Container or Saline Line.

D.5 System Messages

The following is a list of messages the Intego PET Infusion System could display and possible corrections:

NOTE: Follow the error instructions on the display and/or contact Bayer for further support.

D.5.1 Critical Error Messages

NOTE: These errors require the system to be restarted.

Code	Error Message	Explanation
E1000	The measured concentration does not match the expected value based on the <Selected RP> assay information.	◆ Confirm that the assay information was entered correctly. Re-enter if needed. ◆ Confirm that the assay information provided by the pharmacist is correct. Re-enter if needed. ◆ If the assay corrections do not correct the problem, replace the Multi-Dose Vial. ◆ If problem persists, recalibrate the Dose Calibrator. ◆ If the problem continues, record the error message and contact Bayer.
E1055	Activity detected while in Training Mode.	◆ Remove activity from the Dose Calibrator and restart the Intego PET Infusion System in Training Mode. ◆ If the problem continues, record the error message and contact Bayer.

Code	Error Message	Explanation
Any other E code	N/A	◆ Restart the system. ◆ If the message reappears, record the error message and contact Bayer.

D.5.2 Recoverable Errors / Messages

Error Message	Explanation
A dose cannot be prepared from this vial. Replace vial.	◆ No activity is available to extract. ◆ Replace vial.
Activity concentration is less than the acceptable minimum.	◆ Concentration may have decayed below minimum. Verify assay information was entered correctly and reprime. ◆ Verify SAS is installed properly and reprime. ◆ Replace vial and reprime.
Activity concentration too high. Check assay information.	◆ Confirm that the assay information was entered correctly. Re-enter if needed. ◆ Confirm that the assay information provided by the pharmacist is correct. Re-enter if needed. ◆ If the assay was entered correctly, wait until the concentration decays to 3.7 GBq/ml (100 mCi/ml).
Additional saline flush required. Press OK to flush to waste.	◆ A prior dose infusion was aborted or high residual activity remains from priming the SAS. ◆ It should not be necessary to re-prime the SAS, but if problem persists, re-prime the SAS. Refer to the "Recovering from Priming Issues" section.
Air detected in SAS. Need to prime.	◆ Confirm that the SAS is correctly inserted into the Air Detector. ◆ Reprime the SAS. Refer to the "Recovering from Priming Issues"and the "Recovering from Activation of Air Detector" sections.
Assay information incomplete. Enter all required data.	◆ Enter all required information. Source lot number, activity, and date/time are all required. Refer to "Performing Daily QC" section and "Dose Calibrator Calibration" section of this manual.
Calibration already in progress.	◆ Wait until the measurement of the source is complete before proceeding to the next step.
Check printer paper.	◆ Check printer and replace paper as needed. Refer to the "Replacing the Printer Paper" section of this manual. ◆ Close the printer door.
Date/Time write failure. Please try again.	◆ Re-enter the date/time. Refer to the "System Configuration" chapter.
Dose Calibrator is detecting high background activity. Check for nearby sources.	◆ Check for nearby sources. ◆ Check for the calibration source in the Dose Calibrator. Remove if still present. ◆ Check for nearby sources of radiation. If necessary, remove the other source(s) or move the cart. ◆ Check for contamination in the room or in the system, including possible radiopharmaceutical leak from the SAS. Clean the system per site guidelines and replace the SAS.
Dose Calibrator serial number mismatch. Recalibration required.	◆ Recalibrate the Dose Calibrator. Refer to the "Dose Calibrator Calibration" section of this manual.
Dose cannot be infused. Discard and prepare new dose.	◆ Discard the current dose and prepare a new dose.
Error reading preferences, resetting system to defaults.	◆ Re-enter system preferences on Configuration screen.

Error Message	Explanation
Infuse or discard current dose before preparing another.	◆ A patient dose is already prepared. Either infuse or discard the current dose before attempting to prepare another dose or perform a test injection.
Invalid date/time entered.	◆ Re-enter the date/time. Refer to the "System Configuration" chapter of this manual.
Measured activity is higher than expected. Check assay information and prime again.	◆ Confirm that the assay information was entered correctly. Re-enter if needed. Refer to the "Entering RP Assay and Saline Information" section. ◆ Confirm that the assay information provided by the pharmacist is correct. Re-enter if needed. Reprime.
Measured activity is low. Ensure vial and SAS are properly installed. Replace if necessary.	◆ Confirm that the assay information was entered correctly. Re-enter if needed. Refer to the "Entering RP Assay and Saline Information" section. ◆ Confirm that the assay information provided by the pharmacist is correct. Re-enter if needed. ◆ Ensure vial and SAS are properly installed. Refer to the "Recovering from Priming Issues" section.
Operation will exceed Waste Container capacity. Replace SAS.	◆ The Waste Container is full. Replace the SAS.
Pinch Valve error: <Valve Error Details> See Operation Manual.	◆ Confirm that the SAS is properly installed into the Pinch Valves. Refer to the "Recovering from Priming Issues" section. ◆ Valve Error Details will be one of the following phrases: • waste valve failed. • patient valve failed. • both valves failed. • both valves open simultaneously. • tubing absent from waste valve. • tubing absent from patient valve. • tubing absent from both valves. • patient valve failure, tubing absent from waste valve. • waste valve failed, tubing absent from patient valve. ◆ If necessary, restart the system. If the message reappears, record the error message and contact Bayer.
Prime PAS.	◆ The PAS has not been primed and the operation cannot be completed. Prime the PAS before proceeding.
Printer error occurred.	◆ Press the OK button and the system should correct the issue. ◆ If the message reappears, restart the system. ◆ If the message persists, record the error message and contact Bayer.
Residual activity is higher than the acceptable background level. See Operation Manual.	◆ High residual activity remains from priming the SAS. ◆ Refer to the "Recovering from Priming Issues" section. ◆ Re-prime the SAS. ◆ If re-priming does not resolve the issue, check for contamination and replace the SAS.
Saline is low. Please replace.	◆ Replace the saline container. Refer to the "Daily Setup" chapter.
Saline prime failed. Ensure SAS is properly installed. replace if necessary.	◆ Confirm that the SAS is correctly installed. Refer to the "Recovering from Priming Issues" section. ◆ Replace the SAS if necessary.

Error Message	Explanation
SAS detected. Press Cancel to use or remove SAS and press OK.	◆ The system requires that the vial and SAS be removed before the Multi-Dose Vial assay information can be reset. ◆ The SAS should only be used with a single Multi-Dose Vial. Change the SAS each time a new Multi-Dose Vial is installed. Refer to the "Entering RP Assay and Saline Information" section. ◆ Press CANCEL if the REMOVE button was pressed by mistake.
<Selected RP> or saline information missing. Enter all required data.	◆ Enter all required information. Multi-Dose Vial Lot Number, Assay (or Reference) Date and Time, Total Activity, RP Volume, and Saline Volume are all required. Refer to the "Entering RP Assay and Saline Information" section.
<Selected RP> prime failed. Check SAS and vial installation, then reprime.	◆ Confirm that the SAS is correctly installed. Refer to the "Recovering from Priming Issues" section. ◆ Replace the SAS if necessary.
Source not detected.	◆ During calibration, the source was not inserted into the Dose Calibrator. Insert the appropriate source to continue with calibration. Refer to the "Dose Calibrator Calibration" section of this manual.
The date and time must be set.	◆ Set the date/time. Refer to the "System Configuration" chapter.
The process has been interrupted and must be restarted.	The STOP button has been pressed. To resume operation, do one of the following: During SAS and PAS priming: ◆ Start the priming sequence again by pressing the PRIME button. During dose extraction: ◆ Press the DISCARD button and then press the PREPARE button to extract a new dose. During test Injection: ◆ Press the RESUME button to resume the test injection OR ◆ Press the STOP button to stop the test injection During patient infusion: ◆ Press the RESUME button to resume the infusion OR ◆ Press the DISCARD button to discard the dose.
Process paused.	During test Injection: ◆ Press the RESUME button to resume the test injection OR ◆ Press the STOP button to stop the test injection During patient infusion: ◆ Press the RESUME button to resume the infusion OR ◆ Press the DISCARD button to discard the dose.
The supplier activity is too high to perform the requested operation.	◆ Activity exceeds 25.9 GBq (700 mCi). ◆ Confirm that the assay information was entered correctly. Re-enter if needed. ◆ Confirm that the assay information provided by the pharmacist is correct. Re-enter if needed. ◆ If the assay was entered correctly, wait until the activity decays to 25.9 GBq (700 mCi).
This action is not allowed in Training Mode.	◆ Dose Calibrator calibration and the Linearity check can not be performed while the system is in Training Mode. Restart the system in Clinical Mode to perform a Dose Calibrator calibration or a Linearity check.

Error Message	Explanation
Unable to detect printer. Press OK to retry.	◆ Press the OK button and the system should correct the issue. ◆ If the message reappears, restart the system. ◆ If the message persists, record the error message and contact Bayer.
Vial activity exceeds safe shielding level. Immediately return vial to protective container.	◆ Activity exceeds 27.75 GBq (750 mCi), which is the maximum safe shielding level for the system. ◆ The Vial Shield should be removed from the system until the total activity in the Multi-Dose Vial decays to ◆ 27.75 GBq (750 mCi).
Waste activity will exceed shielding limit. Replace SAS.	◆ Activity in the Waste Container exceeds ◆ 3.7 GBq (100 mCi). ◆ Replace SAS before continuing operations.
Residual activity exceeds 2% of original dose. Press OK to flush to waste.	◆ High residual activity remains after infusion. ◆ Refer to the "PAS Occlusion Recovery" section.
PAS Occlusion Detected. WARNING PAS is pressurized and may leak if disconnected.	◆ If the infusion is stopped prematurely, do not disconnect disposables before discarding dose, as liquid leakage can occur. ◆ Reference "PAS Occlusion Recovery".

E Vials and Vial Shields

⚠WARNINGS

Patient injury may result from using vials with other than a 20 mm diameter cap.

- Use of vials having other than a 20 mm diameter cap may create a particulate hazard.

No modification of this equipment is allowed.

Biological contamination may occur if aseptic techniques are not followed.

- Follow facility guidelines to disinfect the vial shield.

Radiation Exposure Hazard.

- The Intego PET Infusion System Vial Shields must be transported in a container that complies with applicable transportation regulations.

⚠CAUTIONS

If using an absorbent disc, place no more than one disc in the vial shield before installing an RP vial.

- Additional absorbent discs may result in improper fit and potential damage to the vial.

If using an absorbent disc, do not expose the absorbent disc to cleaning solutions or other fluids.

- Discard and replace the absorbent disc if it has been exposed to any fluid. Failure to replace the absorbent disc that has been exposed to fluids may result in improper fit and potential damage to the vial.

If using an absorbent disc, discard and replace the absorbent disc after each use.

- Failure to replace the absorbent disc may result in improper fit and damage to the vial.

Do not place labels, tags, or other materials on the exterior of the vial shield surface.

- Application of labels, tags, or other materials may cause the vial shield to fit incorrectly in needle insertion device and result in improper fit and functioning of the system.

Ensure that the needle insertion guide is clean and free of debris.

- Debris may cause the vial shield to fit incorrectly in needle insertion device and result in improper fit and functioning of the system.


NOTE: Table D-1 represents the vials and compatible vial shields for use with the MEDRAD® Intego PET Infusion System. For additional vials not specified in this table, please contact Bayer.

Table E - 1: Vial Shield and Vial Specifications Summary

Vial Shield		Vial(s)				
Collar Color	Weight	Volume Range	ISO 8362-1:2003 ISO 8362-4:2003 Size Specification	Body Diameter	Glass Height	Cap Diameter
Silver	8.6 kg (18 lbs)	30-35 ml	30H	36 ± 0.5mm (1.417 ± 0.020”)	62.8 ± 0.7 mm (2.472 ± 0.027”)	20 mm (0.787”)
Blue	8.6 kg (18 lbs)	20-25 ml	20R	30 ± 0.25 mm (1.181 ± 0.009”)	55 ± 0.7 mm (2.165 ± 0.027”)	20 mm (0.787”)
Gold	8.6 kg (18 lbs)	20-25 ml	20H	32 ± 0.45mm (1.259 ± 0.017”)	58 ± 0.6 mm (2.283 ± 0.023”)	20 mm (0.787”)
Green	8.6 kg (18 lbs)	10-15 ml	10R	24 ± 0.2mm (0.944 ± 0.007”)	45 ± 0.5 mm (1.771 ± 0.020”)	20 mm (0.787”)
Brown	8.6 kg (18 lbs)	10-15 ml	15R	24 ± 0.25mm (0.944 ± 0.009”)	60 ± 0.5 mm (2.362 ± 0.020”)	20 mm (0.787”)
Black	8.6 kg (18 lbs)	10-15 ml	--	25 ± 0.25mm (0.984 ± 0.009”)	53.5 ± 0.5 mm (2.106 ± 0.020”)	20 mm (0.787”)
			10H	25.4 ± 0.4mm (1.000 ± 0.015”)		
			--	25.5 ± 0.25mm (1.003 ± 0.009”)		
Violet	8.6 kg (18 lbs)	10 ml	--	21.25 ± 0.2mm (0.837 ± 0.007”)	51 ± 0.5 mm (2.008 ± 0.020”)	20 mm (0.787”)
				21.75 ± 0.25mm (0.856 ± 0.009”)	52 ± 0.5 mm (2.047 ± 0.020”)	20 mm (0.787”)

Table E - 2: Vial Shield and Vial Specifications - Silver

Vial Shield		Vial(s)				
Collar Color	Weight	Volume Range	ISO 8362-1:2003 ISO 8362-4:2003 Size Specification	Body Diameter	Glass Height	Cap Diameter
Silver	8.6 kg (18 lbs)	30-35 ml	30H	36 ± 0.5mm (1.417 ± 0.020")	62.8 ± 0.7 mm (2.472 ± 0.027")	20 mm (0.787")



Table E - 3: Vial Shield and Vial Specifications - Blue

Vial Shield		Vial(s)				
Collar Color	Weight	Volume Range	ISO 8362-1:2003 ISO 8362-4:2003 Size Specification	Body Diameter	Glass Height	Cap Diameter
Blue	8.6 kg (18 lbs)	20-25 ml	30H	30 ± 0.25mm (1.181 ± 0.009”)	55 ± 0.7 mm (2.165 ± 0.027”)	20 mm (0.787”)



Table E - 4: Vial Shield and Vial Specifications - Gold

Vial Shield		Vial(s)				
Collar Color	Weight	Volume Range	ISO 8362-1:2003 ISO 8362-4:2003 Size Specification	Body Diameter	Glass Height	Cap Diameter
Gold	8.6 kg (18 lbs)	20-25 ml	20H	32 ± 0.45mm (1.259 ± 0.017”)	58 ± 0.6 mm (2.283 ± 0.023”)	20 mm (0.787”)



Table E - 5: Vial Shield and Vial Specifications - Green

Vial Shield		Vial(s)				
Collar Color	Weight	Volume Range	ISO 8362-1:2003 ISO 8362-4:2003 Size Specification	Body Diameter	Glass Height	Cap Diameter
Green	8.6 kg (18 lbs)	10-15 ml	10R	24 ± 0.2mm (0.944 ± 0.007")	45 ± 0.5 mm (1.771 ± 0.020")	20 mm (0.787")



Table E - 6: Vial Shield and Vial Specifications - Brown

Vial Shield		Vial(s)				
Collar Color	Weight	Volume Range	ISO 8362-1:2003 ISO 8362-4:2003 Size Specification	Body Diameter	Glass Height	Cap Diameter
Brown	8.6 kg (18 lbs)	10-15 ml	15R	24 ± 0.25mm (0.944 ± 0.009")	60 ± 0.5 mm (2.362 ± 0.020")	20 mm (0.787")


Table E - 7: Vial Shield and Vial Specifications - Black

Vial Shield		Vial(s)				
Collar Color	Weight	Volume Range	ISO 8362-1:2003 ISO 8362-4:2003 Size Specification	Body Diameter	Glass Height	Cap Diameter
Black	8.6 kg (18 lbs)	10-15 ml	--	25 ± 0.25mm (0.984 ± 0.009")	53.5 ± 0.5 mm (2.106 ± 0.020")	20 mm (0.787")
			10H	25.4 ± 0.4mm (1.000 ± 0.015")		
			--	25.5 ± 0.25mm (1.003 ± 0.009")		



Table E - 8: Vial Shield and Vial Specifications - Violet

Vial Shield		Vial(s)				
Collar Color	Weight	Volume Range	ISO 8362-1:2003 ISO 8362-4:2003 Size Specification	Body Diameter	Glass Height	Cap Diameter
Violet	8.6 kg (18 lbs)	10 ml	--	21.25 ± 0.2mm (0.837 ± 0.007”)	51 ± 0.5 mm (2.008 ± 0.020”)	20 mm (0.787”)
				21.75 ± 0.25mm (0.856 ± 0.009”)	52 ± 0.5 mm (2.047 ± 0.020”)	20 mm (0.787”)

F Components and Catalog Numbers

	Catalog Number
System	
MEDRAD® Intego PET Infusion System	INT SYS 200
Accessories	
Power Cord - North America	EPC 700 A
Power Cord - International	EPC 700 C
Power Cord - China	CPC 700
Power Cord - Brazil	BRZ-PWR-CRD
Power Cord - UK/Saudi Arabia	84575445
Calibration Source Holder	3017468
Roll - Printer Labels - Thermal	60937964
Pinch Valve Covers	INT PVC
Disposables	
Source Administration Set (SAS)	INT CSS
Patient Administration Set (PAS)	INT CPS
Manuals	
Multilingual CD (Includes Operation Manual and Quick SAS Installation Guide)	86473038
Service Manual - English Only	86465116
Vial Shields	
Silver	INT VS SLV
Blue	INT VS BLU
Gold	INT VS GLD
Green	INT VS GRN
Brown	INT VS BRN
Black	INT VS BLK
Violet	INT VS VLT

G Country-Specific Information

The following information is country-specific. The information within each section pertains only to that country.

G.1 Australia

Australian Communications and Media Authority (ACMA) Regulatory Compliance Marking (RCM)

G.2 New Zealand

R-NZ New Zealand Radio Label Regulatory Compliance Mark (RCM)

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G.8 Singapore

Complies with IDA Standards DA107248

Infocomm Development Authority of Singapore (IDA) Compliance Label

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OMAN - TRA
R/3560/16
D090024

Oman Compliance Label

A
Access Cap 19, 57, 77, 115
Activity Bar 67, 68
Activity Remaining Graph 70
Air Detector 6, 13, 15, 56, 77, 113, 114, 118
Air Detector Holder 15, 55
B
Battery Maintenance 96
Brake Release 26
C
Cable Storage 17, 26, 27, 28
Calibrating the Pinch Valves 93
Calibration Source 34, 46, 90, 118
Calibration Source Holder 19, 47, 87, 91, 92
Carrying Handle 19, 57, 77, 78
Caster Brakes 17, 26, 27, 28, 105
Catalog Numbers 3
Cleaning 34, 83, 109
Cleaning Guidelines 83, 84
Clinical Mode 39, 40, 80, 81, 120
Components 15, 16, 17, 18, 39, 50, 51, 56, 59, 84, 133
Components and Catalog Numbers 133
Configuration Settings 72
Configuring the System 29
Confluence 17, 54, 115
Critical Error Messages 117
D
Daily QC 24, 25, 34, 45, 46, 47, 79, 85
Device Lock 21, 22, 26, 35, 39
Discarding the Dose 69
Disposal of Equipment 97
Dose Calibrator 7, 9, 13, 15, 17, 19, 21, 25, 79, 85, 86, 87, 90, 92, 93, 97
Dose Calibrator Calibration 79, 85, 90, 104
Drive Controller 17, 27, 85
Drive Engage Switch 17, 27, 85
Drive Override 17, 26, 27, 28
Drive Speed Switch 6, 17, 27
E
Entering Assay Information 48
G
Geometry Check 93
I
INFUSE button 70
Installing the Multi-Dose Vial 48, 51, 57, 59
Installing the PAS 58
Installing the SAS 48, 50, 51, 57, 59

L
Linearity Check 79, 85, 86, 87, 88
Lock Code 22, 25, 35, 36, 39
Locking the Device 22
Locking the Shielded Chamber 23
Luer Connector 18, 59
M
Maintenance Schedule 84, 85
Maintenance Settings 34
Mobility- INT SYS 200 105
Moving the Intego PET Infusion System - INT SYS 200 26
N
Needle Cartridge 6, 15, 17, 53, 54, 77, 113
Needle Insertion Device 15, 53, 54, 57, 58
O
Occlusions 70, 71
On/Shutdown/Standby Button 16, 20, 21, 22
On/Shutdown/Standby button 4
One-Way Check Valve 14, 18
Operator Position 27, 28, 100
P
PAS Compartment 15, 16, 58, 59
PAS Compartment Cover 50, 51, 57, 58, 59
PAS Occlusion Recovery 116, 117
PAS Pinch Valve 15, 54, 55, 77, 114, 115
PAS Tubing 65, 68
Patient Information 23, 33, 43, 61, 62, 63, 71, 72
Patient Tube 17, 54, 55, 60, 77
Personalized Dose 66, 67
Post-Coil Tube 17, 54, 60
Power Switch 16, 20
Pre-Coil Tube 17, 54, 60
Preparing an RP Dose 68
Priming the PAS 63
Priming the SAS 59
Printer 16, 118, 119, 121
R
Recommended Maintenance Schedule 84
Recoverable Errors / Messages 118
Recovering a Dose Due to a System Failure 114
Recovering from Activation of Air Detector 113
Recovering from Priming Issues 113
Removable Prime Tube 59, 64
Replacing the Printer Paper 95
RP Assay 24, 48, 51, 57, 59, 60
RP Pump 5, 13, 15, 17, 56, 77, 113, 115

RP Tube 17, 54, 56, 60
S
Safety Features 13
Saline Disconnect Luer 17, 52, 113, 114
Saline Hook 7, 16, 50, 51
Saline Pump 5, 13, 16, 17, 52, 53, 77, 113, 114
Saline Spike 5, 17, 51
Saline Tube 15, 17, 52, 53, 60, 77
Saline Tube Holder 15, 53
SAS Coil 17, 54, 60, 78
SAS Components 57
SAS Confluence Holder 15, 53, 54
SAS Track 15, 17, 53, 54
Securing the Intego PET Infusion System 21
Shielded Chamber 7, 13, 14, 15, 16, 21, 23
Shielded Chamber Lid 9, 14, 15, 16, 26, 27, 47, 50, 57, 58, 59
Shielded Chamber Lid Latch 16, 27, 47, 51, 58, 77
Shut Down 9
Specifications 3, 99
Specifications, Electrical 102
Specifications, Environmental 101
Specifications, Fluid Delivery 103
Specifications, Mechanical 99
Standby 16, 21, 22, 26, 27
Status Indicators 9, 17
Stop Button 4, 16, 120
Swabbable Valve 15, 17, 18, 55, 59, 77, 115, 116
Swabbable Valve Holder 15, 55, 77
Symbols 4
Symbols and Icons 4
System Configuration 29
System Messages 117
System Overview 13
T
T-Connector 6, 15, 17, 54
T-Connector Holder 15, 54, 55
Test Inject 14, 32, 65, 114
Training Mode 40, 79, 80, 81
Troubleshooting 113
U
Unlocking the Device 22
USB Port 16, 34, 35, 41, 42, 45
V
Vial Cap 19
Vial Monitoring 73, 74
Vial Shield 13, 15, 19, 48, 49, 57, 58, 60

Vial Shield Compartment 15, 53, 56, 57, 58, 115
W
Waste Container 13, 15, 17, 24, 27, 54, 60, 69, 70, 105, 117, 119, 121
Waste Pinch Valve 15, 54, 55, 77
Waste Storage 15, 54, 78
Waste Tube 17, 54, 60, 77

Bayer reserves the right to modify the specifications and features describe herein or to discontinue any product or service identified in this publication at any time without prior notice or obligation. Please contact your authorized Bayer representative for the most current information.

All patient data that appear in this document are fictitious. No actual patient information is shown.

Bayer, the Bayer Cross, MEDRAD, and MEDRAD Intego are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries. Other trademarks and company names mentioned herein are properties of their respective owners and are used herein solely for informational purposes. No relationship or endorsement should be inferred or implied.

© 2014-2018 Bayer. This material may not be reproduced, displayed, modified or distributed without the express prior written consent of Bayer.

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To provide feedback or request support, please use the contact form provided on radiology.bayer.com/contact

CE 0086

Manufacturer
Bayer Medical Care Inc.
1 Bayer Drive
Indianola, PA 15051
U.S.A.
Phone: +1-412-767-2400
+1-800-633-7231
Fax: +1-412-767-4120

EC REP **Authorized European Representative**
Bayer Medical Care B.V.
Horsterweg 24
6199 Maastricht Airport
The Netherlands
Phone: +31(0)43-3585601
Fax: +31(0)43-3656598

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©2007-2019 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer.

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Bayer, the Bayer Cross, MEDRAD, and MEDRAD Intego are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries.

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A glossary of the symbols used on the MEDRAD® Intego administration sets can be found on page 2 of these Instructions for Use.

Bayer Medical Care Inc.
1 Bayer Drive
Indianola, PA 15051-0780
U.S.A.
Phone: +1 (412) 767-2400
+1 (800) 633-7231
Fax: +1 (412) 767-4120

INT CSS:
AvailMed S.A. De C.V.
C. Industrial Lt.001 Mz. 105
No. 20905 Int. A Col. Cd. Industrial
Tijuana, Baja California
Mexico 22444

EC REP **Bayer Medical Care B.V.**
Horsterweg 24
6199 AC Maastricht Airport
The Netherlands
Phone: +31 (0) 43-3585601
Fax: +31 (0) 43-3656598

Importer BMC B.V.:
Bayer Medical Care B.V.
Avenue Céramique 27
6221 KV Maastricht
The Netherlands

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ENGLISH

Introduction: Read the information contained within this document completely before installing or using the MEDRAD *Intego PET Infusion System* Administration Sets. Understanding this information will assist the user in safely operating the MEDRAD Intego System.

Intended Use: The MEDRAD *Intego PET Infusion System* is intended to deliver accurate doses of ^{18}F Fluorodeoxyglucose (FDG) and ^{18}F Sodium Fluoride (NaF) radiopharmaceuticals and commonly used flushing solutions to patients during molecular imaging (nuclear medicine) diagnostic procedures. The MEDRAD *Intego PET Infusion System* is also intended to provide effective radiation shielding to medical personnel from Fluorine-18 (^{18}F) radiation exposure during nuclear medicine diagnostic procedures.

Indications for Use: The MEDRAD *Intego PET Infusion System* is indicated for the administration of ^{18}F Fluorodeoxyglucose (FDG), ^{18}F Sodium Fluoride (NaF), and commonly used flushing solutions to patients during molecular imaging (nuclear medicine) procedures.
NOTE: MEDRAD Intego disposables labeled for ^{18}F FDG or ^{18}F NaF are intended for the delivery of NaF with software version CR10.1 and later. Please contact the Bayer representative for additional information regarding administration of ^{18}F NaF.

Restricted Sale: The United States Food and Drug Administration (FDA) restricts sale of this device to physicians or those with written authorization from a physician.

Important Safety Notice: This device is intended to be used by trained personnel with experience in diagnostic imaging studies.

Report any serious incident that has occurred in relation to this device to Bayer (radiology.bayer.com/contact) and to your local European competent authority (or, where applicable, to the appropriate regulatory authority of the country in which the incident has occurred).

⚠ Warnings

For devices labeled for single use, please note: This product is intended for single use only. Do not resterilize, reprocess or reuse. The disposable devices have been designed and validated for single use only. Re-use of the single use disposable devices pose risks of device failure and risks to the patient. Potential device failure includes significant component deterioration with extended use, component malfunction, and system failure. Potential risks to the patient include injury due to device malfunction or infection as the device has not been validated to be cleaned or re-sterilized.

Biological contamination can result from reusing the Patient Administration Set or failure to follow sterile technique. DO NOT reuse or attempt to resterilize. Properly discard disposable items after use. If there is any possibility that contamination may have occurred during set-up or use, disassemble and set-up a new sterile product.

Severe operator injury can occur if care is not taken when opening or closing the lids over the shielded chamber or the Patient Administration Set compartment. Be certain that all fingers are clear before opening or closing the shielded chamber lid or the Patient Administration Set compartment lid.

Only the fluid path of the set is sterile and non pyrogenic. Do not use in a sterile or aseptic area without proper precautions.

Air embolism or improper injection dosage may occur if saline line is damaged. Make certain not to pinch the saline tube within the saline pump as this may damage the administration set.

Biological contamination may occur if aseptic techniques are not followed. Follow facility guidelines to disinfect the vial septum.

Radiation Exposure Hazard: Injector may contain radioactive materials. Do not reach under cart while system contains activity. Do not remove access panels while system contains activity.

Do not use if sterile package is opened or damaged. Patient or operator injury may result if package is opened or damaged, or if damaged components are used. Visually inspect contents and package before each use.

Puncture hazard. The SAS contains sharps. To prevent injury, take care when handling the SAS to control the sharps. Dispose of the sharps per individual facility guidelines.

DO NOT use the Source Administration Set with more than one multi-dose vial. Using the Source Administration Set with more than one multi-dose vial could result in injury to patient or operator.

Radiation exposure hazard. Using the lid of the Shielded Chamber as a work space may expose the operator to excess radiation. Do not use the lid to the Shielded Chamber as a work space while there are radioactive materials present within the system. Using the lid as a work space or placing items on the lid while it is open may damage the lid, hampering its ability to provide an adequate radiation shield.

Radiation exposure hazard. Radiopharmaceuticals emit gamma radiation. Operator should exercise extreme caution when handling radiopharmaceuticals. Doses may exceed maximum permissible dose if appropriate precautions are not properly exercised. Follow all individual facility guidelines for the handling and disposal of radioactive materials. DO NOT open the Shielded Chamber without first taking radiation precautions as outlined in facility procedures.

⚠ Warnings

DO NOT use any radiopharmaceutical products other than ^{18}F Fluorodeoxyglucose (FDG) or ^{18}F Sodium Fluoride (NaF) with the MEDRAD Intego System. Use of radiopharmaceuticals other than ^{18}F FDG or ^{18}F NaF may damage the MEDRAD Intego System and may cause inaccurate dose measurement and/or excessive radioactive exposure.

Improper injection dosage or radiation spillage may result if the Shielded Chamber lid is closed onto the tubing set. Ensure that the tubing set and all administration set components are properly installed before closing the Shielded Chamber lid.

Patient infection may result from the use of non-sterile components. Maintain sterility of all disposable components using aseptic techniques.

The PVC tubing in the PAS contains DEHP, a chemical known to the State of California to cause birth defects or other reproductive harm.

Extremely high doses of phthalates such as DEHP have been reported to cause birth defects in research animals. Although limited scientific data is available on the effect of phthalates in humans, health authorities recommend that exposure to devices containing these compounds be minimized as much as possible for children and nursing or pregnant women.

DEHP is also a substance on the candidate list of Substances of Very High Concern (SVHC) according to Article 59 of REGULATION (EC) no. 1907/2006 concerning the Registration, Evaluation, Authorization and Restriction of Chemicals (REACH).
Substance Name: Bis (2-ethyl(hexyl)phthalate) (DEHP)
CAS Number: 117-81-7
Ec Number: 204-211-0

INSTALLING THE SOURCE ADMINISTRATION SET (SAS)

Always install a new SAS when installing a new multi-dose vial. Install the multi-dose vial and the SAS in the Shielded Chamber using the following procedure:
1. Release the latch securing the Chamber Lid and slide the lid towards the front of the cart to expose the Shielded Chamber.
2. Open the disposable tray and remove the SAS from the packaging.
3. Place the saline bag onto the hanger on the rear of the cart near the handle and insert spike at the end of the SAS into the saline bag.
4. Insert the Saline Tube into the Saline Pump so that the connector is between the pump and the saline bag, with the connector close to the pump.
5. Insert the tubing into the hold-down clip and route the Saline Tube along the SAS track, starting at the Saline Pump and ending at the Vial Shield Compartment where the Vial Shield will be installed.

NOTE: Make certain that the Saline Tube is centered on the pump rollers. The tubing must not fall below the rollers. When closing the pump door, make certain that the Saline Tube is centered in the pump door so that the tubing is not pinched.

6. Insert the Sharps Cartridge into the clip on the shield lid. The tube exiting the Sharps Cartridge should be towards the RP Pump.

NOTE: Do not remove the protective sheath from the needles of the Sharps Cartridge.

7. Insert Confluence into the SAS Confluence holder.
8. Verify dose calibrator is empty before installing coil.
9. Insert the SAS Coil into the Dose Calibrator so that the loose ends of the tubing are towards the top.

NOTE: Be sure that the SAS Coil is placed at the bottom of the Dose Calibrator.

10. Unroll the Waste Bag and insert it into the shielded Waste Storage.
11. Insert the T-connector into the retainer.
12. Install the section of tubing between the T-connector and the Waste Bag into the waste pinch valve.

NOTE: The pinch valve is opened by pushing on the button on the top of the valve. To install the tubing, slide the tubing into the valve while pushing down on the button.

13. Install the pinch valve cover on the waste pinch valve.
14. Install the section of tubing between the T-connector and the Swabbable Valve into the PAS pinch valve.

NOTE: The pinch valve is opened by pushing on the button on the top of the valve. To install the tubing, slide the tubing into the valve while pushing down on the button.

15. Install the pinch valve cover on the PAS pinch valve.
16. Insert the free end of the SAS (exiting the PAS pinch valve) into the Air Detector.
17. Using the Vial Shield Compartment handle, lift the shield lid straight up and then rotate the lid counter-clockwise.
18. The shield lid is hinged. Lift the lid to expose the Sharps Cartridge Clip.
19. Lower the shield lid over the Sharps Cartridge.
20. Insert the Vial Shield containing a full vial of the radiopharmaceutical indicated for use into the Vial Shield Compartment. DO NOT remove the vial access cap of the Vial Shield. Leave the carrying handle on the Vial Shield.

NOTE: The Vial Shield is keyed so that it can only be properly inserted into the Vial Shield Compartment with the Vial Shield Compartment posts resting in the keyed areas located next to the flat surfaces of the Vial Shield.

21. Remove the cap from the Vial Shield by pushing down on the carrying handle, rotating the vial access cap counter-clockwise, and then lifting the cap straight off.
22. Follow aseptic techniques as outlined in individual facility guidelines to disinfect the vial septum.
23. Insert the sharps into the multi-dose vial by pivoting the shield lid back over the well and pressing the needles straight down into the vial.
24. Insert the tubing from the Sharps Cartridge into the RP pump. Ensure that the tube collar rests immediately outside the RP pump, on the side closest to the Vial Shield Compartment. Close the pump lid.
25. Verify that all of the SAS is lying flat within the channel, that the Vial Shield Compartment handle is down and that the door of the RP pump is closed with the tubing completely installed. Close the chamber lid.

NOTE: Follow facility procedures to ensure that the radiopharmaceutical indicated for use is properly secured within the MEDRAD *Intego PET Infusion System*.
NOTE: The MEDRAD *Intego PET Infusion System* is designed to be used with standard locks. Bayer does not supply a lock with the system.

26. Enter multi-dose vial assay information by following the procedure outlined in the Entering Multi-dose Vial Assay Information section of the MEDRAD Intego Operation Manual.

INSTALLING THE PATIENT ADMINISTRATION SET (PAS)

Always install a new PAS for each patient. Install the PAS using the following procedure:
1. Open the shield cover over the PAS area of the Shielded Chamber.
2. Clean the Swabbable Valve at the distal end of the SAS following aseptic techniques described in individual facility guidelines.
3. Attach the PAS to the Swabbable Valve at the distal end of the SAS.
4. Depending on which side of the patient the MEDRAD *Intego PET Infusion System* is located, route the PAS out either the left or right side of the shielding.
5. Ensure that the PAS will not be pinched, and close the shield cover.

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

FONTS ON THIS VERSION HAVE ALL BEEN CONVERTED TO OUTLINE/CURVES

SPECIFICATIONS

1.0 REQUIREMENTS:

1.1 MATERIAL TYPE: N/A

1.2 ADHESIVE : N/A

1.3 COLOR

BASE: N/A

NOMENCLATURE: BLACK

1.4 FINISH: N/A

1.5 LABEL SIZE: 6.25” X 8.75”

1.6 PACKAGE TYPE:
N/A

2.0 SPECIAL CONSIDERATIONS:

3.0 CERTIFICATION: NOT REQUIRED

DRAWN BY: ELIZA SHEARER		DATE 11/14/18
DO NOT SCALE DRAWING	UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ON: 2 PL DECIMALS ± .01 3 PL DECIMALS ± .005 FRACTIONS 0" - 6" ±1/32 FRACTIONS 6" - UP ±1/16 ANGULAR ± 0°30'	All technical information, equipment design, and application data contained within this document is the CONFIDENTIAL property of Bayer HealthCare. No use or disclosure thereof may be made without written permission from Bayer HealthCare.

Bayer HealthCare

TITLE

Label, Tyvek, Intego, SAS, Flex, MC, WLD

SIZE	DWG NUMBER	REV
C	DN-258049	A

SCALE: FULL | DATE COPIED: | SHEET 1 OF 1