

Surgical Gown

Instructions for use

Intended use

This is a gown made from non-woven and worn by the operating room personnel. It may be fluid resistant or impervious to fluids. A surgical gown is used during surgical procedures to provide a sterile barrier between the patient and the clinician.

A surgical gown is designed to isolate potential contaminants (e.g., microorganisms, particulate matter) present on the operating room personnel from the surgical site.

The forearms of the surgeon are in frequent contact with the patient during surgery. An additional surgical sleeve can be added over the surgical gown to provide a more efficient barrier at the level of forearms.

These products are intended to be used in sterile conditions.

Direction for use


Wear the surgical gown such as normal garment. Put the arms into the surgical gown and arrange to fit to the body.

Warnings

There's NO known warning explicit on labelling.

Cautions

There's NO known warning explicit on labelling.



Shelf life

The products duration of use is the same as the duration of the surgery, in average few hours.

The shelf-life of the products is five (5) years.

Storage condition

General warehouse conditions, keep dry.

Material

Part Name	Material Description
OPS Advanced Surgical Gowns	SMS material (PP), with soft treatment
OPS Ultimate Surgical Gowns	Trilaminate material (3 layers of material: SPP+PE+SPP), breathable
Eclipse Surgical Gowns	Spunlace (blend of PET (45%) + woodpulp (55%))

Legal Manufacturer

Medline International France SAS

5, rue Charles Lindbergh Chateaubriant 44110 France