



**STADIOMETER 100
MANUAL BOOK**



445 mm

210 mm




Technical Data

Brand	: InnoQ
Product Code	: INQPMC101003
Model / Type	: Stadiometer 100
Measurement Range	: 0 - 200 cm
Weight	: 2.4 ± 0.5 kg
Dimension Installed	: 530 x 351 x 2142 mm
Dimension Folded	: 596 x 351 x 171 mm
Graduation	: 1 mm
User	: Child - Adult
Feature	: Foldable, Easy to carry

Warranty

The warranty period for this product is 24 months upon delivery or installation data. This warranty extends to defects in materials or workmanship and includes shipping cost. Warranty claims will be accepted if the original product packaging is complete, therefore all packaging must be maintained. Warranty claims cannot be accepted if the product has been disassembled by people who do not have the authority to do so. Contact local sales agents in case of warranty issues.

Service Center





www.idsmmed.com

**INFANTOMETER 100
MANUAL BOOK**


210 mm

445 mm


Install



Infantometer Hinge System




Hinges Open



Folded Hinges

1



2


Technical Data

Brand	: InnoQ
Product Code	: INQPMC101002
Model / Type	: Infantometer 100
Measurement Range	: 0 - 100 cm
Weight	: 1.6 ± 0.5 kg
Dimension (L x W x H)	
When Used	: 1108 x 322 x 117 mm
When Folded	: 586 x 322 x 117 mm
Graduation	: 1 mm
Used For	: Baby, Age 0 - 24 Months
Feature	: Foldable and Lightweight

Warranty

The warranty period for this product is 24 months upon delivery or installation data. This warranty extends to defects in materials or workmanship and includes shipping cost. Warranty claims will be accepted if the original product packaging is complete, therefore all packaging must be maintained. Warranty claims cannot be accepted if the product has been disassembled by people who do not have the authority to do so. Contact local sales agents in case of warranty issues.

Service Center



3