

HOGY Surgical Drape

Reuse Prohibited

Product name: Surgical Drape**Intended use:**

This Drape is intended for use by medical personnel to protect both the wearer and the patient from the transfer of microorganisms, fluids and blood. The Drape is appropriate for use in major procedures.

Performance and Standards:

The product complies with EN 13795-1:2019 Standard Performance

Warnings & Contraindications:

- 1) Do not use this product if you are allergic to polypropylene, polyester or cotton.
- 2) Do not use this product after the expiry date has passed.
- 3) Do not use this product if the packaging is opened, wet or visibly damaged.

Shelf Life

Three years from the date of manufacture.

Handling and Storage Conditions

Store in a dry, cool, ventilated environment away from flammable and corrosive materials. Protect from direct sunlight and moisture.

Transport conditions

During transport, avoid damage to the outer packaging, direct sunlight and rain, snow and other forms of precipitation.

Disposal Instructions:

Please dispose of Drape as clinically contaminated waste as per your organization's guidelines for waste disposal.

Sterilization:

Product is sterilized by EO.

<Contact information (etc.)>

[Manufacturer/distributor and manufacturer name (etc.)]

<Manufacturer/distributor>

PT. Hogy Medical
Address: Infinia Park No. 45 Blok A
57, Jalan Dr. Saharjo, Manggarai,
Tebet, Kota Jakarta Selatan,
Daerah Khusus Ibukota Jakarta
Phone: (021) 83797001

Department<Manufacturer>

HOGY MEDICAL CO.,LTD.