

金石型矽膠管組說明書

(English Version)

Silicone Tubing Assembly

Product No.: GH – RST – S

Product Description

- WITH 22mm OD or 22mm ID Connector
- Length: Various length available upon request.
- With bore straight or curving inner surface.
- Our silicone circuits are designed for use with ventilators.

Disinfection/Sterilization Reusable Silicone Tubing by commonly used methods:

- Boiling.
- Pasteurizing.
- Liquid Chemical Soaking.
- ETO Sterilizing.
- Autoclaving up to 134 degrees Centigrade

Manufacturer;

Cheen Hounq Enterprise Co., Ltd..

23, Alley 11, Lane 65, San Dreen St.,
Shulin (238) Taipei, TAIWAN R.O.C.

Sale Representative;

Golden Hygeia International Co., Ltd.

2FL-1, 58 NANKING EAST ROAD, SEC.5,
TAIPEI, TAIWAN R.O.C.

Manufacture Date :

Expiration Date :

LOT CODE :