

Reinforced Universal Pack I (PP/PE)

Instructions For Use

Product components

Product description	Components
INQGNS202002: Reinforced Universal Pack I (PP/PE)	Hand towels * 4 (60 cm x 50 cm) OP tape * 2 (50 cm x 9 cm) Reinforced side drapes with adhesive + pouch * 2 (117 cm x 200 cm) Reinforced bottom drape with adhesive * 1 (190 cm x 190 cm) Reinforced top drape with adhesive * 1 (240 cm x 150 cm) Table cover * 1 (180 cm x 210 cm) Mayo stand cover * 1 (75 cm x 145 cm) Suture bag * 1 (16 cm x 30 cm)

Intended user

Medical staff and patients

Intended Use

It protects the health of medical staff while isolating the risk of cross-infection of the substance brought by the former to the patient.

Direction for use

1. Please use the products immediately after opening package.
2. Products have been sterilized, do not use if the packing is torn or broken.
3. Put the products in clean, dry, ventilated area and far away from fire.

Labeling symbols



Shelf life

5 years

Sterilization

Sterilized by Ethylene Oxide

Legal manufacturer

Global SCP Medical Systems Pte. Ltd., Singapore