金石型 T 型引流器/集水瓶说明书

(English Version)

Water Trap Assembly

Reusable Type Product No.: GH-RWT-012
Single Use Type Product No.: GH-RWT-013

Product Description

- WITH 22mm OD x 22mm OD Wye
- Spring-loaded plunger maintains circuit integrity while bottle emptied.

Disinfection/Sterilization Reusable Water Trap by commonly used methods:

- Boiling.
- Pasteurizing.
- Liquid Chemical Soaking.
- ETO Sterilizing.
- Autoclaving up to 132 degrees Centigrade

Manufacturer;
Cheen Houng Enterprise Co., Ltd..
23, Alley 11, Lane 65, San Dreen St.,
Shulin (238) Taipei, TAIWAN R.O.C.

Sale Representative;
Golden Hygeia International Co., Ltd.
2FL-1, 58 NANKING EAST ROAD, SEC.5,
TAIPEI, TAIWAN R.O.C.

Manufacture Date：
Expiration Date：

LOT CODE :