

gloveon Innova

Latex Exam Gloves
Powder Free, Standard Cuff

Innova is a light weight powder free latex gloves with a unique textured surface for better grip in both wet and dry environment. It is a versatile glove that has high tensile strength and superb barrier protection properties which are suitable for any demanding task in any market segment.

NON STERILE | LATEX | 2 | CE | MD | EN ISO 374-1:2016/Type B KPT | EN ISO 374-5:2016 VIRUS | Tested in accordance with EN ISO 374-5:2016 Resistance to Bacteria and Fungi - pass Resistance to Viruses - pass | 1068-12/2019

Length (mm)	
	≥ 230
Thickness Measurements (mm)	
Palm (centre of Palm)	0.11 ± 0.03
Finger (13mm ± 3mm from tip)	0.13 ± 0.03

Physical Properties	Before Ageing	After Ageing
Tensile Strength (MPa)	≥ 18	≥ 16
Elongation (%)	≥ 650	≥ 500

Inspection Levels & AQL	Inspection Level	AQL
Watertightness	G1	1.5
Physical Dimensions	S2	4.0
Physical Properties	S2	4.0
Visual Inspection (Major)	S4	2.5
Visual Inspection (Minor)	S4	4.0
Particulate Residue	N = 5	≤ 2mg/glove
Protein Content	N = 3	≤ 50µg/g

*Full test reports on Sensitisation Testing, Primary Skin Irritation and Residual Proteins Testing are available upon request.

FEATURES

- Fully textured
- Powder free
- Made with natural rubber latex
- Protein level ≤ 50µg/g
- Ambidextrous
- Standard cuff
- Natural colour

PACKAGING

100 gloves per box (XS-XL)
10 boxes per carton

REGULATORY COMPLIANCE

ARTG 127794, FDA 510(k),
MDR (EU) 2017/745, REACH, EU 10/2011,
EC 1935/2004, EU 2016/425

STANDARDS

ASTM D5151, ASTM D6124, ASTM D3578,
ASTM D5712, EN ISO 374-1 (Type B),
EN 374 part 2, 4 & 5, EN 420, EN 1186,
EN 13130, CEN/TS 14234, EN 16523-1,
ISO 10993 part 5 & 10

MANUFACTURING ACCREDITATIONS

ISO 9001, ISO 13485, EN ISO 13485

Measured breakthrough time (minutes)	>10	>30	>60	>120	>240	>480
Permeation performance level	1	2	3	4	5	6

Chemical	EN 16523-1:2015 Permeation Level	EN 374-4:2013 Mean Degradation (%)
K 40% Sodium Hydroxide	3	-24.7
P 30% Hydrogen Peroxide	4	19.8
T 37% Formaldehyde	5	-23.1

Hartalega

Hartalega Sdn Bhd
T: +603 6277 1733
W: hartalega.com.my/contact-us/

Protection Always On
Visit us at
gloveonglobal.com

Distributed by:

IDS Medical Systems (Thailand) Co., Ltd.
Lasalle Tower, Floor G/2, Floor 1-3, 10/11 Moo 16 Srinakarin Road, Bangkaew, Bangphi, Samutprakarn 10540, Thailand
Tel: +66 2349 4780 Fax: +66 2758 8883
E-mail: salesth@idsmed.com